

## Cover Sheet for

# ENVIRONMENTAL CHEMISTRY METHOD

***Pesticide Name:*** Diazinon

***MRID #:*** 414904-01

***Matrix:*** Water

***Analysis:*** GC/NPD

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at [flynt.elizabeth@epa.gov](mailto:flynt.elizabeth@epa.gov).

TO THE HONORABLE SPEAKER OF THE HOUSE OF REPRESENTATIVES

10-48943 10-48944

07-1053      98-10-24

10-10-68

SPRINGBORN LABORATORIES, INC.  
ENVIRONMENTAL SCIENCES DIVISION  
CHEMISTRY DEPARTMENT

DIAZINON  
414904-01

**TITLE:** A Gas-Liquid Chromatographic Method for the Measurement of Diazinon Residue in Pond Water

**AUTHOR:** T.Z. Kendall

**ABSTRACT:** A gas-liquid chromatographic (GC) method is described for the analysis of diazinon in pond water. The method employs a dichloromethane partition of the aqueous sample prior to analysis by GC utilizing nitrogen-specific detection. Recoveries of diazinon from pond water fortified in the range of 0.100 - 100 ppb resulted in a mean recovery of 105% with a recovery range of 82.3-140%. Control values were all < 0.0439 ppb.

**DATE:** July 18, 1989

## A Gas-Liquid Chromatographic Method for Measurement of Diazinon Residue in Pond Water

### INTRODUCTION

A method for determining diazinon [ O,O-Diethyl O-(2-isopropyl-6-methyl-4-pyrimidinyl) phosphorothioate ] in pond water is presented.

### PRINCIPLE AND APPLICATION

Diazinon is extracted from pond water (to which is added a saturated solution of sodium chloride) with dichloromethane which is subsequently dried through sodium sulfate, evaporated to dryness and dissolved in hexane. A gas-liquid chromatograph equipped with a thermionic detector operated in the nitrogen-specific mode measures the extracted residue and integrates the resulting peak. Linear regression analysis of diazinon peak areas for samples and reference standards permits calculation of each diazinon sample concentration.

### ANALYTICAL METHOD

#### Reagents

Dichloromethane, Fisher Optima, Reagent grade

Water, Burdick and Jackson, HPLC grade

Hexane, Burdick and Jackson, HPLC grade

Acetone, Burdick and Jackson, HPLC grade

Sodium chloride, Mallinckrodt, Reagent, A.C.S.

Sodium sulfate, Doe and Ingalls, Inc., Analytical Reagent

Diazinon, Lot No. S87-1185-1, 96.7% a.i., supplied by Ciba Geigy

#### Equipment

Balance, Ohaus Galaxy 160, four-place analytical balance

Flasks, volumetric, assorted sizes

Pipets, volumetric, assorted sizes

Separatory funnels, 500 mL

Glass funnels, assorted sizes

Glass wool, prewashed with dichloromethane

Roundbottom flasks, 500 mL

GC vials with crimp caps and Teflon-faced septa

Serum bottles, Wheaton, assorted sizes, with Teflon-lined lids and metal crimp caps

Syringes, Hamilton, assorted sizes

Rotary evaporator, Buchii Model R110, with vacuum pump, 40 °C water bath

## Detailed Procedure

### I. Preparation of Stock Solution

#### A. Diazinon (1 mg/mL)

1. Weigh 100 milligrams (a.L.) of diazinon on an analytical balance.
2. Transfer the diazinon to a 100-mL volumetric flask and dissolve to the mark with acetone.
3. Transfer the stock solution to a 100-mL amber serum vial and seal with a Teflon-lined crimp cap.
4. Store this stock solution in a refrigerator maintained at 4 °C.

### II. Control Sample Fortification

- A. Rinse all glassware with dichloromethane prior to fortification.
- B. To each separatory funnel, add 250 mL of pond water.
- C. Fortify each sample with diazinon by volumetric addition of dilutions of the primary stock solutions.

NOTE: The fortification levels produced in the control pond water samples for the method validation/recovery were 100, 10.0, 1.00, 0.500 and 0.100 ppb (three replicates at each level). An additional three pond water samples were left unfortified and utilized as control samples.

### III. Extraction

#### A. Pond Water

1. Add 10 mL of an aqueous saturated sodium chloride solution to the sample.
2. Add 100 mL of dichloromethane to the 250 mL aqueous sample.
3. Shake the separatory funnel for 90 seconds and allow the phases to separate.
4. Drain the lower dichloromethane layer through approximately 30 grams of sodium sulfate (glass wool plug in glass funnel) into a 500-mL roundbottom flask.
5. Repeat steps 1-4, combining the dichloromethane extracts in the roundbottom flask.
6. Rinse the sodium sulfate with approximately 25 mL of dichloromethane and collect the rinse in the same flask.
7. Evaporate the sample to approximately 0.5 mL on a rotary evaporator.
8. Evaporate the remaining dichloromethane under a gentle stream of nitrogen.
9. Volumetrically pipet the requisite volume of reagent-grade hexane into the flask and swirl the flask in order to ensure complete dissolution of the diazinon residues.

#### NOTE:

- a) For the 0.100 ppb recovery samples a volume of 1 mL would yield a concentration of 25.0 µg/L.
  - b) The volume added for samples of unknown concentration would be such that the final concentration falls in the midrange of the standard curve.
10. Remove an aliquot of each solution and proceed to Section IV, Gas-Liquid Chromatography.

### IV. Gas-Liquid Chromatography

A. Instrumental Conditions: Hewlett-Packard Model 5840 gas-liquid chromatograph equipped with Hewlett-Packard Model 7671A autosampler and thermionic detector.

Column: DB-17. 15 M x 0.53 mm ID with 1.0 µm film thickness (J & W Cat #125-1712)

Column Flowrate: Helium @ 20 mL/minute

Detector: Hydrogen @ 3.3-3.5 mL/minute

Air: @ 80-90 mL/minute  
Injector Temperature: 250 °C  
Column Temperature: 170 °C (isothermal)  
Detector Temperature: N-P @ 275 °C  
Chartspeed: 0.3 cm/minute  
Injection Volume: 4 µL  
Attenuation: 2<sup>5</sup>  
Slope Sensitivity: 0.20

### C. Analysis

1. Prepare standard solutions containing diazinon. Standard solution concentrations used for the recovery study were 25.0, 50.0, 100, 300 and 500 µg/L.
2. Inject 4 µL of the 25.0 µg/L standard solution. Adjust the attenuation so that the peak signal results in at least a fifteen percent deflection from the baseline.
3. Inject 4 µL of each of the standards, document the peak areas, and determine the correlation coefficient of the line. Proceed to step 4 if the correlation coefficient is greater than or equal to 0.985.
4. Inject 4 µL of several samples.
5. After each set of samples, reinject 4 µL of each of the standards and document the peak areas.
6. Construct a standard curve for the analyte (using all standard results) by plotting peak area observed versus the concentration (µg/L) of the standard injected.
7. The standard linear regression analysis for diazinon is used to determine the concentration in each sample.
8. In order to determine the analytical result for each sample, the following equation is used:

$$\text{Analytical Result (ppb)} = A \times \text{D.F.} / d,$$

where:

Analytical Result = concentration of diazinon

A = concentration ( $\mu\text{g/L}$ ) of sample from the regression analysis  
D.F. = dilution factor, ratio of the final volume (mL) of the sample to the initial volume (mL) of sample used  
 $d_w$  = density of water (kg/L)

## RESULTS AND DISCUSSION

Pond water samples were prepared which contained a variety of concentrations of diazinon. To each sample was added 10 mL of an aqueous saturated solution of NaCl. Each sample was then partitioned into chloromethane, dried through sodium sulfate, rotary evaporated to approximately 0.5 mL, evaporated to dryness under nitrogen and dissolved in GC-grade hexane.

Recoveries from pond water fortified with diazinon at 100 ppb ranged from 92.7 to 98.7% with an average recovery of 95.6%; recoveries at 10.0 ppb ranged from 101 to 104% with an average recovery of 102%; recoveries at 1.00 ppb ranged from 82.3 to 103% with an average recovery of 92.8%; recoveries at 0.500 ppb ranged from 106 to 113% with an average recovery of 109%; recoveries at 0.100 ppb ranged from 117 to 140% with an average recovery of 126%. The minimum detectable level for diazinon was determined by linear regression analysis of the calculated value of one-half the mean of the peak areas of the lowest standard, resulting in an M.D.L. of less than 0.0439 ppb.

Linear regression analyses for the response of diazinon is shown in Figure 16A. Representative chromatograms of a diazinon standard, of an extracted diazinon pond water sample and a pond water control sample are shown in Figure 17A.




Table 19A. Analytical results for the recovery of diazinon from pond water.

| Fortified Level (ppb) | Sample Volume (mL) | Level Found (ppb) | Percent Recovery (%) |
|-----------------------|--------------------|-------------------|----------------------|
| 100                   | 250                | 92.7              | 92.7                 |
| 100                   | 250                | 95.4              | 95.4                 |
| 100                   | 250                | 98.7              | 98.7                 |
| 10.0                  | 250                | 10.1              | 101                  |
| 10.0                  | 250                | 10.1              | 101                  |
| 10.0                  | 250                | 10.4              | 104                  |
| 1.00                  | 250                | 0.823             | 82.3                 |
| 1.00                  | 250                | 0.929             | 92.9                 |
| 1.00                  | 250                | 1.03              | 103                  |
| 0.500                 | 250                | 0.563             | 113                  |
| 0.500                 | 250                | 0.537             | 107                  |
| 0.500                 | 250                | 0.530             | 106                  |
| 0.100                 | 250                | 0.117             | 117                  |
| 0.100                 | 250                | 0.120             | 120                  |
| 0.100                 | 250                | 0.140             | 140                  |
| Control               | 250                | < 0.0439          | NA                   |
| Control               | 250                | < 0.0439          | NA                   |
| Control               | 250                | < 0.0439          | NA                   |

Average Recovery:  $105.5\% \pm 13.7$ .

Theoretical minimum detectable concentration is &lt; 0.0439 ppb for a 250-mL sample.

Figure 16A. The linear regression analysis for diazinon standards used in the pond water recovery study.



Springborn Laboratories, Inc.

Figure 17A. Representative GC chromatograms of a diazinon standard (upper), of an extracted pond water sample (middle) and a control pond water sample (lower).




HP RUN # 19 JUN/22/89 TIME 18:26:44  
BOTTLE 19  
AREA %

| RT   | AREA   | AREA %  |
|------|--------|---------|
| 4.38 | 159388 | 100.000 |

Standard  
500 µg/L

DIL FACTOR: 1.0000 E+0




HP RUN # 23 JUN/22/89 TIME 18:58:13  
BOTTLE 23  
AREA %

| RT   | AREA  | AREA % |
|------|-------|--------|
| 1.53 | 215   | 1.284  |
| 3.98 | 1727  | 10.314 |
| 4.84 | 1864  | 11.132 |
| 4.37 | 18298 | 61.451 |
| 5.36 | 2649  | 15.828 |

Recovery  
0.100 µg/L  
A µg/ml

DIL FACTOR: 1.0000 E+0



HP RUN # 27 JUN/22/89 TIME 19:31:09  
BOTTLE 27  
AREA %

| RT   | AREA | AREA % |
|------|------|--------|
| 3.42 | 1456 | 18.383 |
| 3.98 | 2244 | 16.773 |
| 4.38 | 1489 | 18.531 |
| 4.48 | 2786 | 28.226 |
| 5.11 | 1636 | 12.228 |
| 5.35 | 3928 | 29.359 |

Recovery  
Control  
0

DIL FACTOR: 1.0000 E+0

# 6.0 APPENDIX III - METHOD VALIDATION: DIAZINON IN SEDIMENT

NO. 1000

ANALYST: J. L. HARRIS  
 DATE: 10/10/90  
 METHOD: EPA 821-R  
 SAMPLE NO.: 1000  
 ANALYST: J. L. HARRIS

ANALYST: J. L. HARRIS  
 DATE: 10/10/90  
 METHOD: EPA 821-R  
 SAMPLE NO.: 1000  
 ANALYST: J. L. HARRIS

ANALYST: J. L. HARRIS  
 DATE: 10/10/90  
 METHOD: EPA 821-R  
 SAMPLE NO.: 1000  
 ANALYST: J. L. HARRIS