

Cover Sheet for

ENVIRONMENTAL CHEMISTRY METHOD

Pesticide Name: Prometryn

MRID #: 415464-01

Matrix: Soil

Analysis: GC/FPD

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at flynt.elizabeth@epa.gov.

ORGANICS LABORATORY

Minnesota Valley Testing Laboratories, Inc.
1126 North Front Street
New Ulm, MN 56073

Method No.	Edition	Revision
CHR-12111671	03/28/89	1.0
Subject:		
Determination of Prometryn and GS-11354 in Soils		
Submitted By: Daron Decker		
Approved By: Daron Decker		
References: CIBA-GEIGY Corporation		

1.0 SCOPE

This method describes the procedure for extraction, cleanup and quantitation of Prometryn and GS-11354 residues in soil. This method is sensitive to 0.01 ppm.

2.0 PRINCIPLE

Prometryn and GS-11354 are extracted from soil by refluxing with acetonitrile/water (80:20). The extract is concentrated to aqueous and cleaned up by use of an Extrelut column. The analytes are eluted with ethyl acetate/hexane (25:75) and the eluate is dried by rotary evaporation. The analytes are brought up in ethyl acetate/methanol (95:5) and analyzed by capillary gas chromatography with flame thermionic detection.

3.0 CHEMICALS AND SOLUTIONS

Prometryn Standard Reference Material (CIBA-GEIGY Corporation)
GS-11354 Standard Reference Material (CIBA-GEIGY Corporation)
Acetonitrile
Ethyl Acetate

3.0 Chemicals and Solutions (continued)

Glass Wool
Hexane
Methanol
Water, Reverse Osmosis/Deionized
Extrelut

4.0 APPARATUS

4.1 Equipment

Heating Mantle or Hot Plate
Glass Microfibre Filter
Buchner Funnel
Teflon Boiling Chips
Nitrogen Atmosphere Evaporator
Extrelut Reservoir Tube, 75 mL capacity

4.2 Glassware

Graduated Cylinder, 200 mL
Boiling Flask, 500 mL
Liebig Condenser
Erlenmeyer Sidearm Suction Flask, 500 mL
Pipets, Class A Volumetric
Pipettes, Dispo Pastuer
Vial, Solvent Saver (7 mL)
Vial, Auto Sampler (1.5 mL) with screw cap

4.3 Gas Chromatograph

Instrument: Shimadzu GC-9A with Flame
Thermionic Detector
Autosampler: Shimadzu AOC-9
Column: 15 x 0.53 mm ID, 5% Phenyl methyl
silicon, 1 μ m DF
Data System: Maxima Data System, Dynamic
Solutions or equivalent

5.0 ANALYTICAL PROCEDURE

5.1 Extraction

- 5.1.1 Weigh 50 g (to 0.001g) of soil into a 500 mL boiling flask and record the weight on the Raw Data Summary Sheet.
- 5.1.2 Add 250 mL of acetonitrile:water (80:20) and assemble the flask, water cooled Liebig condenser, and heating mantle for refluxing.

5.1 Extraction (continued)

5.1.3 Reflux for one hour.

5.1.4 Remove heat from the flask and allow it to cool to approximately room temperature.

5.1.5 Filter sample solution through Buchner Funnel filter with a microfibre filter. Rinse boiling flask and filtrate with approximately 50 mL of 80/20 acetonitrile/water. Collect all of the filtered extract.

5.1.6 Quantitatively transfer the filtered extract to a 500 mL boiling flask and rotary evaporate until approximately 20 mL aqueous remains.

5.1.7 Add one package of Extrelut to the sample flask and mix well. Allow the sample to adsorb onto the Extrelut for a minimum of 20 minutes.

5.1.8 Place a plug of glass wool in the tip of a polyethylene Extrelut reservoir tube and transfer the adsorbed sample.

5.1.9 Elute the analytes of concern into a 500 mL boiling flask with 225 mL of ethyl acetate:hexane (25:75).

5.1.10 Rotary evaporate the sample to dryness and bring up with approximately 5 mL of ethyl acetate:methanol (95:5).

5.1.11 Transfer the sample to a saver vial and dry on a nitrogen evaporator.

5.1.12 Bring up in 2 mL of ethyl acetate:methanol (95:5) and transfer to an autosampler vial for analysis by capillary gas chromatography.

6.0 GAS CHROMATOGRAPHIC ANALYSIS

6.1 Instrument Set Points

- 6.1.1 Column: 15 m x 0.53 mm ID, 5% Phenyl methyl silicon capillary column, 1 um DF
- 6.1.2 Temperature: Column oven-150 degrees C
Injector/Detector Block-225 degrees C

6.1 Instrument Set Points (continued)

6.1.3 Times: Run Time-12 minutes
Approximate Elution Times:
GS-11354 - 5.6 minutes
Prometryn - 9.0 minutes

6.1.4 Gases: Carrier - Helium at
0.8kg/cm2
Detector - Air and Hydrogen
at appropriate setting to
achieve optimum sensitivity

6.1.5 Miscellaneous: Range 2
Injection Volume - 5 uL

6.2 Calibration and Quantitation

6.2.1 Set up instrument and detector as per
Shimadzu GC-9A I and II Instruction
Manuals. Verify performance and
stability as per manuals.

6.2.2 Inject standards and samples in a ratio
of at least one standard for every three
samples. Be sure standards cover the
entire concentration range of the
samples. Have at least four different
concentrations of standards across the
concentration range.

6.2.3 Use standard peak heights and
concentrations to calculate the linear
regression and interpolate residue
concentration in sample extracts.

6.2.4 Calculation:

$$\text{ppm} = \frac{\text{ng residue found}}{\text{mg sample extracted}}$$