

Cover Sheet for

ENVIRONMENTAL CHEMISTRY METHOD

Pesticide Name: Prometryn

MRID #: 415464-01

Matrix: Soil

Analysis: LC/MS/MS

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at flynt.elizabeth@epa.gov.

Appendix IX

Enseco Thermospray LC/MS/MS Confirmations

Standard Operating Procedure:

"Determination of Prometryn, GS-11354 and GS-11526
in Soil by LC/MS/MS"

Mass Spectra of Prometryn and its metabolites
by LC/MS and LC/MS/MS

STANDARD
OPERATING
PROCEDURE

Subject or Title:
Determination of Prometryn,
GS-11354 and GS-11526 in Soil by LC/MS/MS

Page 1 of 5

SOP No.:
LM-CAL-3045

Revision No.:
Original

Effective Date:
April 23, 1990

Supersedes: None

1. PRINCIPLE

Prometryn and metabolites (GS-11354 and GS-11526) are extracted from soil using a 80:20 (v/v) acetonitrile-water solution. Sample extracts are centrifuged, reduced under nitrogen, and then analyzed by TSP-LC/MS/MS. The detection limit for this method is 10 ug/kg.

2. STANDARDS

2.1 Standards

2.1.1 Prepare a 1 mg/mL individual stock solutions of each analyte in HPLC grade methanol. Replace stock solution every four months or sooner.

2.1.2 Using the 1 mg/mL individual stock solutions from 2.1.1 prepare a 10 ug/mL mixed standard in methanol. Replace mixed analyte solution every four months or sooner.

2.1.3 Using the 10 ug/mL mixed analytical standard solution from 2.1.2 prepare 0.10 ug/mL, 1.0 ug/mL and 10 ug/mL fortification standards in methanol. Replace spiking standards every four months or sooner.

2.1.4 Using the 10 ug/mL mixed standard solution from 2.1.2 prepare calibration standards at 5 ng/mL, 25 ng/mL, 100 ng/mL, 250 ng/mL, and 500 ng/mL in 50/50 Meoh/H₂O. Replace calibration standards every four months or sooner.

Prepared By:
Frank Kenney

Date:
April 23, 1990

Management Approval:

Date:

QA Officer Approval:

Date:

Subject or Title:
Determination of Prometryn,
GS-11354 and GS-11526 in Soil by LC/MS/MS

Page 2 of 5

SOP No.:
LN-CAL-3045

Revision No.:
Original

Effective Date:
April 23, 1990

3. PROCEDURE

3.1 Extraction

3.1.1 Weigh 10.0 grams of soil into a 40-ml vial.

3.1.2 If applicable, prepare fortification samples by adding fortification standard.

3.1.2.1 Add 1.0 ml of the 0.10 ug/ml fortification standard to yield 10.0 ug/kg (ppb).

3.1.2.2 Add 1.0 ml of the 1.0 ug/ml fortification standard to yield 100.0 ug/kg (ppb).

3.1.2.3 Add 1.0 ml of the 10.0 ug/ml fortification standard to yield 1,000 ug/kg (ppb).

3.1.3 Add the acetonitrile-water solution (80:20 v/v).

3.1.3.1 Add 20.0 ml to the sample and method blank.

3.1.3.2 Add 19.0 ml to the fortified samples to give a total volume of 20.0 ml of the extracting solution.

3.1.4 Adjust the "pH" of the extracting solution to 7.2 - 7.8 with acetic acid or NH_4OH .

3.1.5 Shake the extracts on the orbital shake for 3 hours at ca 240 rpm with the vials positioned horizontally.

3.1.6 Centrifuge the extracts up to 1000 rpm for about 5 minutes to obtain a clear supernatant.

3.1.7 Aliquot 5.0 ml of the supernatant and filter it through a 0.45 μm syringe filter into an 8-ml test tube. Filter an additional 0.5 ml of 80:20 acetonitrile-water through the filter into the test tube.

3.1.8 Reduce under nitrogen to less than 1 ml to remove the acetonitrile.

STANDARD
OPERATING
PROCEDURE

Subject or Title:
Determination of Prometryn,
GS-11354 and GS-11526 in Soil by LC/MS/MS

Page 3 of 5

SOP No.:
LN-CAL-3045

Revision No.:
Original

Effective Date:
April 23, 1990

3.1.9 Adjust the final volume to 1.0 mL with methanol.

3.1.10 Transfer the extract to a 4 mL vial. Final concentration of the sample is 2.5 g/1.0 mL.

3.1.11 Organize the sample vials into the vial box, and store it in the refrigerator at 2-6° C until analysis.

4. LIQUID CHROMATOGRAPH OPERATING PARAMETERS

4.1 The LC chromatography system parameters have been optimized for these analytes as follows:

Column: 1) JIP Hypersil ODS 6 cm x 4.6 mm ID (3µm)
No Guard or equivalent

Flow rate: 0.8 mL/min. + 0.4 mL/min.
Post column addition of Buffer

Mobile Phase: ACN/Water
(Water contains 0.1M ammonium acetate)

Gradient: 30/70 ACN/H₂O
Gradient to 65/35 at 2.5 minutes
to 100% ACN at 4 min.
Return to 30/70 at 6.5 minutes
or equivalent conditions which provide
baseline resolution of all analytes.

Retention Times: Prometryn 7:20 minutes (approx.)
GS-11354 2:40 minutes (approx.)
GS-11526 1:50 minutes (approx.)

5. TSP - MASS SPECTROMETER OPERATING PROCEDURES

5.1 Instrument Tuning. - The instrument must be capable of detecting 100 pg of each analyte, (20 µL of 5 ng/mL) with a signal/noise (S/N) ratio of 5:1.

5.1.1 If desired, a mixture of prometryn and the metabolites (1 mg/L) may be used as a tuning solution. This provides fine tuning of the TSP interface parameters.

Subject or Title:
Determination of Prometryn,
GS-11354 and GS-11526 In Soil by LC/MS/MS

Page 4 of 5

SOP No.:
LN-CAL-3045

Revision No.:
Original

Effective Date:
April 23, 1990

5.2 TSP-LCMS Operation Parameters.

5.2.1 The following TSP-LCMS parameters must be optimized in order to achieve an optimum performance.

Vaporizer Control Point - 75 °C - 120 °C
Aerosol Temp - 250 °C - 320 °C
Repeller Voltage - 40 V - 75 V

5.2.2 It has been determined that optimal sensitivity is achieved with filament on ionization mode, however, either mode is acceptable.

5.3 Selected Ion Monitoring.

5.3.1 Samples and standards products are analyzed by selected reaction monitoring (SRM) of products derived from the protonated molecular ions (MH)⁺. OS - Atrazine is included as an internal standard for retention time verification and possible future quantitation.

5.3.2 The appropriate masses to be monitored are listed below:

Analyte	Precursor (MH ⁺)	Product Ion	Dwell Time
Prometryn	m/z 242	m/z 158 ± .3 amu	0.3 seconds/window
GS-11354	m/z 200	m/z 158 ± .3 amu	0.3 seconds/window
GS-11526	m/z 212	m/z 128 ± .3 amu	0.3 seconds/window
OS-Atrazine	m/z 221	m/z 179 ± .3 amu	0.3 seconds/window

The seconds/window time may be changed providing that at least 6 scans are acquired per peak.

5.4 Calibration and Sample Analysis.

5.4.1 During the course of analysis of a batch of samples each of the calibration standards (see Section 2.1.4) is run once. Between 1 and 6 samples is analyzed between these standards.

5.4.2 Calculate the calibration factor (CF) for each standard, defined as the ratio of the standard concentration to the response for each analyte at each standard concentration.

Subject or Title:
Determination of Prometryn,
GS-11354 and GS-11526 in Soil by LC/MS/MS

Page 5 of 5

SOP No.:
LM-CAL-3045

Revision No.:
Original

Effective Date:
April 23, 1990

5.4.3 If the percent relative standard deviation is less than 20% over the working range, linearity can be assumed and the average calibration factor can be used in place of a calibration curve.

$$5.4.3.1 \text{ Average Calibration Factor} = \frac{\text{Concentration (ng/mL)}}{\text{Total Area Counts}} \cdot n$$

Where n = number of calibration injections.

5.4.4 Inject 20-50 µL of extract into the TSP-LCMS system for analysis.


6. CALCULATIONS

6.1 The concentration of each analyte in a sample is calculated as follows:

$$\text{Concentration (ug/Kg, dry weight)} = \frac{(CF) (Vx) (Ax)}{(Ws) (Rw) (1000)}$$


CF = Average Calibration factor (ng/mL)
Vx = Sample extract final volume (mL)
Rw = Dry Weight/Wet Weight Ratio
Ax = Sample Total Area Counts
Ws = Sample Wet Weight (Kg)

TSP LC/MS of Prometryn and Two Degradation Products




20 ul of 0.5 ppb Spike

TSP LC/MS/MS of Prometryn and Two Degradation Products




CHR0: PD9003281 ver 2 on UIC 002002 26-NHR-90 Elapset: 00:05:40.9 241 PROMETHEO + METAS
 Sample: SS-9003108 10HC/IL Start: 12:49:20 363 BY LC/ms/ms
 Coasm: 20UL ACDAUFROM V=39 S=280 EH=120088 FILLOFF CID=2.0MTORR
 Model: TSP + DRU LHR SMITH GAS UP LR
 Gases: FHR
 Peak: 1000.00 amu Label under 1.0 363 Initial: 123 179
 Area: 6, 4.00 Baseline: 15, 3 Masses: 123 179
 Label: -2, 10.00
 Sim traces



CHRO: PD980328J ver 2 on UIC 002002 20-MAR-90 Elapsed: 00:05:40.9 241
 Sample: SS-9003108 1M/C/NL Start: 12:49:20 363
 Cont: 20UL ACIDURPRON V-90 S-200 EH-1200XB FILOFF CID-2.0RTORR
 Model: TSP + DBU LMR SITHH GAS UP LR
 Oper: FHK
 Peak: 1000.00 auu Label: undet 1 > 363 Inlet: 1
 Area: 0, 4.00 Baseline: 15, 3 Masses: 128 > 179
 Label: -2, 10.00

Prometryn and Metabolites
 by LC/MS/MS
 Slimey Traces



QUAN: TEMP.09 ver 1 on UIC 002002

Entries: 4

Data file(s) used

File: PD900328J ver 2 on UIC 002002 28-MAR-90

Samp: SS-9003108 10MG/ML

Start : 12:49:20 338

Comm: 20UL ACDAUPROM V=90 S=280 EM=1200XB FILOFF CID=2.0HTORR

Mode: TSP +DAU LMR SYNTH GAS UP LR

Oper: FHM

Inlet :

Disp: SEAR CHRO


Amount = (area * ref amt) / (ref area * response factor)

Response factors from quantitation library

#	Compound Name	Time	MH+	Area	Amount
1	DS-ATRAZINE	5:32	179	686323.	50.00
2	PROMETRYNE	7:28	138	149523.	10.00
3	CS11526	1:43	128	72230.	10.00
4	CS11354	2:36	138	226703.	10.00

#	Compound Name	RF	Scan	Pred Time	Delta Time	Height	S/N
1	DS-ATRAZINE	0.000073	235	5:32	0.0	97763.	6666
2	PROMETRYNE	0.000067	318	7:28	0.0	29101.	653.8
3	CS11526	0.000138	71	1:43	0.0	7283.	228.0
4	CS11354	0.000044	109	2:36	0.0	27640.	371.3


SPEC1 TS8902203 ver 2 on UIC 2 5 20-FEB-89 DERIVED SPECTRUM 9
 SAAP1 PROMETRYNE Start 1 091051.44 10
 Coma1 IHL OF 10PP3 COFF--20, CID=0.7
 Model EI + DRU 242.2 0 -20eV LMR SYNTH GAS UP LR
 Oper1 ACQDRU619 Inlet 1
 Base1 158.0 Inten 1 103711 Masses: 50 > 250
 Norm1 158.0 RIC 1 324186 # peaks: 123
 Peak1 1000.00 m/z
 Data1 + 68 > 71 - 58




SPEC: PR0800 ver 4 on UIC 002005 1-APR-08 DERIVED SPECTRUM 9
 Samp: GS11526 Start: 17:53:19 10
 Conn: 100HC WITH 55% MECH FULL SCAN DRU V=95
 Model: EI +DRU 212.1 0 -20eV LMR SYNTH GAS UP LR
 Oper: AC03TSP

Base: 128.0 Inlet: 37636 Masses: 60 > 220
 Norm: 128.0 RIC: 105114 # peaks: 116
 Data: / 12 > 15 - 5


Daughter Spectrum of Mass 212
 (GS-11526)




SPEC: PROM81 ver 1 on UIC 002005 1-APR-88 DERIVED SPECTRUM 9
 Samp: CS11354 Start: 17150131 10
 Comp: 100HG WITH 53% MECH FULL SCAN DRU V-95
 Model: EI + DRU 200.1 eV -20eV LMR SYNTH GAS UP LR Daughter Spectrum of Mass 200
 Oper: RCO3TSP Inlet: 1 (CS-11354)
 Base: 157.9 Inlet: 1500678 Masses: 60 > 228
 Name: 157.9 RIC: 13312322 # peaks: 179
 Peak: 1000.00 m/z
 Data: + 9 > 13 - 4



SPEC: PROMETRI ver 3 on UIC 002005 1-FEB-89 DERIVED SPECTRUM 9
 SAMP: PROMETRYNE
 COM: 20UL OF 10UG/ML, .05M, FILAMENT-OFF, JET=250, VAP=230 10 1st Stage Spectrum Prometryn
 Model: EI +03MS LMR UP LR
 Oper: ACQISPLP/R
 Base: 242.2 Inlet: 2225331 Masses: 100 > 300
 Norm: 242.2 RIC: 306532 # peaks: 38
 Peak: 1000.00 muw
 Data: / 9112-418 mw = 241.14



SPEC: GS11354 ver 1 on VIC 002002 24-APR-90 DERIVED SPECTRUM 9
 Sample GS-11354 Start: 14142134 10 Dealkalated Metabolite GS-11354
 Concentration 200.0 OF 1000/mL
 Model: E1 + Q3MS LMR UP LR
 Operator CKT/FHK
 Base: 200.0 Inlet: 1 Masses: 125 > 250
 Noise: 200.0 RIC: 1 3200969 # peaks: 36
 Peak: 1000.00 m/z
 Data: 1 36



432

SPECI CS11526R ver 1 on UIC 002082 24-APR-98 DERIVED SPECTRUM 9 Hydroxy Permetrin CS-11526
 Sample CS-11526
 Concentration 200L OF 10UG/ML
 Model EI +QMS LMR UP LR
 Operator CKT/FHK
 Base 212.0 Inlet 1 4764539 Masses 125 > 250
 Normal 212.0 RIC 1 5330829 0 peaks 71
 Peak 1000.00 muu
 Data 1 44

