

Cover Sheet for

ENVIRONMENTAL CHEMISTRY METHOD

Pesticide Name: Triclopyr

MRID #: 444561-06

Matrix: Water

Analysis: Immunoassay

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at flynt.elizabeth@epa.gov.


GRM.: 95.11
EFFECTIVE: May 20, 1997
SUPERSEDES: New

Determination of Residues of Triclopyr and Trichloropyridinol in Water Using
Magnetic Particle-Based Immunoassay Test Kits


D. L. Young and C. A. Mihaliak
Global Environmental Chemistry Laboratory—Indianapolis Lab
DowElanco
Indianapolis, Indiana 46268-1054

A. Scope

This method is applicable for the quantitative determination of residues of triclopyr and its metabolite, 3,5,6-trichloro-2-pyridinol (trichloropyridinol, TCP) in water. The concentration range for analysis of triclopyr is 0.10 to 300 ng/mL, with a validated limit of quantitation of 0.10 ng/mL. The concentration range for analysis of trichloropyridinol is 0.50 to 300 ng/mL, with a validated limit of quantitation of 0.50 ng/mL.



Triclopyr
CAS Number: 55335-06-3



3,5,6-Trichloro-2-pyridinol
CAS Number: 6515-38-4

Effective Date: May 20, 1997

GRM 95.11

B. Principle

The Triclopyr and Trichloropyridinol RAPID Assay test kits apply the principles of enzyme-linked immunosorbent assay (ELISA) (1, 2) for the determination of residues in water samples. An aliquot of each sample is pipeted into a disposable test tube. Enzyme-conjugated triclopyr or trichloropyridinol and paramagnetic particles coated with specific antibodies are sequentially added to the tube. During an incubation period, the sample residue and the enzyme conjugate compete for antibody sites on the magnetic particles. At the end of the incubation period, a magnetic field is applied to the particles. The sample residue and enzyme conjugate bound to the antibodies on the particles are held in the tube by the magnetic field while the unbound reagents are decanted. After decanting, the particles are washed to remove unbound enzyme conjugate. The presence of triclopyr or trichloropyridinol is detected by adding the enzyme substrate (hydrogen peroxide) and a chromogen (3,3',5,5'-tetramethylbenzidine; TMB), generating a colored product. After another incubation period, the reaction is stopped and stabilized by the addition of acid. Since the enzyme conjugate is in competition with the sample residue for the antibody sites, the level of color development is inversely proportional to the concentration of triclopyr or trichloropyridinol in the sample (i.e., lower concentrations result in greater color development). The absorbance at 450 nm is measured in each tube using the RPA-1 RAPID Analyzer. A calibration curve is generated and the residue concentration in unknown samples is calculated from the regression equation using the preprogrammed software capabilities of the RPA-1 RAPID Analyzer.

C. Safety Precautions

1. Each analyst must be acquainted with the potential hazards of the reagents and products used in this method before commencing laboratory work. SOURCES OF INFORMATION INCLUDE: MATERIAL SAFETY DATA SHEETS, PRODUCT LITERATURE, AND OTHER RELATED DATA. Safety information on non-DowElanco products should be requested from the supplier. Disposal of reagents must be in compliance with local, state, and federal laws and regulations.
2. Avoid contact of the Stopping Solution (0.5% sulfuric acid) with skin and mucous membranes. Wear protective clothing and proper eye protection when working with this material. If this reagent comes in contact with skin, flush the exposed area with water.
3. Volatile and flammable organic solvents such as acetone and methanol must be used in well-ventilated areas away from ignition sources.

D. Equipment (Note N.1.)

1. Balance, analytical, Mettler, Model AE50, Mettler Instrument Corporation, Hightstown, NJ 08520.
2. Magnetic Separator Rack and Base, 60-position, catalog number A00004, Strategic Diagnostics Inc., Newtown, PA 18940.

Effective Date: May 20, 1997

GRM 95.11

3. Mixer, Vortex-Genie, catalog number 12-812, Fisher Scientific, Pittsburgh, PA 15238.
4. Photometer, fixed wavelength spectrophotometer RPA-1 RaPID Analyzer, catalog number A00003, Strategic Diagnostics Inc.
5. Pipetter, Eppendorf, 2-20 μ L, catalog number 21-381-201, Fisher Scientific.
6. Pipetter, Eppendorf, repeater, catalog number 21-380-8, Fisher Scientific.
7. Pipetter, Eppendorf, tri-volume, (100- μ L, 200- μ L, 250- μ L), catalog number 21-278-38, Fisher Scientific.
8. Timer, minutes and seconds with alarm, catalog number 14-649-14, Fisher Scientific.

E. Glassware and Materials (Note N.1.)

1. Bottles, 8-oz. (237-mL) graduated, with PTFE-lined caps, catalog number 03-320-11G, Fisher Scientific.
2. Culture tubes, disposable glass, 16 x 100 mm, catalog number 14-962-10B, Fisher Scientific.
3. Cylinders, mixing, 50-mL, graduated, with stoppers, catalog number 08-56C, Fisher Scientific.
4. Pipet tips, Eppendorf Combitip for repeater pipet, 12.5 mL, catalog number 21-380-8C, Fisher Scientific.
5. Pipet tips, Eppendorf disposable, 0.1 mL-1.0 mL, catalog number 21-372-4, Fisher Scientific.
6. Pipet tips, Eppendorf disposable, 1-100 μ L, catalog number 21-381-303, Fisher Scientific.
7. Pipets, 1-mL disposable, catalog number 13-678-25B, Fisher Scientific.
8. Pipets, 5-mL disposable, catalog number 13-678-25D, Fisher Scientific.
9. Pipets, 10-mL disposable, catalog number 13-678-31J, Fisher Scientific.
10. Vials, 40 mL glass, catalog number 03-339-5C, Fisher Scientific.

Effective Date: May 20, 1997

GRM 93.11

F. Reagents and Prepared Solutions (Note N.1.)

1. Reagents

- a. Acetone, Omnisolv, catalog number AX0116-1, EM Science, Gibbstown, NJ 08027.
- b. Methanol, ChromAR HPLC grade, catalog number 3041-09, Mallinckrodt Specialty Chemicals Company, Paris, KY 40361.
- c. Trichloropyridinol RaPID Assay Test Kit, catalog number A00208, Strategic Diagnostics Inc. Kit contents include (Note N.2.):
697
 - (1) Trichloropyridinol Antibody, coupled to paramagnetic particles
 - (2) Enzyme Conjugate
 - (3) Calibration Standards (0.50 ng/mL, 2.5 ng/mL and 6.0 ng/mL)
 - (4) Quality Control Sample
 - (5) Diluent/Zero Standard
 - (6) Color Solution
 - (7) Stopping Solution
 - (8) Washing Buffer
 - (9) Test Tubes
- d. Triclopyr RaPID Assay Test Kit, catalog number A00171, Strategic Diagnostics Inc. Kit contents include (Note N.2.):
5 + 1
 - (1) Triclopyr Antibody, coupled to paramagnetic particles
 - (2) Enzyme Conjugate
 - (3) Calibration Standards (0.10 ng/mL, 1.0 ng/mL and 3.0 ng/mL)
 - (4) Quality Control Sample
 - (5) Diluent/Zero Standard
 - (6) Color Solution
 - (7) Stopping Solution
 - (8) Washing Buffer
 - (9) Test Tubes
- e. Sample Diluent: trichloropyridinol, catalog number A00210; triclopyr, catalog number A00173, Strategic Diagnostics Inc. (Note N.2.)
12
12
- f. Standard:
Obtain triclopyr and TCP analytical standards from Test Substance Coordinator, DowElanco, 9330 Zionsville Rd., Building 306, Indianapolis, IN 46268-1053.
- g. Water, OmniSolv HPLC grade, catalog number WX001-4, EM Science.

Effective Date: May 20, 1997

GRM 95.11

G. Preparation of Fortification Stock Solutions (Note N.3.)

1. Preparation of Trichloropyridinol Fortification Solutions

- a. Weigh 0.050 g of the TCP analytical standard using an analytical balance. Quantitatively transfer to a 100-mL volumetric flask and dissolve in acetone. Dilute to volume with acetone to obtain a stock solution containing 500 µg/mL of TCP.
- b. Dilute 10.0 mL of the above 500-µg/mL TCP stock solution with acetone in a 100-mL volumetric flask to obtain a 50.0-µg/mL stock solution.
- c. Dilute 10.0 mL of the above 50.0-µg/mL TCP solution with acetone in a 100-mL volumetric flask to obtain a 5.0-µg/mL stock solution.
- d. Dilute 10.0 mL of the above 5.0-µg/mL TCP solution with acetone in a 100-mL volumetric flask to obtain a 0.50-µg/mL stock solution.
- e. Prepare the following fortification stock solutions in either Trichloropyridinol Sample Diluent or distilled water. Use the TCP solutions described in G.1.b., G.1.c., and G.1.d. These solutions should be prepared immediately before use.

Original Concentration µg/mL	Original Volume mL	Diluent Volume mL	Final Volume mL	Final Concentration µg/mL
50.0	6.0	14.0	20.0	15.0
50.0	1.0	19.0	20.0	2.5
5.00	1.0	19.0	20.0	0.25
0.50	1.0	19.0	20.0	0.025
0.50	0.5	19.5	20.0	0.0125

- f. Use the fortification stock solutions from G.1.e. for preparation of the fortified water samples as follows. The samples should be prepared immediately before use.

Stock Solution Concentration µg/mL	Stock Volume mL	Final Volume mL	Final Concentration µg/mL
15.0	1.0	50	0.3
2.5	1.0	50	0.05
0.25	1.0	50	0.005
0.025	1.0	50	0.0005
0.0125	1.0	50	0.00025

Effective Date: May 20, 1997

GRM 95.11

2. Preparation of Triclopyr Fortification Solutions

- a. Weigh 0.050 g of the triclopyr analytical standard using an analytical balance. Quantitatively transfer to a 100-mL volumetric flask and dissolve in methanol. Dilute to volume with methanol to obtain a stock solution containing 500 µg/mL of triclopyr.
- b. Dilute 2.0 mL of the above 500-µg/mL triclopyr solution with methanol in a 100-mL volumetric flask to obtain a 10.0-µg/mL stock solution.
- c. Dilute 10.0 mL of the above 10.0-µg/mL triclopyr solution with methanol in a 100-mL volumetric flask to obtain a 1.0-µg/mL stock solution.
- d. Dilute 10.0 mL of the above 1.0-µg/mL triclopyr solution with methanol in a 100-mL volumetric flask to obtain a 0.10 µg/mL stock solution.
- e. Prepare the following fortification stock solutions in either Triclopyr Sample Diluent or distilled water. Use the triclopyr solutions described in G.2.c. and G.2.d. These solutions should be prepared immediately before use.

Original Concentration µg/mL	Original Volume mL	Diluent Volume mL	Final Volume mL	Final Concentration µg/mL
1.0	3.0	17.0	20.0	0.150
1.0	1.0	19.0	20.0	0.050
0.10	2.5	17.5	20.0	0.0125
0.10	1.0	19.0	20.0	0.005
0.10	0.5	19.5	20.0	0.0025

Effective Date: May 20, 1997

GRM 95.11

f. Use the fortification stock solutions from G.2.b. and G.2.e. for preparation of the fortified water samples as follows. The samples should be prepared daily.

Stock Solution Concentration µg/mL	Stock Volume mL	Final Volume mL	Final Concentration µg/mL
10.0	1.5	50	0.3
10.0	0.75	50	0.15*
10.0	0.375	50	0.075
10.0	0.125	50	0.025
1.0	0.25	50	0.005
0.15	1.0	50	0.003
0.03	2.0	50	0.002
0.05	1.0	50	0.001
0.05	0.5	50	0.0005
0.0125	1.0	50	0.00025
0.005	1.0	50	0.00010
0.0025	1.0	50	0.00005

* Use 0.15 µg/mL as non-specific binding reagent (NSB).

Effective Date: May 20, 1997

GRM 93.11

H. Instrument Settings

1. To obtain results from the Trichloropyridinol RaPID Assay, use the following parameter settings on the RPA-1 RaPID Analyzer:

Parameter	RPA-1 Abbreviation	Setting
Protocol Name	Protocol Name	TCP
Data Reduction	Data Reduct	Linear Regression
Transformation	Xformation	Ln/LgtB
Number of Calibrators	# of Calibrator	4
Number of Calibrator Replicates	# of Reps:	2
Calibrator #1 Concentration	Cal #1 Conc	0.00
Calibrator #2 Concentration	Cal #2 Conc	0.50
Calibrator #3 Concentration	Cal #3 Conc	2.50
Calibrator #4 Concentration	Cal #4 Conc	6.00
Minimum Correlation	Correlation Flag	0.990
Maximum Concentration (pg/mL)	Normal Range Hi	6.00
Minimum Concentration (pg/mL)	Normal Range Low	0.25
Number of Controls	# of Controls	1
Number of Control Replicates	Ctrl Replicates	2
Number of Reagent Blanks	# Rgt Blk	0
Wavelength	Wavelength	450 nm
Read Mode	Read Mode	Absorbance
Units	Units	ng/mL
Precision of Calibrators	Rep %CV Flag	100

Effective Date: May 20, 1997

GRM 95.11

2. To obtain results from the Triclopyr RaPID Assay, use the following parameter settings on the RPA-1 RaPID Analyzer:

Parameter	RPA-1 Abbreviation	Setting
Protocol Name	Protocol Name	Triclopr
Data Reduction	Data Reduct	Linear Regression
Transformation	Xformation	Ln/LgtB
Number of Calibrators	# of Calibrator	4
Number of Calibrator Replicates	# of Reps:	2
Calibrator #1 Concentration	Cal #1 Conc	0.00
Calibrator #2 Concentration	Cal #2 Conc	0.10
Calibrator #3 Concentration	Cal #3 Conc	1.00
Calibrator #4 Concentration	Cal #4 Conc	3.00
Minimum Correlation	Correlation Flag	0.990
Maximum Concentration (pg/mL)	Normal Range Hi	3.00
Minimum Concentration (pg/mL)	Normal Range Low	0.05
Number of Controls	# of Controls	1
Number of Control Replicates	Ctrl Replicates	2
Number of Reagent Blanks	# Rgt Blk	1
Wavelength	Wavelength	450 nm
Read Mode	Read Mode	Absorbance
Units	Units	ng/mL
Precision of Calibrators	Rep %CV Flag	100

I. Determination of Recovery of Triclopyr and Trichloropyridinol in Water

1. Preparation of Recovery Samples

- Allow the samples to warm to room temperature. Measure 50.0-mL portions of the control water samples into 8-oz glass bottles using volumetric graduated cylinders. For laboratory recovery samples, add the appropriate spiking solution from Section G. An unfortified control sample, if available, should be carried through the method with each sample set.
- Vortex mix the samples well before removing the analysis aliquot.
- Assay each sample according to the procedure described in Section 1.2.
- If the sample contains more than 3.0 ng/mL of triclopyr or 6.0 ng/mL of trichloropyridinol, perform an additional dilution of the sample from Step 1.1.b. prior to assay (e.g., for a 1:10 dilution, pipet 4.5 mL of Sample Diluent or distilled water into a culture tube, add 0.5 mL of the sample from Step 1.1.b. to the tube and vortex to mix). After vortexing, wait at least 5 minutes before proceeding with the assay. Assay the diluted aliquot as described in Section 1.2.

Effective Date: May 20, 1997

GRM 95.11

2. Assay Procedure

Conduct each test in an individual test tube. The average of duplicate analyses of a sample or standard constitute a single result. A standard curve and the appropriate control and recovery samples must be included in each analytical batch. For further details, consult the Triclopyr or Trichloropyridinol RaPID Assay Kit Inserts.

Remove all kit reagents from refrigerated storage and allow them to equilibrate to room temperature prior to use. A minimum of 30 minutes is recommended for warming.

- a. Turn on the RPA-I Photoanalyzer at least 30 minutes prior to measuring absorbance in the completed assay.
- b. Label test tubes for standards, controls, and samples. Place the tubes in the proper rack position. Be sure that the rack is removed from the Magnetic Separator.
- c. Critical step: Use the indicated pipeting technique in this step to obtain accurate and precise data. Add the standard, quality control solution, non-specific binding (NSB) reagent (see G.2.f.) or sample to each test tube using an Eppendorf pipetter (250 µL with the TCP kit, 200 µL with the triclopyr kit). Pipet each sample or standard directly to the bottom of the tube; avoid liquid adhering to the sides of the test tube. Use a fresh pipet tip for each standard and sample.
- d. Using an Eppendorf repeater pipet equipped with a 12.5-mL Combitip, add 0.25 mL (Dial Setting = 1) of Enzyme Conjugate down the inside wall of each tube.
- e. Before use, thoroughly mix the Antibody Coupled Paramagnetic Particles by swirling the bottle. Avoid vigorous shaking and foaming.
- f. Using a repeater pipet equipped with a 12.5-mL Combitip, add 0.50 mL (Dial Setting = 2) of the Antibody Coupled Paramagnetic Particles down the inside wall of each tube.
- g. When dispensing of the magnetic particles has been completed, mix the samples by gently vortexing (Vortex setting = 3-4) each tube for 1-2 seconds.
- h. Incubate at room temperature (20 minutes for analysis of TCP, 30 minutes for analysis of triclopyr).
- i. After the incubation period, combine the rack and the magnetic base. Seat all tubes by pressing them into the base. Allow 2 minutes for the particles to separate.
- j. Do not separate the tube rack from the magnetic base. Using a smooth motion, invert the combined rack assembly over a collection container and pour out the tube contents. Keep the rack inverted and gently blot the test tube rims on several layers of paper towels. Do not shake or bump the rack as the magnetic particles may fall out of the tubes.

Effective Date: May 20, 1997

GRM 93.11

k. Using a repeater pipet equipped with a 12.5-mL. Combitip, add 1 mL (Dial Setting = 4) of Washing Buffer to each tube. Wait 2 minutes. Do not separate the tube rack from the magnetic base. Using a smooth motion, invert the combined rack assembly over a collection container and pour out the tube contents. Keep the rack inverted and gently blot the test tube rims on several layers of paper towels.

l. Repeat the wash step described in Step 1.2.k. an additional time.

m. Remove the tube rack from the magnetic separator and then add 0.50 mL of Color Solution to each tube using a repeater pipet equipped with a 12.5-mL. Combitip (Dial Setting = 2).

n. Gently vortex each tube for 1-2 seconds.

o. Incubate for 20 minutes at room temperature. During this incubation, pour approximately 1 mL of Washing Buffer into a clean tube for use as an instrument blank.

p. At the end of the incubation period, add 0.5 mL of Stopping Solution to each tube using a repeater pipet equipped with a 12.5-mL. Combitip (Dial Setting = 2).

q. Analyze each tube using the RPA-1 RaPID Analyzer within 15 minutes after adding the Stopping Solution.

3. Operating Procedure for the RPA-1 RaPID Analyzer

The RPA-1 RaPID Analyzer is pre-programmed with the protocols for several RaPID Assay procedures. The following steps describe how to set up and run the analyzer to measure absorbance in the tubes for the Triclopyr or Trichloropyridinol RaPID Assay.

a. Switch on the instrument and allow it to warm up at least 30 minutes prior to use. The RPA-1 RaPID Analyzer will perform a self test. If all parameters are satisfactory, the "SELECT COMMAND" prompt will appear.

b. At the "SELECT COMMAND" prompt, press "RUN".

c. At the "RUN PROTOCOL" prompt, scroll through the protocols using the arrow keys until "TCP" appears for trichloropyridinol analysis or "Triclopyr" appears for triclopyr analysis. Press "ENTER".

d. At the "SPL. REPLICATES" (sample replicates) prompt, press "2" to indicate the number of replicates for each sample, then press "ENTER".

e. At the "BLANK TUBE/INSERT TUBE" prompt, insert the tube containing approximately 1 mL of Washing Buffer. The display will briefly read "EVALUATING TUBE" then "REMOVE TUBE" and the instrument will produce

Effective Date: May 20, 1997

GRM 95.11

an audible beep indicating that the absorbance has been measured. After hearing the beep, remove the tube.

f. At the "CAL. #1 REP. #1/INSERT TUBE" prompt, insert the first replicate of the first standard/calibrator (0.00 ng/mL). Remove the tube after the beep.

g. Follow the prompts on the instrument display until all of the standards have been measured. The tube order is important. The RPA-1 RAPID Analyzer has been programmed to evaluate the standards in ascending order, in duplicate, starting with 0.00 ng/mL.

h. After all of the standards have been evaluated, the instrument will report the equation of the line, the transformed data and the standards data.

i. Insert the quality control tubes at the "CNTRL. #1 REP. #1" and "CNTRL. #1 REP. #2" prompts. The instrument will report the calculated values for each replicate of the quality control sample.

j. Evaluate the results for the standard curve and the quality control sample. At the "EDIT CALIBRATORS YES/NO" prompt, press "NO".

k. At the "SPL. #1 REP #1/INSERT TUBE" prompt, insert the first sample tube. Remove the tube after the beep.

l. Continue sample analysis following the prompts on the instrument display. Press "STOP" after all the samples have been evaluated and the results have been reported by the RPA-1 RAPID Analyzer.

I. Calculations

1. Calibration Curve

The RPA-1 RAPID Analyzer contains preprogrammed data reduction capabilities which calculate a calibration curve for each analytical batch using the absorbances of the standards supplied with the kit. The calibration curve is constructed by linear regression after performing a ln/Logit data transformation of the concentration and absorbance values, respectively.

The regression equation is :

$$\text{Logit } \frac{B}{B_0} = [\text{slope} \times \ln(\text{Conc})] + Y \text{ intercept}$$

Where:

$$\text{Logit } \frac{B}{B_0} = \ln \frac{B/B_0}{1 - B/B_0}$$

Effective Date: May 20, 1997

GRM 93.11

B = the absorbance measured at a specific concentration
B₀ = the absorbance measured for the 0.00 ng/mL standard
Conc = the concentration of the standard

An example of a calibration curve is presented in Figure 1.

2. Calculation of Trichloropyridinol or Triclopyr in Samples

The RPA-1 RaPID Analyzer will calculate the concentration of trichloropyridinol or triclopyr in each sample using the preprogrammed data reduction parameters (Note N.4.). It will report the absorbance value and calculated concentration for each sample tube as well as the mean absorbance, the mean sample concentration and the percent coefficient of variation (%CV) of the duplicate measurements for each sample. The mean values are the final result for each sample.

To calculate sample concentration, use the following equation:

$$\text{Measured Concentration} = e^a$$

Where:

$$a = \frac{\left(\text{Logit } \frac{B}{B_0} - Y \text{ intercept} \right)}{\text{slope}}$$

Example:

Water sample fortified at 0.10 ng/mL, rack positions 35 and 36 (Figure 2):

Mean absorbance (450 nm) (B) = 0.905
Mean absorbance, 0 calibrator (B₀) = 1.213

$$\begin{aligned} \frac{B}{B_0} &= \frac{0.905}{1.213} \\ &= 0.746 \end{aligned}$$

Therefore:

$$\begin{aligned} \text{Logit } \frac{B}{B_0} &= \ln \frac{0.746}{1 - 0.746} \\ &= \ln (2.937) \\ &= 1.077 \end{aligned}$$

Effective Date: May 20, 1997

$$\begin{aligned}\text{Measured Concentration} &= e^{\left[\frac{1.077 - (-0.482)}{-0.655} \right]} \\ &= 0.093 \text{ ng/mL} \\ &= 0.000093 \text{ } \mu\text{g/mL}\end{aligned}$$

3. Calculation of the Method Factor and Gross Sample Concentration

Method Factor = Dilution Factor

Thus, for a sample receiving a 1:10 additional dilution:

$$\text{Gross Sample Concentration} = 0.000093 \text{ } \mu\text{g/mL} \times 10$$

$$\text{Gross Sample Concentration} = 0.00093 \text{ } \mu\text{g/mL}$$

4. Calculation of Percent Recovery

The percent recovery for the fortified samples is calculated as follows:

$$\text{Recovery} = \frac{\text{Concentration Found}}{\text{Concentration Added}} \times 100\%$$

From J.4.:

$$\text{Recovery} = \frac{0.000093 \text{ } \mu\text{g/mL}}{0.00010 \text{ } \mu\text{g/mL}} \times 100\%$$

$$\text{Recovery} = 93\%$$

5. Correction for Recovery

For those analyses that require correction for method recovery, use the average recovery of all the recovery samples in a given analytical batch to correct for method efficiency as follows:

$$\text{Sample Concentration (corrected } \mu\text{g/mL)} = \frac{\text{Sample Concentration (gross } \mu\text{g/mL)} \times 100}{\text{Average \% Recovery}}$$

K. Determination of Triclopyr or Trichloropyridinol in Water

1. Prepare treated samples, a system (reagent) blank, fortified recovery samples, and an untreated control as described in Section I.1.

Effective Date: May 20, 1997

GRM 95.11

2. Use the RPA-1 RaPID Analyzer to calculate a standard curve and determine the gross concentration in each sample as described in Section J.

L. Quality Control

1. Analytical Batch Definition

An analytical batch of samples is defined as a group of 60 tubes. The size of the batch is based on the capacity of the magnetic separator rack. An analytical batch of less than 60 tubes can be analyzed. The first 10 tubes (rack positions 1-10) are used for duplicate analysis of the four standards and the quality control solution. Following the quality control solution, up to 25 samples (recovery or study samples) may be analyzed in duplicate (2 tubes). If more samples are to be analyzed than can be accommodated in one rack, the remaining samples should be analyzed as a different analytical batch with a new standard curve, control, and recovery samples.

2. Quality Control Solution

A quality control solution containing 3.0 ng/mL of trichloropyridinol or 2.0 ng/mL of triclopyr (supplied with the kit; Section F.I.c.(4) or F.I.d.(4), respectively, should be analyzed at the beginning of every batch of samples (rack positions 9 and 10). The quality control solution should be assayed in the same manner as all other samples. Additional fortified matrix recovery samples should be analyzed to further ensure proper execution of the method.

3. Study Samples

Study samples should be assayed in duplicate. If the concentration of the sample exceeds the range of the assay, dilute with Sample Diluent (typically a 10-fold dilution is performed) and then assay the diluted sample aliquot. Multiply the result by the appropriate method factor to obtain the final result.

4. Criteria for Acceptance of an Analytical Batch

The correlation coefficient for the linear regression of the calibration curve should be greater than 0.990. The absorbance %CV should be less than 10% for each duplicate pair of standards. The concentration %CV should be less than 20% for the quality control sample, and the recovery value should be within $\pm 20\%$ of the expected value (i.e., 2.4-3.6 ng/mL for trichloropyridinol or 1.6-2.4 ng/mL for triclopyr). If the data fail to meet these performance criteria, the analyst should evaluate the results, determine the potential source of the variation, and repeat the analysis if necessary. An example of calibration and quality control data are graphed and summarized in Figure 1.

Effective Date: May 20, 1997

GRM 95.11

5. Interferences

Thirty-seven pesticides, seventeen organic/inorganic compounds and four solvents were tested for the potential to interfere with conjugate binding in the trichloropyridinol assay (3) (Table I). Forty-five pesticides, seventeen organic/inorganic compounds and four solvents were tested for the potential to interfere with conjugate binding in the triclopyr assay (4) (Table II). None of the pesticides exhibited an I_{50} concentration below 10 ng/mL in either test kit. The I_{50} concentration is the concentration which results in a 50% inhibition of conjugate binding to the available antibodies, and it is a commonly used reference value for expressing cross reactivity and determining the extent of interference. In comparison, the I_{50} for triclopyr in the Triclopyr RaPID Assay test kit is approximately 0.78 ng/mL and 2.31 ng/mL for trichloropyridinol in the Trichloropyridinol RaPID Assay test kit.

6. Specificity/Sensitivity

The metabolites of triclopyr, 2-methoxy-3,5,6-trichloropyridine and 3,5,6-trichloro-2-pyridinol, have been tested to determine whether the Triclopyr RaPID Assay test kit will detect their presence in a water sample (3). The Triclopyr RaPID Assay test kit is sensitive to the methoxypyridine metabolite ($I_{50} = 4$ ng/mL) but not to the trichloropyridinol metabolite ($I_{50} > 10,000$ ng/mL). The Trichloropyridinol RaPID Assay test kit is not sensitive to either triclopyr or methoxypyridine metabolite at an I_{50} concentration below 10 μ g/mL.

7. Modifications and Uses

Modifications to the assay procedure are not recommended. This procedure is valid only when using reagents manufactured by SDI. This procedure is for use on water samples. Gross particulate matter should be removed prior to analysis by allowing the sample to settle, by centrifugation, or by filtration through a glass fiber filter. Validation of the method for analysis of other sample matrices would be required prior to implementing this method for sample analysis.

M. Results and Discussion

1. Method Validation

a. Recovery Levels and Precision

A method validation study was conducted to determine the recovery levels and the precision of the residue method for trichloropyridinol and triclopyr in water. The results are summarized by sample type and fortification level in Tables III-IV. Average recoveries of trichloropyridinol at all fortification levels ranged from 104% to 120% and the corresponding relative standard deviations ranged from 3% to 8%. Average recoveries of triclopyr at all fortification levels ranged from 92% to 118% and the corresponding relative standard deviations ranged from 5% to 20%.

Effective Date: May 20, 1997

GRM 95.11

b. Standard Curve Linearity

The correlation coefficient (r^2) values for the linear regression equations describing the absorbance as a function of standard concentration ranged from 0.9990 to 0.9998 for trichloropyridinol and from 0.9984 to 0.9999 for triclopyr during the method validation.

c. Limits of Detection and Quantitation

Following a published technique (5), the limit of detection (LOD) and limit of quantitation (LOQ) were calculated using the standard deviation from the results of the recovery samples fortified at the targeted LOQ, which was 0.50 ng/mL for trichloropyridinol and 0.10 ng/mL for triclopyr. The LOD was calculated as three times the standard deviation (3s), and the LOQ was calculated as ten times the standard deviation (10s). The calculated LOD was 0.141 ng/mL for trichloropyridinol and 0.06 ng/mL for triclopyr. The calculated LOQ was 0.47 ng/mL for trichloropyridinol and 0.19 ng/mL for triclopyr (Tables III-IV). The calculated LOQ values generally support the validated LOQ of 0.50 ng/mL for trichloropyridinol and 0.1 ng/mL for triclopyr. For residues that are above the LOD but below the validated LOQ, numerical values should be reported as less than the LOQ (< 0.50 ng/mL for TCP or < 0.10 ng/mL for triclopyr) because these residue levels cannot be reliably quantitated.

2. False Positive/False Negative Rate

Unfortified control samples (matrix blanks) and samples fortified at the LOD were analyzed during the study to determine the false positive and false negative rate. A false positive result occurs when residue at or above the established LOD is found in a control sample known to be free of analyte (6). A false negative occurs when no residue is detected in a sample fortified at the LOD. There were no false positives from the unfortified control samples and no false negatives reported from the LOD fortified samples analyzed during this study.

3. Ruggedness

Individual recovery values presented in Tables III-IV are generally within the range of 70-120% with the following exceptions: 14806704 well water samples fortified at 0.50 ng/mL with TCP, six samples, 126%, 122%, 132%, 130%, 132%, and 124%; 12713212 pond water samples fortified at 0.50 ng/mL with TCP, two samples, 122% and 122%; 14806701 tap water samples fortified at 0.10 ng/mL with triclopyr, three samples, 130%, 130%, and 140%; 14806706 creek surface water samples fortified at 0.10 ng/mL with triclopyr, two samples, 60% and 60%; 14806701 tap water samples fortified at 0.50 ng/mL with triclopyr, three samples, 124%, 130%, and 122%; 12713212 pond water samples fortified at 0.50 ng/mL with triclopyr, one sample, 128%; 12713212 pond water samples fortified at 1.25 ng/mL with triclopyr, one sample, 125%; 14806706 creek surface water samples fortified at 5.0 ng/mL with triclopyr, one sample, 126%; 14806706 creek surface water samples fortified at

Effective Date: May 20, 1997

GRM 95.11

75 ng/mL with triclopyr, two samples, 124% and 121%; 14806706 creek surface water samples fortified at 150 ng/mL with triclopyr, one sample, 132%. These results include data generated by multiple analysts on multiple days (7).

4. Confirmatory Method

The detection and/or quantitation of trichloropyridinol and triclopyr can be confirmed by capillary gas chromatography with mass selective detection, using method GRM 95.18, "Determination of Residues of Triclopyr, 3,5,6-Trichloro-2-pyridinol, and 2-Methoxy-3,5,6-trichloropyridine in Water by Capillary Gas Chromatography with Mass Selective Detection" (8).

5. Assay Time

The time required to analyze a typical analytical batch (25 samples or recoveries, four standards and the quality control sample, in duplicate) is from 2 to 3 hours.

N. Notes

1. Equipment, glassware, materials, reagents, and chemicals equivalent to those specified may be substituted with the understanding that their performance must be confirmed by appropriate tests. Common laboratory supplies are assumed to be readily available and are not listed here. Other immunochemical reagents should not be substituted.
2. Refrigerate all kit reagents and sample diluent at 2-8 °C. Do not freeze. Reagents may be used until the expiration date printed on the labels. The test tubes require no special storage conditions and may be stored separately from the kit reagents.
3. Standard solutions can be prepared at other concentrations by making appropriate dilutions. Store the solutions in a location that is protected from light.
4. The following information may appear as part of the raw data report:

"x.xx nd" (e.g., 0.80nd) indicates that the calculated concentration is below the minimum concentration programmed in the RPA-1 RaPID Analyzer parameters (20 pg/mL). The result should be reported "not detected" and the minimum concentration value should be noted.

"nd" indicates the absorbance measured is greater than or equal to the absorbance of the 0.00 ng/mL standard; therefore, a concentration cannot be calculated. The result should be reported as "not detected".

"x.xx HI" (e.g., 3.82HI) indicates that the calculated concentration exceeds the maximum concentration programmed in the RPA-1 RaPID Analyzer parameters (3.6 ng/mL). The sample should be diluted (e.g., 1:10) with Sample Diluent and then re-analyzed.

Effective Date: May 20, 1997

GRM 95.11

O. References

1. Tijssen, T. *Practice and Theory of Enzyme Immunoassays*, Elsevier Science: Amsterdam, 1985.
2. Meulenbergh, E. P.; Mulder, W. H.; Stoks, P. G., *Environ. Sci. Technol.*, 1995, 29, 553-561.
3. Mihaliak, C. A. "Trichloropyridinol Immunoassay Kit Development Report", 1995, unpublished report of DowElanco.
4. Mihaliak, C. A. "Triclopyr Immunoassay Kit Development Report", 1994, unpublished report of DowElanco.
5. Keith, L. H.; Crummett, W.; Deegan, J., Jr.; Libby, R. A.; Taylor, J. K.; Wentler, G., *Anal. Chem.* 1983, 55, 2210-2218.
6. Keith, L. H. *Environmental Sampling and Analysis, A Practical Guide*; Lewis: USA, 1991.
7. Mihaliak, C. A.; Berberich, S. A. In *Immunoanalysis of Agrochemicals*, ACS Symposium Series 586; Nelson, J. O.; Karu, A. E.; Wong, R. B., Eds.; American Chemical Society: Washington, DC, 1995; pp 288-300.
8. Olberding, E. L. "Determination of Residues of Triclopyr, 3,5,6-Trichloro-2-pyridinol, and 2-Methoxy-3,5,6-trichloropyridine in Water by Capillary Gas Chromatography with Mass Selective Detection", GRM 95.18, 1996, unpublished residue method of DowElanco.

The information herein is presented in good faith, but no warranty, express or implied, is given nor is freedom from any patent owned by DowElanco or by others to be inferred. In the hands of qualified personnel, the procedures are expected to yield results of sufficient accuracy for their intended purposes, but recipients are cautioned to confirm the reliability of their techniques, equipment, and standards by appropriate tests. Anyone wishing to reproduce or publish the material in whole or in part should request written permission from DowElanco.

Effective Date: May 20, 1997

GRM 93.11

Table 1. Compounds Tested for the Potential to Interfere in the Trichloropyridinol RaPID Assay Test Kit

Pesticides	Solvents (Max Tolerated)	Organic/Inorganic compounds
Alachlor	Acetone (15%)	Calcium (chloride dihydrate)
Aldicarb	Acetonitrile (5%)	Copper (chloride)
Ametryn	DMF (>20%)	Humic Acid
Atrazine	Methanol (10%)	Iron (chloride hexahydrate)
Azinphos methyl		Magnesium (chloride hexahydrate)
Benomyl		Manganese (chloride)
Carbaryl		Mercuric (chloride)
Carbendazim		Nickel (sulfate hexahydrate)
Carbofuran		Nitrate (sodium)
Chlorpyrifos		Peroxide (hydrogen)
Chlorpyrifos methyl		Phosphate (sodium, heptahydrate)
Clopyralid		Silicates (sodium meta-)
2,4-D		Sodium chloride
Diazinon		Sulfate (sodium)
Dinoseb		Sulfite (sodium)
Fenitrothion		Thiosulfate (sodium, pentahydrate)
Fluroxypyr		Zinc (chloride)
Glyphosate		
Lindane		
Malathion		
MCPA		
Methamidophos		
Metolachlor		
Methomyl		
Oxamyl		
Parathion		
Parathion methyl		
Phosmet		
Picloram		
Pirimicarb		
Pirimiphos-ethyl		
Pirimiphos-methyl		
Profenfos		
Propachlor		
Terbufos		
Thiophanate-methyl		
Triclopyr		

Effective Date: May 20, 1997

GRM 95.11

Table II. Compounds Tested for the Potential to Interfere in the Triclopyr RAPID Assay Test Kit

Pesticides	Solvents (Max Tolerated)	Organic/Inorganic compounds
Alachlor	Acetonitrile (2%)	Calcium (chloride dihydrate)
Aldicarb	Acetone (1%)	Copper (chloride)
Aldicarb Sulfoxide	DMF (0.5%)	Humic Acid
Aldicarb Sulfone	Methanol (10%)	Iron (chloride hexahydrate)
Atrazine		Magnesium (chloride hexahydrate)
Benomyl		Manganese (chloride)
Butachlor		Mercuric (chloride)
Butylate		Nickel (sulfate hexahydrate)
Captaf		Nitrate (sodium)
Carbaryl		Peroxide (hydrogen)
Carbendazim		Phosphate (sodium, heptahydrate)
Carbofuran		Silicates (sodium meta-)
Chlorpyrifos		Sodium chloride
Chlorpyrifos methyl		Sulfate (sodium)
Clopyralid		Sulfite (sodium)
2,4-D		Thiosulfate (sodium, pentahydrate)
2,4,5-T		Zinc (chloride)
Diazinon		
Dicamba		
Dichloropropene		
Dinoseb		
Disulfoton		
Fenitrothion		
Fenoxaprop		
Fluroxypyr		
Glyphosate		
Hexazinone		
Imazapyr		
Isopropyl		
Lindane		
MCPA		
Methamidophos		
Metolachlor		
Metribuzin		
Metsulfuron		
Pentachlorophenol		
Picloram		
Pirimicarb		
Pirimiphos-ethyl		
Pirimiphos-methyl		
Profenfos		
Propachlor		
Terbufos		
Thiabendazole		
Thiophanate-methyl		

Effective Date: May 20, 1997

GRM 95.11

Table III. Recovery of Trichloropyridinol from Water

Control Sample Number	Date of Analysis	TCP, ng/mL		Percent Recovery	Statistical Calculations
		Added	Found		
14806701	24-May-1995	0	ND ^a	--	--
14806704	24-May-1995	0	ND	--	--
14806706	24-May-1995	0	ND	--	--
12713212	24-May-1995	0	ND	--	--
14806701	24-May-1995	0.25	NA ^b	--	--
14806704	24-May-1995	0.25	NA	--	--
14806706	24-May-1995	0.25	NA	--	--
14806706	24-May-1995	0.25	NA	--	--
12713212	24-May-1995	0.25	NA	--	--
12713212	24-May-1995	0.25	NA	--	--
14806701	24-May-1995	0.50	0.58	116	
14806701	24-May-1995	0.50	0.59	118	
14806701	24-May-1995	0.50	0.57	114	
14806704	24-May-1995	0.50	0.63	126	
14806704	24-May-1995	0.50	0.61	122	
14806704	24-May-1995	0.50	0.66	132	
14806704	24-May-1995	0.50	0.65	130	
14806704	24-May-1995	0.50	0.66	132	
14806704	24-May-1995	0.50	0.62	124	
14806706	24-May-1995	0.50	0.57	114	
14806706	24-May-1995	0.50	0.46	92	$\bar{x} = 0.60$
14806706	24-May-1995	0.50	0.58	116	$s = 0.047$
12713212	24-May-1995	0.50	0.61	122	$(3s)^c = 0.141$
12713212	24-May-1995	0.50	0.60	120	$(10s)^d = 0.47$
12713212	24-May-1995	0.50	0.61	122	RSD = 8%

Table III. (Cont.) Recovery of Trichloropyridinol from Water

a. None detected at a detection limit of 0.25 ng/mL.
b. Not applicable (analyte was detected, but it was below the 0.50 ng/mL limit of quantitation).
c. Calculated limit of detection.
d. Calculated limit of quantitation.

Effective Date: May 20, 1997

GRM 95.11

Table IV. Recovery of Triclopyr from Water

Control Sample Number	Date of Analysis	Triclopyr, ng/mL		Percent Recovery	Statistical Calculations
		Added	Found		
14806701	26-Jan-1995	0	ND*	--	--
14806701	26-Jan-1995	0	ND	--	--
14806701	26-Jan-1995	0	ND	--	--
14806701	26-Jan-1995	0	ND	--	--
14806701	26-Jan-1995	0	ND	--	--
14806701	26-Jan-1995	0	ND	--	--
14806701	26-Jan-1995	0	ND	--	--
14806701	26-Jan-1995	0	ND	--	--
14806701	26-Jan-1995	0	ND	--	--
14806701	27-Jan-1995	0	ND	--	--
14806701	27-Jan-1995	0	ND	--	--
14806701	27-Jan-1995	0	ND	--	--
14806701	27-Jan-1995	0	ND	--	--
14806701	27-Jan-1995	0	ND	--	--
14806704	03-Feb-1995	0	ND	--	--
14806704	03-Feb-1995	0	ND	--	--
14806704	03-Feb-1995	0	ND	--	--
14806704	03-Feb-1995	0	ND	--	--
14806706	22-Mar-1995	0	ND	--	--
14806706	22-Mar-1995	0	ND	--	--
14806706	22-Mar-1995	0	ND	--	--
14806706	22-Mar-1995	0	ND	--	--
12713212	22-May-1995	0	ND	--	--
12713212	22-May-1995	0	ND	--	--
12713212	22-May-1995	0	ND	--	--
12713212	22-May-1995	0	ND	--	--
12713212	22-May-1995	0	ND	--	--

Effective Date: May 20, 1997

GRM 95.11

Table IV. (Cont.) Recovery of Triclopyr from Water

Control Sample Number	Date of Analysis	Triclopyr, ng/mL		Percent Recovery	Statistical Calculations
		Added	Found		
14806701	26-Jan-1995	0.05	NA ^b	--	--
14806701	26-Jan-1995	0.05	NA	--	--
14806701	26-Jan-1995	0.05	NA	--	--
14806701	26-Jan-1995	0.05	NA	--	--
14806701	26-Jan-1995	0.05	NA	--	--
14806704	03-Feb-1995	0.05	NA	--	--
14806704	03-Feb-1995	0.05	NA	--	--
14806704	03-Feb-1995	0.05	NA	--	--
14806704	03-Feb-1995	0.05	NA	--	--
14806704	03-Feb-1995	0.05	NA	--	--
14806706	22-Mar-1995	0.05	NA	--	--
14806706	22-Mar-1995	0.05	NA	--	--
14806706	22-Mar-1995	0.05	NA	--	--
14806706	22-Mar-1995	0.05	NA	--	--
14806706	22-Mar-1995	0.05	NA	--	--
12713212	22-May-1995	0.05	NA	--	--
12713212	22-May-1995	0.05	NA	--	--
12713212	22-May-1995	0.05	NA	--	--
12713212	22-May-1995	0.05	NA	--	--
12713212	22-May-1995	0.05	NA	--	--
14806701	26-Jan-1995	0.10	0.09	90	
14806701	26-Jan-1995	0.10	0.10	100	
14806701	26-Jan-1995	0.10	0.09	90	
14806701	26-Jan-1995	0.10	0.10	100	
14806701	26-Jan-1995	0.10	0.09	90	
14806701	26-Jan-1995	0.10	0.09	90	
14806701	26-Jan-1995	0.10	0.10	100	
14806701	26-Jan-1995	0.10	0.11	110	
14806701	26-Jan-1995	0.10	0.11	110	
14806701	26-Jan-1995	0.10	0.10	100	
14806701	26-Jan-1995	0.10	0.12	120	
14806701	26-Jan-1995	0.10	0.09	90	
14806701	26-Jan-1995	0.10	0.10	100	
14806701	26-Jan-1995	0.10	0.09	90	
14806701	26-Jan-1995	0.10	0.12	120	
14806701	26-Jan-1995	0.10	0.10	100	
14806701	26-Jan-1995	0.10	0.12	120	
14806701	26-Jan-1995	0.10	0.10	100	
14806701	26-Jan-1995	0.10	0.13	130	
14806701	26-Jan-1995	0.10	0.11	110	

Effective Date: May 20, 1997

GRM 95.11

Table IV. (Cont.) Recovery of Triclopyr from Water

Control Sample Number	Date of Analysis	Triclopyr, ng/mL		Percent Recovery	Statistical Calculations
		Added	Found		
14806701	27-Jan-1995	0.10	0.11	110	
14806701	27-Jan-1995	0.10	0.10	100	
14806701	27-Jan-1995	0.10	0.11	110	
14806701	27-Jan-1995	0.10	0.13	130	
14806701	27-Jan-1995	0.10	0.11	110	
14806701	27-Jan-1995	0.10	0.11	110	
14806701	27-Jan-1995	0.10	0.11	110	
14806701	27-Jan-1995	0.10	0.12	120	
14806701	27-Jan-1995	0.10	0.14	140	
14806701	27-Jan-1995	0.10	0.12	120	
14806704	03-Feb-1995	0.10	0.08	80	
14806704	03-Feb-1995	0.10	0.11	110	
14806704	03-Feb-1995	0.10	0.09	90	
14806704	03-Feb-1995	0.10	0.07	70	
14806704	03-Feb-1995	0.10	0.07	70	
14806704	03-Feb-1995	0.10	0.07	70	
14806704	03-Feb-1995	0.10	0.07	70	
14806704	03-Feb-1995	0.10	0.09	90	
14806704	03-Feb-1995	0.10	0.07	70	
14806704	03-Feb-1995	0.10	0.07	70	
14806706	22-Mar-1995	0.10	0.07	70	
14806706	22-Mar-1995	0.10	0.08	80	
14806706	22-Mar-1995	0.10	0.08	80	
14806706	22-Mar-1995	0.10	0.06	60	
14806706	22-Mar-1995	0.10	0.06	60	
14806706	22-Mar-1995	0.10	0.08	80	
14806706	22-Mar-1995	0.10	0.08	80	
14806706	22-Mar-1995	0.10	0.08	80	
14806706	22-Mar-1995	0.10	0.08	80	
14806706	22-Mar-1995	0.10	0.07	70	

Effective Date: May 20, 1997

Table IV. (Cont.) Recovery of Triclopyr from Water

Control Sample Number	Date of Analysis	Triclopyr, ng/mL		Percent Recovery	Statistical Calculations
		Added	Found		
12713212	22-May-1995	0.10	0.10	100	
12713212	22-May-1995	0.10	0.10	100	
12713212	22-May-1995	0.10	0.11	110	
12713212	22-May-1995	0.10	0.11	110	
12713212	22-May-1995	0.10	0.12	120	
12713212	22-May-1995	0.10	0.11	110	$\bar{x} = 0.097$
12713212	22-May-1995	0.10	0.12	120	$s = 0.019$
12713212	22-May-1995	0.10	0.11	110	$(3s)^c = 0.057$
12713212	22-May-1995	0.10	0.11	110	$(10s)^d = 0.19$
12713212	22-May-1995	0.10	0.10	100	RSD = 20%
14806701	26-Jan-1995	0.25	0.25	100	
14806701	26-Jan-1995	0.25	0.28	112	
14806701	26-Jan-1995	0.25	0.26	104	
14806704	03-Feb-1995	0.25	0.20	80	
14806704	03-Feb-1995	0.25	0.19	76	
14806704	03-Feb-1995	0.25	0.23	92	
12713212	22-May-1995	0.25	0.25	100	$\bar{x} = 0.246$
12713212	22-May-1995	0.25	0.27	108	$s = 0.031$
12713212	22-May-1995	0.25	0.28	112	RSD = 13%
14806701	26-Jan-1995	0.50	0.60	120	
14806701	26-Jan-1995	0.50	0.62	124	
14806701	26-Jan-1995	0.50	0.55	110	
14806701	26-Jan-1995	0.50	0.65	130	
14806701	26-Jan-1995	0.50	0.57	114	
14806701	26-Jan-1995	0.50	0.61	122	
14806701	27-Jan-1995	0.50	0.57	114	
14806701	27-Jan-1995	0.50	0.57	114	
14806701	27-Jan-1995	0.50	0.56	112	
14806704	03-Feb-1995	0.50	0.47	94	
14806704	03-Feb-1995	0.50	0.48	96	
14806704	03-Feb-1995	0.50	0.48	96	
14806706	22-Mar-1995	0.50	0.58	116	
14806706	22-Mar-1995	0.50	0.59	118	
14806706	22-Mar-1995	0.50	0.54	108	
12713212	22-May-1995	0.50	0.59	118	$\bar{x} = 0.57$
12713212	22-May-1995	0.50	0.64	128	$s = 0.05$
12713212	22-May-1995	0.50	0.60	120	RSD = 9%

Effective Date: May 20, 1997

Table IV. (Cont.) Recovery of Triclopyr from Water.

Control Sample Number	Date of Analysis	Triclopyr, ng/mL		Percent Recovery	Statistical Calculations
		Added	Found		
14806701	26-Jan-1995	1.0	1.08	108	
14806701	26-Jan-1995	1.0	1.11	111	
14806701	26-Jan-1995	1.0	1.14	114	
14806701	26-Jan-1995	1.0	1.13	113	
14806701	26-Jan-1995	1.0	1.08	108	
14806701	26-Jan-1995	1.0	1.09	109	
14806701	27-Jan-1995	1.0	1.10	110	
14806701	27-Jan-1995	1.0	1.09	109	
14806701	27-Jan-1995	1.0	1.08	108	
14806704	03-Feb-1995	1.0	0.95	95	
14806704	03-Feb-1995	1.0	0.95	95	
14806704	03-Feb-1995	1.0	1.07	107	
14806706	22-Mar-1995	1.0	1.01	101	
14806706	22-Mar-1995	1.0	1.12	112	
14806706	22-Mar-1995	1.0	1.07	107	
12713212	22-May-1995	1.0	1.25	125	$\bar{x} = 1.08$
12713212	22-May-1995	1.0	1.03	103	$s = 0.07$
12713212	22-May-1995	1.0	1.10	110	RSD = 6%
14806701	26-Jan-1995	2.0	2.02	101	
14806701	26-Jan-1995	2.0	1.75	88	
14806701	26-Jan-1995	2.0	1.99	100	
14806701	26-Jan-1995	2.0	1.92	96	
14806701	26-Jan-1995	2.0	1.95	98	
14806701	26-Jan-1995	2.0	1.98	99	
14806701	27-Jan-1995	2.0	1.90	95	
14806701	27-Jan-1995	2.0	1.94	97	
14806701	27-Jan-1995	2.0	2.00	100	
14806704	03-Feb-1995	2.0	1.82	91	
14806704	03-Feb-1995	2.0	1.74	87	
14806704	03-Feb-1995	2.0	1.92	96	
14806706	22-Mar-1995	2.0	1.90	95	
14806706	22-Mar-1995	2.0	2.03	102	
14806706	22-Mar-1995	2.0	1.83	92	
12713212	22-May-1995	2.0	2.16	108	$\bar{x} = 1.95$
12713212	22-May-1995	2.0	1.89	95	$s = 0.13$
12713212	22-May-1995	2.0	2.27	114	RSD = 7%

Effective Date: May 20, 1997


Table IV. (Cont.) Recovery of Triclopyr from Water

Control Sample Number	Date of Analysis	Triclopyr, ng/mL		Percent Recovery	Statistical Calculations
		Added	Found		
14806701	26-Jan-1995	3.0	2.82	94	
14806701	26-Jan-1995	3.0	2.75	92	
14806704	03-Feb-1995	3.0	2.48	83	
14806704	03-Feb-1995	3.0	2.56	85	
14806704	03-Feb-1995	3.0	2.60	87	
14806706	22-Mar-1995	3.0	2.77	92	
14806706	22-Mar-1995	3.0	2.60	87	
14806706	22-Mar-1995	3.0	2.67	89	
12713212	22-May-1995	3.0	2.76	92	$\bar{x} = 2.75$
12713212	22-May-1995	3.0	3.12	104	$s = 0.19$
12713212	22-May-1995	3.0	3.09	103	RSD = 7%
14806706	22-Mar-1995	5.0	5.5	110	
14806706	22-Mar-1995	5.0	5.8	116	
14806706	22-Mar-1995	5.0	5.9	118	$\bar{x} = 5.9$
14806706	22-Mar-1995	5.0	6.3	126	$s = 0.3$
14806706	22-Mar-1995	5.0	5.9	118	RSD = 5%
14806706	22-Mar-1995	25	25	100	
14806706	22-Mar-1995	25	27	108	
14806706	22-Mar-1995	25	24	96	$\bar{x} = 25$
14806706	22-Mar-1995	25	22	88	$s = 1.9$
14806706	22-Mar-1995	25	27	108	RSD = 8%
14806706	22-Mar-1995	75	93	124	
14806706	22-Mar-1995	75	91	121	
14806706	22-Mar-1995	75	86	115	$\bar{x} = 88$
14806706	22-Mar-1995	75	83	111	$s = 4.0$
14806706	22-Mar-1995	75	87	116	RSD = 5%

Control Sample Number	Date of Analysis	Triclopyr, ng/mL		Percent Recovery	Statistical Calculations
		Added	Found		
14806706	22-Mar-1995	150	169	113	
14806706	22-Mar-1995	150	156	104	
14806706	22-Mar-1995	150	162	108	$\bar{x} = 171$
14806706	22-Mar-1995	150	198	132	$s = 14$
14806706	22-Mar-1995	150	168	112	RSD = 8%
14806706	22-Mar-1995	300	276	92	
14806706	22-Mar-1995	300	293	98	
14806706	22-Mar-1995	300	290	97	$\bar{x} = 300$
14806706	22-Mar-1995	300	286	95	$s = 28$
14806706	22-Mar-1995	300	356	119	RSD = 9%
				$\bar{x} = 102$	
				$s = 15$	
				$n = 159$	

- Page 30 of 34

Effective Date: May 20, 1997



Concentration ng/mL	ln (Conc)	Absorbance 450 nm	B/B ₀	Logit (B/B ₀)
0	NA*	1.213	NA	NA
0.10	-2.30	0.896	0.7387	1.041
1.0	0.00	0.450	0.3710	-0.526
3.0	1.10	0.287	0.2366	-1.171

* NA = Not Applicable.

Slope = -0.655
Y-intercept = -0.482
Correlation coefficient (r^2) = 0.9994
Quality Control Solution = 2.22 ng/mL

Figure 1. A Typical Calibration Curve for Triclopyr

Effective Date: May 20, 1997

DowElanco Protocol Number RES94069
Analysis DEBRA YOUNG
Analysis Date Time
Start 2/3/95 3:15 PM
End 2/3/95 3:39 PM
Matrix WATER
Analysis Method TRICLOPYR
Method TRICLOPYR

Batch Information:				Instrument Settings:	
Batch ID #	020795A.DLY			Instrument protocol	TRICLOPR
Data Filename	020795A.DLY			Data Reduction	Lin Regression
				Transformation	Ln/Lg18
				Read Mode	Absorbance
				Wavelength	450 nm
				Units	NG/ML
				Quantitative range (ng/mL)	0.1-3.0

NSB Blank Absorbance = 0.147

Calibrator Data:						
Concentration	Absorbance	Difference	Observed concentration	% difference	Absorbance	% CV
0.00	1.215					
0.00	1.211					
Mean	0.00	1.213				0.2
	0.10	0.928	-0.071	0.08	-26.30	
	0.10	0.865	0.019	0.12	16.30	
Mean	0.10	0.896	-0.002	0.10	-2.20	5.0
	1.00	0.451	0.066	1.07	6.20	
	1.00	0.430	0.075	1.08	7.00	
Mean	1.00	0.450	0.071	1.07	6.60	0.3
	3.00	0.279	0.040	3.04	1.30	
	3.00	0.296	-0.297	2.70	-11.00	
Mean	3.00	0.287	-0.155	2.86	-4.70	4.1

Transformed data:		
Concentration	Transformed concentration	Transformed absorbance
0.10	-2.30	1.041
1.00	0.00	-0.536
3.00	1.10	-1.171

Equation of line:	
Slope	-0.655
Intercept	-0.482
Correlation (r)	0.9994

Quality control samples:						
Run position	Absorbance	Concentration	Mean Absorbance	Mean Concentration	% Recovery	% CV
9	0.315	2.37				
10	0.337	2.06	0.326	2.22	111.0	9.9
Comments:						
14806704 WATER VALIDATION						

Figure 2. Representative Data for the Determination of Triclopyr in Water

Effective Date: May 20, 1997

DowElanco Protocol#			RES94069		Matrix	WATER				
Analyst			DEBRA YOUNG		Analyte	TRICLOPYR				
Batch ID #			020395ADLY							
Data Filename			020395A.DLY							
Analysis Date			2/3/95							
Sample ID	Rack position	Rep #	Absorbance	Conc. (ng/mL)	Mean absorbance	Mean conc. (ng/mL)	Conc. %CV	Corr F&C	Final F&C (ng/mL)	
0A	11	1	1.239	nd						
	12	2	1.244	nd						
	mean									
0B	13	1	1.270	nd				1.0		
	14	2	1.241	nd						
	mean									
0C	15	1	1.274	nd				1.0		
	16	2	1.293	nd						
	mean									
0D	17	1	1.272	nd				1.0		
	18	2	1.296	nd						
	mean									
0E	19	1	1.292	nd				1.0		
	20	2	1.289	nd						
	mean									
0.05A	21	1	1.021	0.04nd	1.025	0.04nd	4.6	1.0	0.04nd	
	22	2	1.028	0.03nd						
	mean									
0.05B	23	1	1.047	0.03nd	1.057	0.03nd	14.7	1.0	0.03nd	
	24	2	1.066	0.02nd						
	mean									
0.05C	25	1	1.080	0.02nd	1.090	0.02nd	18.9	1.0	0.02nd	
	26	2	1.100	0.01nd						
	mean									
0.05D	27	1	1.045	0.03nd	1.061	0.02nd	26.5	1.0	0.02nd	
	28	2	1.078	0.02nd						
	mean									
0.05E	29	1	1.091	0.02nd	1.106	0.01nd	31.3	1.0	0.01nd	
	30	2	1.120	0.01nd						
	mean									
0.1A	31	1	0.946	0.070	0.934	0.08	12.7	1.0	0.08	
	32	2	0.921	0.060						
	mean									
0.1B	33	1	0.892	0.100	0.877	0.11	13.1	1.0	0.11	
	34	2	0.862	0.120						
	mean									
0.1C	35	1	0.924	0.080	0.905	0.09	17.3	1.0	0.09	
	36	2	0.887	0.100						
	mean									

Figure 2. (Cont.) Representative Data for the Determination of Triclopyr in Water

Effective Date: May 20, 1997

DowElanco Protocol#			RES94069		Matrix		WATER		
Analyst			DEBRA YOUNG		Analyte		TRICLOPYR		
Batch ID #			020395ADLY						
Data Filename			020395A.DLY						
Analysis Date			2/3/95						
Sample ID	Rack position	Rep #	Absorbance	Conc. (ng/mL)	Mean absorbance	Mean conc. (ng/mL)	Conc. RCV	Corr Factor	Final Conc. (ng/mL)
0.1D	37	1	0.947	0.070	0.950	0.07	7.6	1.0	0.07
	38	2	0.952	0.070					
	mean								
0.1E	39	1	0.955	0.070	0.945	0.07	10.0	1.0	0.07
	40	2	0.933	0.070					
	mean								
0.1F	41	1	0.944	0.070	0.949	0.07	5.3	1.0	0.07
	42	2	0.954	0.070					
	mean								
0.1G	43	1	0.948	0.070	0.943	0.07	4.8	1.0	0.07
	44	2	0.939	0.070					
	mean								
0.1H	45	1	0.885	0.110	0.906	0.09	19.3	1.0	0.09
	46	2	0.926	0.080					
	mean								
0.1I	47	1	0.934	0.080	0.941	0.07	4.9	1.0	0.07
	48	2	0.947	0.070					
	mean								
0.1J	49	1	0.936	0.070	0.946	0.07	10.5	1.0	0.07
	50	2	0.957	0.060					
	mean								
0	51	1	Empty	Empty	Empty	Empty	Empty	1.0	Empty
	52	2	Empty	Empty					
	mean								
0	53	1	Empty	Empty	Empty	Empty	Empty	1.0	Empty
	54	2	Empty	Empty					
	mean								
0	55	1	Empty	Empty	Empty	Empty	Empty	1.0	Empty
	56	2	Empty	Empty					
	mean								
0	57	1	Empty	Empty	Empty	Empty	Empty	1.0	Empty
	58	2	Empty	Empty					
	mean								
0	59	1	Empty	Empty	Empty	Empty	Empty	1.0	Empty
	60	2	Empty	Empty					
	mean								

Figure 2. (Cont.) Representative Data for the Determination of Triclopyr in Water