

Cover Sheet for

ENVIRONMENTAL CHEMISTRY METHOD

Pesticide Name: Tepraloxym

MRID #: 444673-03

Matrix: Soil

Analysis: HPLC/UV

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at flynt.elizabeth@epa.gov.

BASF CORPORATION
AGRICULTURAL PRODUCTS GROUP
AGRICULTURAL PRODUCTS CENTER
P.O. Box 13528
Research Triangle Park, NC 27709-3528

Study Title:

VALIDATION OF BASF METHOD No. D9517: ANALYTICAL METHOD FOR THE DETERMINATION
OF RESIDUES OF BAS 620 H METABOLITES (DP-1 AND DP-2) IN SOIL USING HPLC/UV
DETECTOR

Method No. D9517
Study No. 96020

Data Requirement:

Guideline 164-1 Terrestrial Field Dissipation

Study Completion Date:

November 22, 1996

Authors:

Manasi Saha
Leonard Collins

Performing Laboratory:

BASF CORPORATION
AGRICULTURAL PRODUCTS GROUP
P.O. Box 13528
Research Triangle Park, NC 27709

BASF Registration Document No.96/5254

This report consists of 75 pages.

This page intentionally left blank.


GOOD LABORATORY PRACTICES COMPLIANCE STATEMENT

This study meets the requirements for 40 CFR 160, Good Laboratory Practices, with the exceptions given in Appendix A, and the following: not all data entered in the study notebook was signed and dated on the day of entry.

Study Director: _____

Manasi Saha

Sponsor: _____



Submitter: _____

BASF CORPORATION
AGRICULTURAL PRODUCTS GROUP
AGRICULTURAL PRODUCTS CENTER
P.O. Box 13528
Research Triangle Park, NC 27709-3528

**ANALYTICAL METHOD FOR THE DETERMINATION OF RESIDUES OF BAS 620 H METABOLITES
(DP-1 AND DP-2) IN SOIL USING HPLC/UV DETECTOR**

AUTHORS: Manasi G. Saha and Leonard Collins
STUDY DIRECTOR: Manasi G. Saha (919) 547-2232
SUPERVISORY PERSONNEL: Arno Krotzky (919) 547-2213
WORK DONE BY: Leonard Collins and Manasi Saha
BASF Corporation Agricultural Product Center
P. O. Box 13528
Research Triangle Park, NC 27709

Method No. D9517

Report Date: November, 1996

ABSTRACT:

Analytical Method No. D9517 was developed to determine trace residues of 3-hydroxy-2-(1-iminopropyl)-5-(tetrahydropyran-4-yl)cyclohex-2-en-1-one (DP-1) and 2-ethyl-6-(tetrahydropyran-4-yl)-4, 5, 6, 7-tetrahydrobenzoxazol-4-one (DP-2) in soil. The purpose of the study is to determine recovery efficiency of the above analytes in soil. BAS 620 H is the active ingredient and DP-1 and DP-2 are two major metabolites found in several environmental fate studies (Reference 1). Method development was conducted at BASF Corporation. Soil samples were extracted with methanol-water (75:25, v/v) using a homogenizer. The extracts were concentrated to dryness under vacuum. The residues were chromatographed by RP C18 SPE and determined by HPLC using a UV detector. This study shows that Analytical Method No. D9517 is suitable for measuring residues of DP-1 and DP-2 down to levels of 0.01 ppm in soil. Over a fortification span from 0.01 to 1.0 ppm, recoveries ranged from 77 to 119% for DP-1 with a mean recovery of $92 \pm 10\%$ (N = 29). Over the same span, recoveries for DP-2 ranged from 72 to 155% with a mean recovery of $93 \pm 19\%$ (N = 29).

PAGES OF REPORT: 75

Study Initiation Date:	January 30, 1996
Experimental Initiation Date:	February 5, 1996
Experimental Completion Date:	May 2, 1996

BASF CORPORATION AGRICULTURAL PRODUCTS GROUP

Agricultural Research Center, Research Triangle Park, N.C. 27709

STATEMENT OF THE QUALITY ASSURANCE UNIT

Method Number: D9517

BASF Study Number: 96020

Name/Number of Test Substance:

3-hydroxy-2-(1-iminopropyl)-5-(tetrahydropyran-4-yl)cyclohex-2-en-1-one (BH 620-DP1),
Lot number 00345-268, 97.9 %.

2-ethyl-6-(tetrahydropyran-4-yl)-4,5,6,7-tetrahydrobenzoxazol-4-one (BH 620-DP2),
Lot number 00665-27, 98.3% purity.

Type of Study: Method Validation

Study Initiation Date: January 30, 1996

The quality assurance unit of the testing facility at the APC has audited this study, the raw data, and the final report and reported its findings to the study director and to management.


Date of Inspection	Report to Study Director and to Management
February 6, 1996	February 6, 1996
February 13, 1996	February 13, 1996
June 19, 1996	June 19, 1996
 Signature of QAU	

TABLE OF CONTENTS

ABSTRACT	4
1. INTRODUCTION AND SUMMARY	9
1.1 Scope and Source of the Method	
1.1.1 Scope	9
1.1.2 Source	9
1.2 Description of Test and Reference Substances	9
1.3 Principle of the Method	10
2. MATERIALS AND METHODS	10
2.1 Equipment	10
2.2 Reagents and Chemicals	12
2.3 Standard Substances and Solutions	12
2.3.1 Standard Solutions for Fortifications	13
2.3.2 Standard Solutions for HPLC-UV Analysis	14
3. ANALYTICAL PROCEDURE	15
3.1 Parent and Metabolite Isolation and Cleanup	15
3.1.1 Sample Preparation	15
3.1.2 Fortification of Procedural Recovery Samples	15
3.1.3 Extraction	15
3.1.4 Sample Clean-up: Silica Gel Solid Phase Extraction	16
3.1.4.1 Column Preparation	16
3.1.4.2. Column Conditioning	16
3.1.4.3. Loading, washing and elution	17
3.1.5 Preparation for Sample Analysis	17
3.1.6 Moisture Determination	17
3.2 Instrumentation	18
3.2.1 Description of Equipment	18
3.2.2 Typical Operating Conditions	18
3.2.3 Calibration Procedures	18
3.2.4 Sample Analysis	19
3.3 Interferences	19
3.3.1 Sample Matrices	19
3.3.2 Other Sources	19
3.4 Confirmatory Techniques	19
3.5 Time Required for Analysis	19
3.6 Potential Problems	20
4. METHODS OF CALCULATION	20
4.1 Calibration	20
4.2 Analyte in Sample	20
4.3 Calculation of Procedural Recoveries	21
5. VALIDATION	22
5.1 Description of Protocol	22
5.2 Solution Stability	23
5.3 Protocol Changes	23

TABLE OF CONTENTS, Continued

6.	RESULTS AND DISCUSSION	23
6.1	General	23
6.2	Accuracy and Precision	24
6.3	Determination Limit	24
6.4	Ruggedness Testing	24
6.5	Limitations	24
7.	CONCLUSIONS	24
8.	QUALITY ASSURANCE PROCEDURES	25
9.	REFERENCES	25
10.	SAFETY AND HEALTH CONSIDERATIONS	25
10.1	General	25
10.2	Solvents, Reagents and Standards	25
11.	SPECIAL NOTES	26
12.	CERTIFICATION	27

FIGURES

Figure 1.	Flow Chart for Analytical Method D9517	28
Figure 2.	Typical Calculation for the HPLC-UV Quantitation of DP-1 Residues in Soil.	29
Figure 3.	Typical Calculation for the HPLC-UV Quantitation of DP-2 Residues in Soil.	30

TABLES

TABLE I.	Recoveries of DP-1 from Soil	31
TABLE II.	Recoveries of DP-2 from soil.	34
TABLE III.	Soil Characterization Data	37
TABLE IV A.	DP-1 Standards used to calculate samples 1044-06-10 to 28.	38
TABLE IV B.	DP-2 Standards used to calculate samples 1044-06-10 to 28.	38
TABLE V A.	DP-1 Standards used to calculate samples 1044-14-09 to 27.	39
TABLE V B.	DP-2 Standards used to calculate samples 1044-14-09 to 27.	39
TABLE VI A.	DP-1 Standards used to calculate samples 1044-22-09 to 27.	40

TABLE OF CONTENTS, Continued

TABLE VI B.	DP-2 Standards used to calculate samples 1044-22-09 to 27.	40
TABLE VII A.	DP-1 Standards used to calculate samples 1044-24-09 to 27.	41
TABLE VII B.	DP-2 Standards used to calculate samples 1044-24-09 to 27.	41
TABLE VIII A.	DP-1 Standards used to calculate samples 1044-34-09 to 27.	42
TABLE VIII B.	DP-2 Standards used to calculate samples 1044-34-09 to 27.	42

APPENDICES

APPENDIX A

Changes to Protocol Number 96020	43-44
--	-------

APPENDIX B

Typical Raw Data for Analyses	45-75
-------------------------------------	-------

1. INTRODUCTION AND SUMMARY

1.1 Scope and Source of the Method

1.1.1 Scope

This method is used to determine the trace residues of 3-hydroxy-2-(1-iminopropyl)-5-(tetrahydropyran-4-yl)cyclohex-2-en-1-one (DP-1) and 2-ethyl-6-(tetrahydropyran-4-yl)-4, 5, 6, 7-tetrahydrobenzoxazol-4-one (DP-2) in soil. Caloxydim is the active ingredient and DP-1 and DP-2 are two major metabolites found in several environmental fate studies (Reference 1).

1.1.2 Source

This method was developed at BASF Corporation, RTP, NC.


1.2 Description of Test and Reference Substances

Fortification and HPLC standards

BASF Code: DP-1

Chemical name: 3-hydroxy-2-(1-iminopropyl)-5-(tetrahydropyran-4-yl)cyclohex-2-en-1-one

Structural formula:



Empirical formula: $C_{14}H_{21}NO_3$

Molecular weight: 251.33


Appearance: White solid

1.2 Description of Test and Reference Substances, continued

BASF Code: DP-2

Chemical name: 2-ethyl-6-(tetrahydropyran-4-yl)-4, 5, 6, 7-tetrahydrobenzoxazol-4-one

Structural formula:

Empirical formula: $C_{14}H_{19}NO_3$

Molecular weight: 249.31

Appearance: White solid

1.3 Principle of the Method

Residues of DP-1 and DP-2 are extracted from soil with methanol-water (75:25, v/v) using a homogenizer. The extracts are then filtered through celite and concentrated to dryness. The dry residues are further cleaned up by C18 SPE. Final determination is made by HPLC using UV detector. See Figure 1 for a flow chart of the analytical method. The limit of quantitation is 0.01 ppm for both analytes.

2. MATERIALS AND METHODS

2.1 Equipment

Suggested Sizes/Manufacturer

Centrifuge

IEC Refrigerated Centrifuge
Model PR 7000
Beckmann Refrigerated
Centrifuge
Model CS-6KR

Centrifuge Bottles, Teflon-lined screw cap

Fisher Scientific Co.
150 mL

Empty glass SPE-column (8 mL)
(with Teflon frits)

Baker Chemical Co.
Cat No. 7308-06

Filtering flask

Fisher Scientific Co.
1000 mL

MATERIALS AND METHODS, continued

Homogenizer and accessories	Omni International Omni 5000 or Omni GLH 115, Generator (G20-195 ST, 20 mm i.d X 198 mm l)
Inert Sampling Adapter, B & J	Burdick & Jackson Cat No. 9473A
Laboratory Shaker	Janke and Kunkel Model HS501-D
Rotary evaporator	Büchi Rotovapor Model RE 111, 114
Standard taper Flat-bottom flasks	VWR Scientific Co. 500 mL, 300 and 125 mL
Solid Phase Extraction Manifold	Baxter Healthcare Corporation
Standard funnels	Fisher Scientific Co.
Fritted (medium porous, about 40-60 mm) funnels, 60 mL	Fisher Scientific Co.
Ultrasonic Bath	Fisher Scientific Co., Model FS-14
Vacuum distilling Vacubrand vacuum pump/controller	Elnik Systems, Inc.
Vacuum distilling adapter	Aldrich Chemical Co. ACS p Cat. No. Z17,067-4
Other general laboratory glassware and supplies as needed.	

NOTE: The equipment listed in this section was used in the development of this method. Equipment with equivalent performance may be used as required.

MATERIALS AND METHODS, continued

2.2 Reagents and Chemicals

Source/Preparation

Acetonitrile, CAS 75-05-8

Baxter Healthcare Corporation,
B & J Brand

Celite 545®

J. T Baker
Cat. No. 3371-01

Formic acid (88%), CAS 64-18-6

Aldrich Chemical Co. ACS
Cat. No. 39, 938-8Gelman Nylon acrodisc
(Membrane disc 0.45 µm)

VWR/ Catalog No 28143-948)

Methanol, CAS 67-56-1

Baxter Healthcare Corporation,
B & J BrandReversed phase silica gel
(C18 bulk pack, 40 µm)J.T. Baker
Cat No. 7025-00

Water, CAS 7732-18-5

Baxter Healthcare Corporation,
B & J Brand

NOTE: Equivalent reagents and chemicals from other suppliers may be substituted.

2.3 Standard Substances and Solutions

The standard compounds shown in the table below were used for method development and validation.

Compound	Code	Lot Number	Purity, %
3-hydroxy-2-(1-iminopropyl)-5-(tetrahydropyran-4-yl)cyclohex-2-en-1-one	BH 620-DP1	00345-268	97.9 %
2-ethyl-6-(tetrahydropyran-4-yl)-4,5,6,7-tetrahydrobenzoxazol-4-one	BH 620-DP2	00665-27	98.3 %

Standard supplied by:

Dr. Rita Laschober
BASF Aktiengesellschaft, APS/UP
Agricultural Research Center
D-67114 Limburgerhof, West Germany
Telephone: 06236/68/2103

MATERIALS AND METHODS, continued

Solid 3-hydroxy-2-(1-iminopropyl)-5-(tetrahydropyran-4-yl)cyclohex-2-en-1-one (DP-1), 2-ethyl-6-(tetrahydropyran-4-yl)-4, 5, 6, 7-tetrahydrobenzoxazol-4-one (DP-2) were maintained frozen ($<-5^{\circ}\text{C}$) until their use in this study. These substances were characterized as required by 40 CFR part 160, FIFRA Good Laboratory Practices. Data on the synthesis and subsequent characterization of these substances are available to BASF and are located at Landwirtschaftliche Versuchsstation der BASF, Limburgerhof, Germany..

Solutions of DP-1 and DP-2 were refrigerated ($+4^{\circ}\text{C}$) during their use in this study. Stock solutions (1 mg/mL) of DP-1 and DP-2 were made fresh every three months and dilutions of the stock solution were made monthly.

Solution Stability Note: Methanol solutions of DP1 and DP-2 are kept in the dark at 4°C were stable for at least 103 days. A statement concerning the stabilities of solutions used in the method validation study for this report is given in Section 5.2.

2.3.1 Standard Solutions for Fortifications

NOTE: These standard concentrations are suggested. A different concentration scheme may be used and additional standards may be prepared as needed.

2.3.1.1 3-hydroxy-2-(1-iminopropyl)-5-(tetrahydropyran-4-yl)cyclohex-2-en-1-one (DP-1): 1.0 mg/mL in Methanol

Prepare a 1.0 mg/mL stock solution by weighing an appropriate amount of DP-1 into a volumetric flask and dissolving it with an appropriate amount of methanol. For example, to prepare a 10 mL stock solution, dissolve 10.0 mg of DP-1 in 10 mL of methanol in 10 mL volumetric flask.

2.3.1.2 2-ethyl-6-(tetrahydropyran-4-yl)-4,5,6,7-tetrahydrobenzoxazol-4-one (DP2): 1.0 mg/mL in methanol

Prepare a 1.0 mg/mL stock solution by weighing an appropriate amount of DP-2 into a volumetric flask and dissolving it with an appropriate amount of methanol. For example, to prepare a 10 mL stock solution, dissolve 10.0 mg of DP-2 in 10 mL of methanol in a 10 mL volumetric flask.

2.3 Standard Substances and Solutions, Continued

2.3.1.3 Mix Standard Solutions of 3-hydroxy-2-(1-iminopropyl)-5-(tetrahydropyran-4-yl)cyclohex-2-en-1-one (DP-1) and 2-ethyl-6-(tetrahydropyran-4-yl)-4,5,6,7-tetrahydrobenzoxazol-4-one (DP2): 50 µg/mL in methanol.

Prepare a 50 µg/mL of mixed standard solution by transferring an appropriate amount of each of the 1.0 mg/mL stock solutions (2.3.1.1 and 2.3.1.2) with a volumetric pipette into a volumetric flask (typically 5.0 mL of each of the 1.0 mg/mL stock solutions in a 100 mL volumetric flask). Dilute to the mark with methanol.

Prepare a 5.0 and 0.5 µg/mL of mixed standard solution by making appropriate dilutions from the 1.0 µg/mL of fortification standard solution with methanol.

Transfer each stock and standard solution to an amber bottle fitted with a Teflon-lined screw cap and store at refrigerator. Replace stock solution 90 days after preparation. Replace standard solutions 30 days after preparation.

2.3.1.4 Injection Standard Solutions of DP-1 and DP-2 for HPLC-UV Analysis in 1:1 Acetonitrile-Water: 1.0, 0.5, 0.25 and 0.125 µg/mL

Prepare a 1.0 µg/mL of injection standard solutions by transferring an appropriate amount of the 50.0 µg/mL mixed standard solution with a volumetric pipette into a volumetric flask. Typically add 1.0 mL of the 50.0 µg/mL standard solution into a 50 mL volumetric flask. Add 1.0 mL of water and then dilute to the mark with 1:1 acetonitrile-water.

Prepare a 0.5, 0.25 and 0.125 µg/mL of injection standard solutions by making dilutions from the 1.0 µg/mL of injection standard solution with 1:1 acetonitrile-water.

For HPLC analysis, inject 50.0 µL of the mixed DP-1 and DP-2 standard solutions of concentrations 0.125, 0.25, 0.5 and 1.0 µg/mL to construct a standard curve. This is a suggested calibration scheme and may be altered as needed.

3. ANALYTICAL PROCEDURE (See Figure 1, Flow chart for Analytical Method D9517)

3.1 Parent and Metabolite Isolation and Cleanup

3.1.1 Sample Preparation

Bulk soil samples received from the field are homogenized using a blender or mill. Homogenized soil samples are stored frozen (-5°C) before analysis. Weigh a 50 g or to the nearest tenth of a gram aliquot of the soil sample into a 150 mL centrifuge bottle.

3.1.2 Fortification of Procedural Recovery Sample

It is recommended to analyze at least two procedural recovery samples and one untreated sample (control) with each analysis set to monitor method efficiency. Typically, run one procedural recovery sample at the limit of quantitation (0.01 ppm) along with one procedural recovery sample at the expected residue level.

For each fortification, pipette an appropriate amount of standard DP-1 and DP-2 fortification solutions prepared in 2.3.1.3 to control soil samples. For example, 1.0 mL of the 5.0 $\mu\text{g/mL}$ standard added to 50 g soil results in a fortification level of 0.1 ppm.

3.1.3 Extraction

Add 80 mL methanol-water (75:25) to the soil samples and extract with homogenizer for 5 minutes (about 6000 RPM). Stop the homogenizer and raise the blade slightly (just enough that the blade still remain in the solution) and turn on the homogenizer for a minute to rinse the main soil clogging in the blade. Raise the blade completely out of the solution and carefully rinse the homogenizer blade with about 20 to 30 mL methanol-water (75:25).

NOTE: Soil clogs in the generators. It is important to rinse this thoroughly and have to be done as soon as polytron is stopped. Rinse in and out of the generators as much as possible with methanol water (75:25).

Cap the centrifuge bottle and vortex for 2 to 3 minutes in a vortex mixer and then shake at 300 rpm for 10 minutes. Cap and centrifuge at 3000 rpm for 10 min. (0°C). Transfer the supernatant into a 500 mL flat bottom flask by decantation.

Add 50 mL methanol-water (75:25) to the soil marc and repeat above extraction. Transfer the supernatant into the above 500 mL flat bottom flask by decantation.

3. ANALYTICAL PROCEDURE, Continued

Add another aliquot of 50 mL methanol-water (75:25) to the soil marc and repeat above extraction. Transfer the supernatant into the above 500 mL flat bottom flask by decantation.

Note: In repeat extractions marc settles very tightly at the bottom of the centrifuge bottle due to centrifugation. Lower down the polytron generator all the way until it hits the marc (analyst will feel a resistance and perhaps it is touching the bottom of the centrifuge bottle, but generator will just penetrate the soil marc) and continue the extraction. Soil Extract after third extraction looks cloudy (specially for clay soil)

Attach a vacuum filtration adapter to a 500 mL flat bottom flask and connect a fritted filtering funnel (60 mL size) to the top of the vacuum filtering adapter. Add enough celite to cover one third of the fritted filtering funnel. Suction filter the combined extracts through the bed of celite. Rinse the flask and celite bed thoroughly with methanol (4 X 20 to 30 mL). Release the suction and rinse the stem of the fritted filtering funnel and vacuum filtering adapter with methanol.

NOTE: Filtering flask with side arm can also be used for filtration. It is recommended not to use Buchner funnel with filter paper for the filtration due to escape of the silt through the celite.

Concentrate the combined extract to about 20 to 30 mL using a rotary evaporator with the water bath temperature set at approximately 60°C (set vacuum initially at about 200 mbar until the removal of methanol and then gradually decrease to about 35 to 40 mbar). Transfer the contents in to a 125 mL round bottom flask (use about 10 to 20 mL methanol to rinse) and evaporate nearly to dryness with a bath temperature maintained at 60°C (vacuum at about 35 to 45 mbar).

NOTE: The extracts needs to evaporate nearly to dryness with few drops (100 ul) to maximum of 1 ml left in the flask.

3.1.4 Sample Clean-up: RP C18 Solid Phase Extraction

3.1.4.1 Column Preparation

Add 1.0 g of C18 silica gel to an empty glass column fitted with a Teflon frit. Pack glass wool (enough to leave half of the column space empty) tightly on the top of the C18 silica gel and add

3. ANALYTICAL PROCEDURE, Continued.

celite on the top (enough to leave 1.0 cm space of the column empty). Add a sampling adapter on the top of the column to connect an additional empty glass reservoir.

3.1.4.2 Column conditioning

NOTE: Vacuum is used for chromatography using a solid phase extraction manifold.

Condition C18 silica gel by passing through approximately 10 mL acetonitrile followed by approximately 25 mL water. Do not allow the column to go dry. The applied vacuum should be sufficient to permit flow rates no more than 2-5 mL/min. Vacuum readings of about 10 kPa (2-5 inches of mercury) have been adequate.

3.1.4.3 Loading, Washing and Elution

Add 1 mL acetonitrile to the sample from section 3.1.3, swirl, sonicate and vortex thoroughly to dissolve the residue from the side of the flask as much as possible before addition of 50 mL of water. Add 50 mL water, sonicate thoroughly [Duration of sonication depends on soil type] to ensure the dissolution of the residues from the side of the flask and quickly transfer (load) to the top of the conditioned C18 silica gel column. Allow the solvent to pass through the column using vacuum. Collect all eluant in a suitable container (250 mL beaker) and discard. Do not allow the column to go dry.

Wash the flask with 50 mL water and add this wash to the column.

Wash the flask with 25 mL 2% acetonitrile-water, v/v and add this wash to the column to further wash the column. Discard all the acetonitrile-water wash.

Add 50 mL of 60:40 acetonitrile-water, v/v to the flask and add to the top of the column. Collect the solvent in a 50 mL brown Teflon-lined screw cap bottle and transfer to a 125 mL flat bottom flask.

NOTE: We found during method development, elution solvent with less percentage of acetonitrile e.g 40:60;v/v acetonitrile-water will also work. But 60:40 v/v acetonitrile-water was required for clay type soil to get high recovery of DP-2.

Evaporate the solution to dryness using a rotary evaporator with the water bath temperature set at approximately 60°C (set vacuum

3. ANALYTICAL PROCEDURE, Continued

initially set at about 200 mbar to remove acetonitrile and then gradually decrease to about 35 to 40 mbar). Remove the sample immediately after evaporation.

3.1.5 Preparation for Sample Analysis

For HPLC-UV determination, dissolve each sample with an appropriate amount of mobile phase (30:70 v/v, acetonitrile-water containing 0.1 % formic acid) just before the analysis. Typical volumes used are 2.0 mL for control and 0.01 ppm fortifications; 10 mL for 0.1 fortifications; and 50 mL for 1 ppm fortifications. For example for samples at the limit of quantitation (0.01 ppm), the final volume is 2.0 mL. Add mobile phase (30:70 v/v, acetonitrile-water containing 0.1 % formic acid, 2.0 mL), swirl, sonicate and vortex to ensure complete dissolution of sample from the side of the flask. Filter the solution through a membrane syringe filter (a 0.45 micron membrane disc fitted to a 1.0 mL disposable plastic syringe and samples were transferred with an aid of glass disposable pipette to the plastic disposable syringe) into the injection vial.

3.1.6 Moisture Determination

Soil analysis results are reported on a "dry weight" basis. Therefore soil sample weights must be corrected for moisture content by any method the laboratory customarily uses. See section 12, Note 2 for an example of a moisture determination procedure.

3.2 Instrumentation

3.2.1 Description of Equipment

Liquid Chromatograph: Varian 9010 pump, Varian 9050 Detector, Varian 9100 Autosampler and VG Data System Multichrom 2 data processing System.

HPLC Column: Nucleosil 5, C18, 250 mm X 3.2 mm

NOTE: The equipment listed was used for method development. Equivalent equipment may be used.

3.2.2 Typical Operating Conditions

Injection Volume: 50 µl

3. ANALYTICAL PROCEDURE, Continued

Isocratic Mobile Phase: 30 % Acetonitrile + 70 % water containing
0.1 % formic acid
Flow Rate: 0.5 ml/min.
Wavelength (λ_{max}): 254 nm.
DP-1 Retention Time: 6.5 min.
DP-2 Retention Time: 9.4 min.

NOTE: The preceding specifications are suggested and may be altered as needed. Actual use conditions and any changes must be documented.

3.2.3 Calibration Procedures

Inject two or more mixed standards of DP-1 and DP-2 until stable responses are observed. Calculation of results is based on peak area measurements using a calibration curve. The calibration curve is obtained by injecting various amounts of the standard solution (e.g. 0.125, 0.25, 0.5 and 1.0 ng/ μ L) concurrently with sample analysis. Different standard concentrations may be used as appropriate. Injection of control in the beginning of the analysis is necessary. Plot detector response (peak area) versus weight of standard injected. The standards should give a linear response.

3.2.4 Sample Analysis

Inject 50 μ L of sample into the HPLC-UV.

Directly compare the response (peak area) of unknown samples injected with the calibration curve to obtain ng of DP-1 and DP-2 injected. Bracket every 2-3 samples with standards to check for shifts in sensitivity or retention time. To do this, an injection sequence including standards and samples must be planned.

If the peak response of the unknown is larger than the highest standard, dilute the unknown appropriately and re-inject.

3.3 Interferences

3.3.1 Sample Matrices

If interfering peaks from the matrix occur in the chromatogram, change the HPLC operating conditions (see 3.2.2) or use an alternative HPLC column. It is desirable to clean the chromatographic system periodically by injecting with the solvent.

3.3.2 Other Sources

3. ANALYTICAL PROCEDURE, Continued

Other Pesticides: None known to date.
Solvents: None known to date.
Labware: None known to date.

3.4 Confirmatory Techniques

HPLC coupled with diode-array UV detection confirms the presence of DP-1 and DP-2 (Figures B.23. to B.27.). GC-MSD (EI Full MS Scan) confirms the structural identity of the analytes (Figures B.28. to B.29.). No problems with interferences or questionable peak identity have been encountered to date.

3.5 Time Required for Analysis

Analysis of a set of 7 soil samples requires 1.5 working days, including sample work-up and HPLC-UV analysis.

3.6 Potential Problems

Potential technical problems have been described at the appropriate points in the method.

4. METHODS OF CALCULATION (See Figure 2 for an example calculation.)

4.1 Calibration

Construct a linear least squares calibration curve in the form $y = bx + c$ from the standards by plotting peak area versus weight of standard injected for both analytes.

4.2 Analyte in Sample

Calculation of results is based on peak area measurements. Using the peak area measurements for DP-1 and DP-2 in the samples, the amount of the analyte in ng from the appropriate least squares calibration curve is determined. See Figure 3.

Calculate ppm values by the equation below.

$$\text{ppm} = \frac{A}{B}$$

where: A = ng value interpolated from calibration curve
(Obtained automatically or could be calculated using Fig. 3)

4. METHODS OF CALCULATION, continued

$$B = \text{mg Sample Injected} = \frac{\text{Dry Sample Wt. (g)} \times \mu\text{L Injected}}{\text{Final dilution volume (mL)}}$$

The "final dilution volume" includes any dilutions which have been made.

Using the Dry Sample Weight in the "mg Sample Injected" calculation will yield a ppm value on a "dry weight" basis. The Dry Sample Weight is obtained after determining the moisture content of the original sample (see Section 3.1.9).

4.3 Calculation of Procedural Recoveries

Correct fortification results for residues found in the control sample as follows:

$$\text{ppm (corrected)} = \text{ppm in fortified control} - \text{ppm in control}$$

Determine percent recovery from the fortification experiments as follows:

$$\% \text{ Recovery} = \frac{\text{ppm (corrected)} \times 100}{\text{ppm DP-1 or DP-2 added}}$$

Only results for procedural recovery samples should be corrected for residues in the controls. Do not correct treated sample results for either control residues or recoveries.

5. VALIDATION:

5.1 Description of Protocol

The validation was carried out as BASF Study 96020. Control soil samples were fortified with DP-1 and DP-2 at levels ranging from 0.01 to 1.0 ppm. The fortified controls were analyzed and the results converted to recovery values for evaluation. Percent recoveries for DP-1 are compiled in Table I. Percent recoveries for DP-2 are compiled in Table II.

The Test Substances used to generate validation data are given below:

Compound	Code	Lot Number	Purity, %
3-hydroxy-2-(1-iminopropyl)-5-(tetrahydropyran-4-yl)cyclohex-2-en-1-one	BH 620-DP1	00345-268	97.9 %
2-ethyl-6-(tetrahydropyran-4-yl)-4, 5, 6, 7-tetrahydrobenzoxazol-4-one	BH 620-DP2	00665-27	98.3 %

These materials were supplied by:

BASF Corporation
Agricultural Research Center
P. O. Box 13528
Research Triangle Park, NC 27709

The Test System consisted of untreated soil samples obtained from caloxydim trial sites (BASF Study 95023). Two different soil types (clay and sandy-loam) were used to validate this method. Both types of soil were obtained from California site and identified as BASF RCN 95008. Soil characterization data for these soil samples are summarized in Table III.

The HPLC column used was manufactured by Phenomenex and had serial number 80249.

5.2 Solution Stability

During the course of this study, the stabilities of fortification and HPLC standard solutions were examined. Solutions were stored in a refrigerator at 4°C. Fortification standards (stock solution) made in methanol showed no degradation for DP-1 and DP-2 after 103

5. VALIDATION, continued

days of storage. Solutions kept in the dark at 4°C were stable for at least 103 days. No degradation of HPLC injection standards made in acetonitrile-water (1:1) occurred for periods up to 21 days. The lifetime of working standard fortification solutions was set at one month and HPLC injection standards was set at three weeks. For further information see section 2.3.

5.3 Protocol Changes

No Changes were made to the validation study protocol. Some typographical error found in the interim report (attachment to the study protocol) were corrected in the Final report (Appendix A).

6. RESULTS AND DISCUSSION

6.1 General

The present work describes an analytical method to measure two major metabolites of BAS 620 H, DP-1 and DP-2 residues in clay and sandy loam soil.

Recoveries of DP-1 ranged from 77 to 119% with a mean recovery of $92 \pm 10\%$ ($N = 29$) and are summarized in Table I. Recoveries of DP-2 ranged from 72 to 155% with a mean recovery of $93 \pm 19\%$ ($N = 29$) and are summarized in Table II. This study has shown that Analytical Method No. D9517 is suitable for measuring residues of DP-1 and DP-2 in soil down to 0.01 ppm.

Two types of soil were used as substrates for the validation study, clay and sandy loam. These soil samples represent the types of soils on which BAS 620 H is customarily used. Soil characterization data are given in Table III.

Standard curves used to calculate the recoveries shown in Tables I and II were generated from standard solutions containing DP-1 and DP-2 injected concurrently with the analysis set. Standard injections bracketed the sample injections. Summaries of all the injection standard data used in this study are given in Tables IV through VIII.

Representative HPLC-UV chromatograms of DP-1 and DP-2 residue analyses of selected control and fortified samples are shown in Figures B.7 through B.22. Typical standard chromatograms are shown in Figures B.1 through B.4. A typical standard curve for DP-1 is given in Figure B.5. A typical standard curve for DP-2 is shown in Figure B.6. A chromatogram and diode-array UV spectra of DP-1 and

6. RESULTS AND DISCUSSION, Continued

DP-2 standards and Sample 96020/1044-22-24 and 96020/1044-34-22 are shown in Figures B.23 through B. 27. Full Scan EI MS spectra of DP-1 and DP-2 standards are presented in Figures B.28 and B.29.

6.2 Accuracy and Precision

Recoveries of DP-1 ranged from 77 to 119% with a mean recovery of $92 \pm 10\%$ ($N = 29$). Recoveries of DP-2 ranged from 72 to 155% with a mean recovery of $93 \pm 19\%$ ($N = 29$). No other statistical analysis was performed on the data.

6.3 Determination Limit

The determination limit for DP-1 and DP-2 residues in soil was 0.01 mg/kg.

6.4 Ruggedness Testing

One analyst executed five sets of extractions and analyses in four different two-day periods. Mean recoveries and estimates of standard deviations for each analysis date were calculated from the data in Tables I and II:

Date of Analysis	DP-1 Recovery	DP-2 Recovery
2/7/96	90 ± 7	91 ± 10
2/9/96	95 ± 14	113 ± 30
2/13/96	95 ± 12	97 ± 16
2/15/96	95 ± 11	83 ± 13
2/22/96	87 ± 4	92 ± 11

¹Six fortified samples and one control sample were analyzed.

The means of the analyses for each set were in the range 70-120%. Chromatograms consisted only of two peaks: one peak at the retention time of each analyte.

6.5 Limitations

None known to date.

7. CONCLUSIONS

This study has shown that Analytical Method No. D9517 is suitable for measuring residues of DP-1 and DP-2 in soil down to 0.01 ppm.

8. QUALITY ASSURANCE PROCEDURES

The raw data and analytical standards of this study will be stored in the BASF archives at:

BASF Corporation
Agricultural Product Center
26 Davis Drive
Research Triangle Park, NC 27709

9. REFERENCES

1. a) Keller, E., BASF Report No.3633 "Aerobic Soil Metabolism of ¹⁴C-BAS 191819 (BAS 620 H)", October, 1994
- b) Keller, E., BASF Report No.3729 "Degradation Behavior of Reg. No. of ¹⁴C-BAS 191819 (BAS 620 H) in Sterile Soil", December, 1993
- c) Yamasaki, R., Nippon Soda Report No. EC-554 "BAS 620 H-anaerobic Aquatic Metabolism Study", December, 1993.
- d) Followings are on-going study and planned to be completed in 1996:
 - i) Kanji Ishihara, Nippon Soda Report No. EC-737 "BAS 620 H Photodegradation Study in water".
 - ii) Shiotani H., NISSO Report No. EC-440 "BAS 620 H- Hydrolysis Study ".
 - iii) Shiotani H., NISSO Report No. EC-518 "BAS 620 H- Photodegradation in Soil Study.

10. SAFETY AND HEALTH CONSIDERATIONS

10.1 General

Use personal protective equipment such as lab coats, safety glasses and gloves (nitrile/latex gloves are recommended) when performing the operations described in this method. Conduct all transfers, partitions, derivatizations, nitrogen-stream evaporations and SPE procedures in a well-ventilated hood. Guard vacuum equipment such as rotovaps to minimize the possibility of injury caused by flying broken glass. Dispose of hazardous wastes in an environmentally acceptable manner, in compliance with applicable laws and regulations.

10.2 Solvents, Reagents and Standards

It is recommended to review the Material Safety Data Sheets (MSDSs) for all solvents and reagents used in this method. The toxicity of the DP-1 and DP-2 are unknown.

11. SPECIAL NOTES

1. Section 3.1.6: An example procedure for moisture determination is as follows:

Weigh 5 g of wet soil ("Wet weight") accurately into a tarred glass petri dish or other container.

Place into a 150°C oven for 16 hours (overnight).

Remove the petri dish from the oven and allow to cool in a desiccator.

Working quickly, remove the cool petri dish from the desiccator and weigh accurately to obtain "Dry Weight".

Determine the ratio $R = \text{"Dry Weight"} / \text{"Wet Weight"}$

Calculate "Dry Sample Weight" = $R \times \text{"Weight Sample Weight"}$ (See section 3.1.1).

The calculated "Dry Sample Weight" is used in section 4.2 to calculate concentration (ppm) values.


2. Sections 3.1.3 to 3.1.5: In the development of this method, BASF personnel had the use of a vacuum pump and controller on rotary evaporators. At various points in the method it was expedient to vary the pressure in the rotary evaporators. The approximate settings used are shown in the table below.

Solvent	Temperature, C	Vacuum pressure, mbar	Method section
Acetonitrile	60	180	3.1.4
Methanol	60	200	3.1.4
Water	60	35	3.1.4

12. CERTIFICATION

We, the undersigned, hereby declare that this report provides a true and accurate record of the results obtained.

Author: Manasi Saha Date: 11/22/96
Manasi G. Saha

Approved By:  Date: 11/18/96
Technical Center Leader
Residue Environmental Fate-1/MS Lab




Figure 1. Flow Chart for Analytical Method D9517

BASF Sample Number 96020/1044-24-12. Fortified Soil at 0.01 ppm for DP-1 Recovery.

A = ng value interpolated from standard curve

$$\text{Standard curve: ng (DP-1)} = \frac{\text{Peak Area} - 3.41566e+002}{2.63151e+003}$$

Peak Area: 29253

Use full computer/calculator precision in any intermediate calculations. Round only the final value.

$$\begin{aligned} A = \text{ng (DP-1)} &= \frac{29253 - 3.416e+002}{2.632e+003} \\ &= 10.987 \end{aligned}$$

$$B = \text{mg sample injected} = \frac{\text{Sample weight (g)} \times \mu\text{L injected}}{\text{Final dilution volume (mL)}}$$

$$= \frac{50 \times 50.0}{2.0} = 1250 \text{ mg}$$

D = molecular weight conversion factor = 1

$$\begin{aligned} \text{ppm} &= \frac{A \times D}{B} = \frac{10.987 \times 1}{1250} \\ &= 0.008789 \text{ ppm} \end{aligned}$$

The corresponding control sample (BASF Sample 96020/1044-24-09) contained 0.000 ppm of DP-1 residue.

$$\text{Net Recovery} = 0.008789 - 0.0000 = 0.008789 \text{ ppm (88\%)}$$

Figure 2. Typical Calculation for the HPLC-UV Quantitation of DP-1 Residues in Soil.

BASF Sample Number 96020/1044-24-12. Fortified Soil at 0.01 ppm for DP-2 Recovery.

A = ng value interpolated from standard curve

$$\text{Standard curve: ng (DP-2)} = \frac{\text{Peak Area} - 2.823856e+002}{1.51176e+003}$$

Peak Area: 22206

Use full computer/calculator precision in any intermediate calculations. Round only the final value.

$$\begin{aligned} A = \text{ng (DP-1)} &= \frac{22206 - 2.823856e+002}{1.51176e+003} \\ &= 14.383 \text{ ng} \end{aligned}$$

B = mg sample injected = $\frac{\text{Sample weight (g)} \times \mu\text{L injected}}{\text{Final dilution volume (mL)}}$

$$= \frac{50 \times 50}{2} = 1250 \text{ mg}$$

D = molecular weight conversion factor = 1

$$\text{ppm} = \frac{A \times D}{B} = \frac{14.383 \times 1.0}{1250}$$

$$= 0.011506 \text{ ppm}$$

The corresponding control sample (BASF Sample 96020/1044-24-09) contained 0.003923 ppm of DP-2 residue.

$$\text{Net Recovery} = 0.011506 - 0.003923 = 0.007583 \text{ ppm (76\%)}$$

Figure 3. Typical Calculation for the HPLC-UV Quantitation of DP-2 Residues in Soil.