

Cover Sheet for

ENVIRONMENTAL CHEMISTRY METHOD

Pesticide Name: Dichloropropene Degradate

MRID #: 445365-04

Matrix: Soil

Analysis: GC/MS

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at flynt.elizabeth@epa.gov.

Appendix B

DowElanco Method GRM 94.17


GRM: 94.17
EFFECTIVE: July 26, 1995
SUPERSEDES: New

Determination of Residues of *cis*- and *trans*-3-Chloroacrylic Acid
in Soil by
Capillary Gas Chromatography with Mass Selective Detection


D. O. Duebelbeis and S. E. Fisher
North American Environmental Chemistry Laboratory
DowElanco
Indianapolis, Indiana 46268-1053

A. Scope

This method is applicable for the quantitation of residues of the 1,3-dichloropropene metabolites, *cis*- and *trans*-3-chloroacrylic acid (CAAC) in soil. The method was validated over the concentration range 0.20 ng/g to 2.0 µg/g with a limit of quantitation of 0.20 ng/g.



cis-3-Chloroacrylic Acid
CAS 1609-93-4



trans-3-Chloroacrylic Acid
CAS 2345-61-1

B. Principle

Residues of CAAC in soil are extracted with acidified acetone. The extract is partially evaporated to remove acetone, diluted with deionized (DI) water, and adjusted to a neutral pH. CAAC residues are concentrated using an ion-exchange solid-phase extraction column (SPE). The CAAC is eluted from the SPE in 0.1 N hydrochloric acid. The eluent is further acidified, saturated with sodium chloride and CAAC residues are partitioned into methyl-*t*-butyl ether (MTBE). The MTBE is passed through a silica gel SPE column to remove water and particulate. Isooctane is added and the MTBE is evaporated. CAAC residues in isooctane are derivatized with N-methyl-N-(*t*-butyldimethylsilyl)trifluoroacetamide (MTBSTFA) to their corresponding *t*-butyldi-methylsilyl esters (CAAC TBDMSE) and analyzed by capillary gas chromatography with mass selective detection (GC/MSD). Soils containing levels of CAAC above 20 ng/g are diluted 100-fold with isooctane, rederivatized and reanalyzed.

Effective Date: July 26, 1995

C. Safety Precautions

1. Each analyst must be acquainted with the potential hazards of the reagents, products, and solvents used in this method before commencing laboratory work. SOURCES OF INFORMATION INCLUDE: MATERIAL SAFETY DATA SHEETS, LITERATURE, AND OTHER RELATED DATA. Safety information on non-DowElanco products should be obtained from the container label or from the supplier. Disposal of reagents, reactants, and solvents must be in compliance with local, state, and federal laws and regulations.
2. Acetic acid, acetone, isooctane, methanol and MTBE are flammable and should be used in well-ventilated areas away from ignition sources.
3. Concentrated acetic acid is corrosive and can cause severe burns. It is imperative that proper eye and personal protection equipment be used when handling this reagent.
4. *cis*- and *trans*-3-Chloroacrylic Acid are corrosive and lachrymators. It is imperative that proper eye and personal protection equipment be used when handling these reagents.
5. MTBSTFA is irritating to eyes, respiratory system and skin. It is imperative that proper eye and personal protection equipment be used when handling this reagent.

D. Equipment (Note N.1.)

1. Automatic sampler, Model 7673, Hewlett-Packard, Wilmington, DE 19808.
2. Balance, analytical, Model AE200, Mettler Instrument Corporation, Hightstown, NJ 08520.
3. Balance, pan, Model BB2440, Mettler Instrument Corporation.
4. Centrifuge, with rotor to accommodate 12-mL vials, Model Centra-8, International Equipment Company, Needham Heights, MA 02194.
5. Evaporator, N-Evap, Model 111, Organomation Associates, Inc., South Berlin, MA 01549.
6. Gas chromatograph, Model 5890 Series II, Hewlett-Packard.
7. Mass selective detector, Model 5971A, Hewlett-Packard, Palo Alto, CA 94304.
8. Mass selective detector data system, Model G1034B, Hewlett-Packard.
9. Shaker, variable speed reciprocating with box carrier, Model 6000, Eberbach Corporation, Ann Arbor, MI 48106.
10. Ultrasonic bath, Model 1200, Branson Cleaning Equipment Company, Shelton, CT 06484.
11. Vacuum manifold box, Model spe-21, J.T. Baker, Inc., Phillipsburg, NJ 08865.
12. Vial Crimper, catalog number 8710-0979, Hewlett-Packard, Wilmington, DE 19808.

Effective Date: July 26, 1995

13. Vortex mixer, Model G-560, Scientific Industries, Inc., Bohemia, NY 11716.
14. Water purification system, Model Milli-Q UV Plus, Millipore Corporation, Milford, MA 01757.

E. Glassware and Materials (Note N.1.)

1. Column, capillary gas chromatography, Durabond-5 liquid phase, 30 m x 0.25 mm i.d., 0.25 μ m film thickness, catalog number 122-5032, J&W Scientific, Folsom, CA 95630.
2. Column, silica gel SPE, catalog number 7086-07, J.T. Baker, Inc.
3. Column, strong anion-exchange (quaternary amine) SPE, catalog number 71225-6013, Varian Sample Preparation Products, Harbor City, CA 90710.
4. Column adapter, PTFE, catalog number 120-1100, Jones Chromatography, Inc., Lakewood, CO 80228.
5. Column inlet liner, deactivated, catalog number 5181-3315, Hewlett-Packard.
6. Column reservoir, 75 mL, catalog number 7120-03, J.T. Baker, Inc.
7. Filter, charcoal, catalog number 7972, Chrompack, Inc., Raritan, NJ 08869. (Note N.2.)
8. Filter, moisture, catalog number 7971, Chrompack, Inc. (Note N.2.)
9. Filter, oxygen, catalog number 7970, Chrompack, Inc. (Note N.2.)
10. Gas, helium, 99.995% purity, Airco, Murray Hill, NJ 07974.
11. Gas, nitrogen, 99.99% purity, Airco.
12. Indicator strips, pH range 0.0 to 6.0, 0.5 pH gradation, product number 4391-01, J.T. Baker, Inc.
13. Indicator strips, pH range 4.5 to 10.0, 0.5 pH gradation, product number 4395-01, J.T. Baker, Inc.
14. Microdispenser, 25 μ L, Drummond Dialamatic Microdispenser, catalog number 3000225, Drummond Scientific Company, Broomall, PA 19008.
15. Microdispenser replacement bore, 25 μ L, catalog number 3000225G, Drummond Scientific Company.
16. Syringes, 100, 250, and 500 μ L capacity, catalog numbers 80600, 80700, and 80800, Hamilton Co., Reno, NV 89520.
17. Tube, Pyrex brand culture tube with threaded end, catalog number 14-957-86D, Fisher Scientific, Pittsburgh, PA 15219.
18. Vial, autosampler, 2 mL, catalog number C4011-1, National Scientific Co., Lawrenceville, GA 30243.

Effective Date: July 26, 1995

19. Vial, 11-mL with PTFE-lined screw cap, catalog number 2504T, Qorpak, Pittsburgh, PA 15205.
20. Vial, 45-mL with PTFE-lined screw cap, catalog number 60958A-11, Kimble Glass, Vineland, NJ 08360.
21. Vial seal, catalog number C4011-1A, National Scientific Company.

F. Reagents and Chemicals (Note N.1.)

1. Reagents

- a. Acetic acid, HPLC grade, catalog number A35-500, Fisher Scientific, Pittsburgh, PA 15219.
- b. Acetone, Optima grade, catalog number A929-4, Fisher Scientific.
- c. Hydrochloric acid, 0.1 N, ACS reagent grade, certified concentration, catalog number SA54-4, Fisher Scientific.
- d. Hydrochloric acid, 2.0 N, ACS reagent grade, certified concentration, catalog number SA431-500, Fisher Scientific.
- e. Isooctane, Optima grade, catalog number O301-4, Fisher Scientific.
- f. Methanol, Optima grade, catalog number A454-4, Fisher Scientific.
- g. MTBE, methyl-*t*-butyl ether, HPLC grade, catalog number E127-4, Fisher Scientific.
- h. MTBSTFA, N-methyl-N-(*t*-butyldimethylsilyl)trifluoroacetamide, catalog number 48920, Pierce, Rockford, IL 61105.
- i. Sodium chloride, ACS reagent grade, catalog number S271-1, Fisher Scientific.
- j. Sodium hydroxide, 0.1 N, ACS reagent grade, certified concentration, catalog number SS276-4, Fisher Scientific.
- k. Sodium sulfate (anhydrous), certified ACS grade, catalog number S421-500, Fisher Scientific.

l. Standards

(1) *cis*-3-Chloroacrylic acid

The *cis*-CAAC standard, TSN100370, used for this study was originally obtained from Aldrich Chemical Co., Milwaukee, WI 53233, catalog number 17,740-7, lot number 04926EW, with a purity of >98% (1).

(2) *trans*-3-Chloroacrylic acid

The *trans*-CAAC standard, TSN100371, used for this study was originally obtained from Aldrich Chemical Co., catalog number C2,235-0, lot number 09114TW, with a purity of >99% (2).

Obtain from Test Substance Coordinator, DowElanco, Indianapolis, IN 46268-1053.

2. Prepared Solutions

- a. 90% acetone/10% 0.1 N hydrochloric acid solution

Pipet 200 mL of 0.1 N hydrochloric acid into a 2000-mL volumetric flask containing approximately 1000-mL of acetone. Swirl the flask and allow to equilibrate to room temperature. Dilute to volume with acetone.

- b. 0.025% Acetic acid in MTBE

Deliver 250 μ L of acetic acid into a 1000-mL volumetric flask containing approximately 500-mL of MTBE. Swirl the flask and dilute to volume with MTBE.

G. Preparation of Standards

All solutions prepared in Section G should be stored in amber bottles and sealed with PTFE-lined caps.

1. Preparation of *cis*- and *trans*-CAAC Stock Solutions

- a. Weigh 0.1000 g of *cis*-CAAC analytical standard and quantitatively transfer to a 100-mL volumetric flask. Dilute to volume with acetone to obtain a 1000 μ g/mL stock solution.
- b. Weigh 0.1000 g of *trans*-CAAC analytical standard and quantitatively transfer to a 100-mL volumetric flask. Dilute to volume with acetone to obtain a 1000 μ g/mL stock solution.

2. Preparation of *cis*- and *trans*-CAAC Spiking Solutions

- a. Transfer 1.0 mL of each of the stock solutions in Sections G.1.a. and b. to a 100-mL volumetric flask and bring to volume with acetone to obtain an initial solution of 10.0 μ g/mL for each *cis*- and *trans*-CAAC. This solution is used for spiking and preparation of diluted spiking and calibration solutions.
- b. Solutions for spiking soil samples are prepared by adding approximately 10 mL of acetone and 10 μ L of acetic acid to a 100-mL volumetric flask. The flask is agitated to allow acetic acid to contact the glass surface. The appropriate aliquot of the initial solution from Section G.2.a. is then added and diluted to volume with acetone as follows:

Effective Date: July 26, 1995

Aliquot of 10.0 µg/mL Soln. mL	Final Soln. Volume mL	Spiking Soln. Final Conc. ng/mL
0.100	100	10.0
0.200	100	20.0
0.500	100	50.0
1.00	100	100.
2.00	100	200.
10.0	100	1000.

- c. A 2.00-ng/mL spiking solution is prepared by adding approximately 10 mL of acetone and 10 µL of acetic acid to a 100-mL volumetric flask. The flask is agitated to allow acetic acid to contact the glass surface. A 1.0 mL aliquot of the 200 ng/mL spiking solution from Section G.2.b. above is added and diluted to volume with acetone. A 1.0 mL aliquot of this solution will fortify a 10-g soil sample at the limit of quantitation, 0.20 ng/g.
- d. A summarization of the spiking solutions prepared, aliquots to be delivered to a 10-g soil sample and the resulting sample concentration is shown below.

Spiking Soln. Concentration	Control Soil Mass g	Aliquot to Deliver mL	Equivalent Sample Conc.* ng/g
2.0 ng/mL	10.0	1.0	0.20
10.0 ng/mL	10.0	1.0	1.0
20.0 ng/mL	10.0	1.0	2.0
50.0 ng/mL	10.0	1.0	5.0
100.0 ng/mL	10.0	1.0	10.0
200.0 ng/mL	10.0	1.0	20.0
1000.0 ng/mL	10.0	1.0	100.
10.0 µg/mL	10.0	0.5	500.
10.0 µg/mL	10.0	2.0	2000.

* The equivalent sample concentration is based on fortifying a 10-g soil sample with the specified aliquot of spiking solution.

3. Preparation of *cis*- and *trans*-CAAC Calibration Solutions

- a. Solutions for calibration are prepared by adding approximately 10 mL of isooctane and 100 µL of acetic acid to a 100-mL volumetric flask. The flask is agitated to allow acetic acid to contact the glass surface. The appropriate aliquot of the initial solution from Section G.2.a. and acetone is then added and diluted to volume with isooctane as follows:

Effective Date: July 26, 1995

GRM 94.17

Aliquot of 10.0 µg/mL Soln. mL	Aliquot of Acetone mL	Final Soln. Volume mL	Calibrn. Soln. Final Conc. ng/mL	Equivalent Sample Conc. ^a ng/g
0.200	4.8	100	20.0	1.00
0.500	4.5	100	50.0	2.50
2.00	3.0	100	200.	10.0
5.00	0.0	100	500.	25.0

^a The equivalent sample concentration of the standard is based on taking the 10-g soil extract to a final volume of 0.5 mL.

- b. A 2.00 and 4.00 ng/mL calibration solution are prepared by adding approximately 10 mL of isooctane and 100 µL of acetic acid to a 100-mL volumetric flask. The flask is agitated to allow acetic acid to contact the glass surface. The appropriate aliquot of the 200 ng/mL solution from Section G.3.a. above and acetone is then added and diluted to volume with isooctane as follows:

Aliquot of 200 ng/mL mL	Aliquot of Acetone mL	Final Soln. Volume mL	Calibrn. Soln. Final Conc. ng/mL	Equivalent Sample Conc. ^a ng/g
1.00	5.0	100	2.0	0.100
2.00	5.0	100	4.0	0.200

^a The equivalent sample concentration of the standard is based on taking the 10-g soil extract to a final volume of 0.5 mL.

H. Gas Chromatography/Mass Spectrometry

1. Column

Install the splitless column inlet liner (Section E.5.) and the capillary column (Section E.1.) in the split/splitless injection port of the GC/MSD following the manufacturer's recommended procedure.

2. Typical Operating Conditions

Instrumentation: Hewlett-Packard Model 5890 (II) Gas Chromatograph
Hewlett-Packard Model 5971A Mass Selective Detector
Hewlett-Packard Model G1034C Data System Software

Column: J&W Scientific fused silica capillary
Durabond-5 liquid phase
30 m x 0.25 mm i.d.
0.25 µm film thickness

Temperatures:

Column 45 °C for 1.0 min
45 °C to 220 °C at 10 °C/min

Injector 230 °C
Interface 300 °C

Effective Date: July 26, 1995

Carrier Gas:	helium
Head Pressure	50 kPa
Linear Velocity	approximately 40 cm/sec at an oven temperature of 140 °C
Injection Mode:	splitless
Purge Delay	0.5 min
Splitter Flow	50 mL/min
Septum Purge	1.0 mL/min
Injection Volume:	1 µL
Detector:	electron impact ionization with selected ion monitoring
Calibration Program	maximum sensitivity autotune (Note N.3.)
Electron Multiplier	1647 volts (tune voltage plus 200)
Ions Monitored:	
<i>cis</i> -CAAC TBDMSE	<i>m/z</i> 163 (quantitation), <i>m/z</i> 165 (confirmation)
<i>trans</i> -CAAC TBDMSE	<i>m/z</i> 163 (quantitation), <i>m/z</i> 165 (confirmation)
Dwell Time	100 msec

Typical mass spectra of *cis*- and *trans*-CAAC TBDMSE are shown in Figures 1 and 2, respectively. Nominal *m/z* 163 and 165 ions monitored result from loss of the *t*-butyl radical (mass 57) and reflect the isotopic contributions of ³⁵Cl and ³⁷Cl, respectively.

3. Calibration Curves

Typical calibration curves for the determination of *cis*- and *trans*-CAAC in soil are shown in Figures 3 and 4, respectively.

4. Typical Chromatograms

Typical chromatograms of a standard, control sample, and a 0.20 ng/g recovery sample for *cis*- and *trans*-CAAC in soil are shown in Figures 5-10.

I. Determination of Recovery of *cis*- and *trans*-CAAC from Soil

To minimize the potential for cross contamination, equipment used to process samples and reusable glassware should be thoroughly rinsed with the 90% acetone/10% 0.1 N hydrochloric acid solution followed by acetone prior to use.

1. Preparation of Recovery Samples

- Weigh 10.0 g of control soil into a series of 45-mL glass vials.
- For preparing fortified samples, use some of the samples as controls and fortify the remaining samples by adding the specified aliquots of the appropriate spiking solutions (Section G.2.d.) in acetone to obtain concentrations ranging from 0.20 to

Effective Date: July 26, 1995

2000 ng/g. A reagent blank, containing no soil sample, should be carried through the method with the samples.

- c. Add 15.0 mL of the 90% acetone/10% 0.1 N hydrochloric acid solution (F.2.a.) to each sample vial and seal with a PTFE-lined cap.
- d. Vortex the samples briefly and sonicate 10-15 seconds.
- e. Shake the samples for a minimum of 30 minutes on a reciprocating shaker at approximately 180 excursions/minute.
- f. Centrifuge each sample for 10 minutes at 2500 rpm.
- g. Carefully decant each extract to a clean 45-mL vial.
- h. Extract each sample a second time by repeating Steps I.1.c., d. and f. Combine the extracts by decanting to the vial in Step I.1.g.
- i. Evaporate the acetone in the extract by placing the vial in an N-Evap evaporator set at 40 °C. Evaporate the sample under nitrogen to approximately 4.0 mL.
- j. Add approximately 20 mL of deionized water to the vial and seal with a PTFE-lined cap.
- k. Vortex and sonicate the samples for 10-15 seconds.
- l. Adjust the sample pH with 0.1 N sodium hydroxide to fall within a range of 6.5 to 8.0 pH units using the following procedure:
 - (1) Determine the pH using an indicator strip (Section E.12.).
 - (2) Add a volume (mL) of 0.1 N sodium hydroxide equal to $240 \times 10^{-\text{pH}}$. For example, if the pH was determined to be three, the amount of 0.1 N sodium hydroxide to add would be 240×10^{-3} or 0.24 mL.
 - (3) Vortex and sonicate the sample for 5 seconds.
 - (4) Determine the pH using an indicator strip (Section E.13.).
 - (5) If the pH falls within the acceptable range proceed with Step I.1.m. If the pH is below the acceptable range add 0.1 N sodium hydroxide dropwise with thorough mixing until an acceptable pH is reached. If the pH is above the acceptable range add 0.1 N hydrochloric acid dropwise with mixing until an acceptable pH is reached. Mixing is critical to ensure that the pH determined represents the total solution.
- m. The samples are then concentrated and purified using the following ion-exchange SPE procedure:
 - (1) Place an ion-exchange (quaternary amine) SPE column (Section E.3.) on the vacuum manifold box.
 - (2) Attach a 75-mL (Section E.6.) reservoir to the top of the column using an SPE column adapter (Section E.4.).
 - (3) Rinse the SPE column and reservoir with approximately 5 mL of methanol. (Do not allow the column bed to dry.)
 - (4) Condition the SPE column with approximately 5 mL of deionized water. (Do not allow the column bed to dry.)

Effective Date: July 26, 1995

- (5) Transfer the sample solution from Step I.1.i(5). to the reservoir and, with the aid of vacuum, pull the sample through the column at a flow rate of approximately 2 mL/min.
 - (6) Rinse the sample vial with approximately 2 mL of deionized water and transfer the rinse to the reservoir. With the aid of vacuum pull the sample through the column at a flow rate of approximately 2 mL/min.
 - (7) Elute the CAAC by passing 5.0 mL of 0.1 N hydrochloric acid solution through the column, collecting the eluent in a 11-mL vial. Discard the SPE column. (Note N.4.)
- n. Add 100 μ L of 2.0 N HCl and 2-3 g of sodium chloride to the 11-mL vial.
 - o. Add 2.5 mL of MTBE to the vial and seal with a PTFE-lined cap. Vortex the vial for 15 seconds and shake the sample for 10 minutes on a reciprocating shaker at approximately 180 excursions/minute.
 - p. Centrifuge the vial for 5 minutes at 2500 rpm.
 - q. The samples are then dried and purified further using the following silica gel SPE procedure:
 - (1) Weigh a 12-mL tube to four significant figures. This tube will be used to collect the eluent at the end of this procedure and the difference between the empty tube weight and the weight at Step I.1.t. will be used to calculate the final sample volume.
 - (2) Place a silica gel SPE column (Section E.2.) on the vacuum manifold box.
 - (3) Rinse the SPE column with approximately 5 mL of MTBE.
 - (4) Condition the SPE column with approximately 5 mL of isooctane. (Do not allow the column bed to dry.)
 - (5) Transfer the MTBE layer (top layer) from the vial in Step I.1.p. to the SPE column and allow the MTBE to pass by gravity flow through the column.
 - (6) Repeat Steps I.1.o. and p. without the shaking procedure and add the MTBE layer to the SPE column. With the aid of vacuum, pull the MTBE through the column at a flow rate of approximately 2 mL/min..
 - (7) Elute the CAAC by passing 10.0 mL of the 0.025% acetic acid in MTBE solution (Section F.2.b.) through the column, collecting the eluent in the 12-mL tube previously weighed. (Note N.4.)
 - r. Add 0.5 mL of isooctane to the eluent, cap the tube with a PTFE-lined cap and vortex the sample for 5 seconds.
 - s. Evaporate the solution at ambient temperature to a volume of approximately 0.25 mL under a gentle flow of nitrogen. (Do not allow sample to evaporate significantly below 0.25 mL.)
 - t. Add 25 μ L of acetone and bring the volume of the sample to approximately 0.5 mL with isooctane by comparison to a set of 12-mL tubes containing a measured volume of 0.5 mL of isooctane. Weigh the tube with sample for use in calculation of final volume. Seal the tube with a PTFE-lined cap after weighing.
 - u. Add approximately 0.1 g of anhydrous sodium sulfate and vortex the sample for 15 seconds.

Effective Date: July 26, 1995

- v. Add 25 μ L of MTBSTFA and vortex the sample for 15 seconds.
- w. Centrifuge the sample for 5 minutes at 2500 rpm.
- x. Transfer the sample to a 2-mL autosampler vial and seal the vial with a cap and crimper.
- y. Transfer 0.5 mL of each of the calibration standards in Section G.3.a. and b. to autosampler vials. Derivatize by adding 25 μ L of MTBSTFA to each vial and seal with a cap and crimper.
- z. Analyze the samples and calibration standards by capillary gas chromatography/mass spectrometry as described in Section H. Samples that demonstrate CAAC levels above 20 ng/g are diluted 100-fold and reanalyzed as follows:
 - (1) Transfer 0.100 mL of the sample to a 10-mL volumetric flask, add 10 μ L of acetic acid, approximately 0.5 mL acetone and dilute to volume with isooctane.
 - (2) Transfer 0.5 mL of the diluted sample to a 2-mL autosampler vial.
 - (3) Add 25 μ L of MTBSTFA and seal the vial with a cap and crimper.
 - (4) Reanalyze as described in Section H.

2. Calculation of Normalized m/z 163 and 165 Sample Peak Area Response

- a. Calculate the final sample weight by subtracting the weight of the empty 12-mL tube (Section I.1.q.) from the weight of the tube and sample (Section I.1.t.).
- b. The final sample volume contains a known volume of 0.025 mL of acetone and an approximate volume of isooctane. The weight contribution of 0.025 mL acetone is subtracted from the sample final weight (Step I.2.a.) giving the weight of isooctane present. The volume of isooctane can then be determined using the density of isooctane. The final sample volume is the addition of the calculated isooctane volume and the 0.025 mL of acetone. Calculate the final sample volume using the following equations:

$$\text{final sample volume (mL)} = \left[\frac{\text{sample final weight} - (0.025 \text{ mL} \times \text{density of acetone})}{\text{density of isooctane}} \right] + 0.025 \text{ mL}$$

$$\text{final sample volume (mL)} = \left[\frac{\text{sample final weight} - (0.025 \text{ mL} \times 0.7899 \text{ g/mL})}{0.687 \text{ g/mL}} \right] + 0.025 \text{ mL}$$

- c. Normalize the sample m/z 163 and 165 peak areas to 0.5 mL using the following equation:

$$\frac{\text{normalized peak area}}{\text{peak area}} = \frac{\text{final sample volume (mL)}}{0.5 \text{ mL}} \times \text{sample peak area}$$

For example, using the data for *cis*-CAAC from Figure 7:

$$\text{d. final sample weight (g)} = 12.22 - 11.74$$

$$\text{final sample weight (g)} = 0.480$$

Effective Date: July 26, 1995

$$e. \text{ final sample volume (mL)} = \left[\frac{0.480 \text{ g} - (0.025 \text{ mL} \times 0.7899 \text{ g/mL})}{0.687 \text{ g/mL}} \right] + 0.025 \text{ mL}$$

$$\text{final sample volume (mL)} = 0.6949$$

$$f. \text{ normalized } m/z \text{ 163 peak area} = \frac{0.6949 \text{ mL}}{0.5 \text{ mL}} \times 415$$

$$\text{normalized } m/z \text{ 163 peak area} = 577$$

$$\text{normalized } m/z \text{ 165 peak area} = \frac{0.6949 \text{ mL}}{0.5 \text{ mL}} \times 152$$

$$\text{normalized } m/z \text{ 165 peak area} = 211$$

3. Calculation of Percent Recovery

- Determine the *m/z* 163 and 165 response areas for both *cis*- and *trans*-CAAC TBDMSE in calibrations standards from Step I.1.y. and z.
- For each standard, calculate the *cis*- and *trans*-CAAC confirmation ratios. The average confirmation ratio for each will be used to confirm the presence of the respective CAAC in the soil samples.

For example, using the data for *cis*-CAAC from Figure 5:

$$\text{Confirmation Ratio} = \frac{\text{peak area of confirmation ion}}{\text{peak area of quantitation ion}}$$

$$\text{Confirmation Ratio} = \frac{\text{peak area at } m/z \text{ 165}}{\text{peak area at } m/z \text{ 163}}$$

$$\text{Confirmation Ratio} = \frac{224}{621}$$

$$\text{Confirmation Ratio} = 0.3607$$

Positive confirmation of the presence of *cis*- and *trans*-CAAC is indicated when the confirmation ratio for the samples is in the range of $\pm 20\%$ of the average found for the respective standards.

- Prepare *cis*- and *trans*-CAAC standard curves by plotting the sample equivalent concentration (ng/g) on the abscissa (x-axis) and the *cis*- and *trans*-CAAC TBDMSE *m/z* 163 peak area on the ordinate (y-axis) as shown in Figures 3 and 4, respectively. Using regression analysis, determine the equation for the curve with respect to the abscissa.

For example, using power regression (3) with the *trans*-CAAC data from Figure 4:

$$Y = \text{constant} \times X^{\text{(exponent)}}$$

$$X = \left(\frac{Y}{\text{constant}} \right)^{1/\text{exponent}}$$

$$\text{trans-CAAC Conc. (ng/g)} = \left(\frac{\text{trans-CAAC TBDMSE peak area}}{\text{constant}} \right)^{1/\text{exponent}}$$

$$\text{trans-CAAC Conc. (ng/g)} = \left(\frac{\text{trans-CAAC TBDMSE peak area}}{3424.0} \right)^{1/1.0147}$$

- d. Determine the net concentration in each recovery sample that does not require dilution (less than or equal to 20 ng/g) by first subtracting the average *cis*- and *trans*-CAAC TBDMSE normalized *m/z* 163 peak area in the control sample from that of the recovery sample. Substitute the peak area obtained into the appropriate equation and solve for the concentration.

For example, using the normalized *m/z* 163 peak areas for *trans*-CAAC data from Figures 9 and 10 and the standard curve equation from Figure 4:

$$\text{trans-CAAC Conc. (ng/g)} = \left(\frac{\text{net trans-CAAC TBDMSE peak area}}{3424.0} \right)^{1/1.0147}$$

$$\text{trans-CAAC Conc. (ng/g)} = \left(\frac{613 - 66}{3424.0} \right)^{1/1.0147}$$

$$\text{trans-CAAC Conc.} = 0.1641 \text{ ng/g}$$

- e. Determine the percent recovery by dividing the net concentration found for each recovery sample by the theoretical concentration added.

$$\text{Recovery} = \frac{\text{Concentration Found}}{\text{Concentration Added}} \times 100\%$$

$$\text{Recovery} = \frac{0.1641 \text{ ng/g}}{0.200 \text{ ng/g}} \times 100\%$$

$$\text{Recovery} = 82\%$$

- f. Determine the diluted concentration in each recovery sample fortified at levels above 20 ng/g by substituting the *cis*- and *trans*-CAAC TBDMSE normalized *m/z* 163 peak area into the appropriate standard curve and solve for the diluted concentration as described in Section I.3.d.

Effective Date: July 26, 1995

- g. Determine the concentration in each diluted recovery sample by multiplying the diluted concentration by the dilution factor of 100.

$$\text{CAAC Conc. (ng/g)} = \text{diluted CAAC Conc. (ng/g)} \times 100$$

- h. Determine the percent recovery for recovery samples fortified above 20 ng/g as described in Section I.3.e.

The average of all the recovery samples in a given analytical set can be used to correct sample results in the set for method efficiency.

J. Determination of *cis*- and *trans*-CAAC in Soil

1. Prepare reagent blank, control, recovery, and treated samples as described in Section I.1.
2. Prepare standard calibration curves for *cis*- and *trans*-CAAC and determine the percentage recovery for each as described in Section I.3.
3. Determine the concentration of *cis*- and *trans*-CAAC in each treated sample by substituting the *cis*- and *trans*-CAAC TBDMSE normalized *m/z* 163 peak areas obtained into the respective equations for the standard calibration curves, and calculate the uncorrected residue results.

For example, using the *cis*-CAAC data from Figures 3 and 7, the uncorrected concentration is calculated as follows:

$$\text{cis-CAAC Conc. (ng/g)} = \left(\frac{\text{cis-CAAC TBDMSE peak area}}{\text{constant}} \right)^{1/\text{exponent}}$$

$$\text{cis-CAAC Conc. (ng/g)} = \left(\frac{577}{3243.3} \right)^{1/1.0244}$$

$$\text{cis-CAAC Conc.} = 0.1854 \text{ ng/g}$$

4. Samples that exceed a *cis*- or *trans*-CAAC concentration of 20 ng/g are diluted and reanalyzed as described in Section I.1.z.
5. Determine the concentration of *cis*- and *trans*-CAAC in each diluted, treated sample as described in Section I.3.f and g.

K. Determination of Soil Moisture

1. Weigh 10.00 g of soil into an aluminum or glass container.
2. Place the sample in an oven at approximately 130 °C and allow to dry for a minimum of 16 hours.

Effective Date: July 26, 1995

3. Remove the sample from the oven, place in a desiccator until the sample has cooled to ambient temperature, and then re-weigh.
4. Calculate the percent moisture on a dry weight basis as follows:

$$\begin{aligned} \text{Percent Moisture} &= \frac{\text{soil moisture weight (g)}}{\text{dehydrated soil weight (g)}} \times 100 \\ &= \frac{(\text{soil weight before drying} - \text{soil weight after drying})}{\text{soil weight after drying}} \times 100 \end{aligned}$$

L. Determination of Corrected *cis*- and *trans*-CAAC in Soil

1. Determine the *cis*- and *trans*-CAAC concentration in the soil samples as described in Section J.
2. Determine the soil moisture as described in Section K.
3. Determine the corrected *cis*- and *trans*-CAAC concentrations in soil samples as follows:

$$\text{CAAC Corrected Conc. (ng/g)} = \left(\frac{\text{CAAC Conc. (ng/g)}}{\% \text{ Recovery}} \right) \left(1 + \frac{\% \text{ Moisture}}{100} \right)$$

M. Results and Discussion

1. Method Validation

a. Recovery Levels and Precision

A method validation study was conducted to determine the recovery levels and the precision of the method for *cis*- and *trans*-CAAC in soil. The results are summarized in Tables I and II.

Recovery values of *cis*-CAAC from soil samples fortified over the concentration range 0.20 to 2000 ng/g averaged 80% with one standard deviation equal to 7%. Recovery values of *trans*-CAAC from soil samples fortified over the concentration range 0.20 to 2000 ng/g averaged 84% with one standard deviation equal to 6%.

b. Standard Curve Fit

The average correlation coefficient (r^2) for the power least squares regression equations describing the detector response as a function of the standard calibration curve concentration was greater than 0.997 for both *cis*- and *trans*-CAAC.

c. Calculated Limits of Quantitation and Detection

Following established guidelines (4), the limits of quantitation (LOQ) and detection (LOD) were calculated using the standard deviation from the 0.20 ng/g

recovery results. The LOQ was calculated as ten times the standard deviation (10s) and the LOD was calculated as three times the standard deviation (3s) of the results of the analysis of nine samples. The results are summarized in Tables III and IV.

For *cis*- and *trans*-CAAC, the calculated statistics support an LOQ of 0.22 and 0.19 ng/g, respectively. Both values are in agreement with the targeted method LOQ of 0.20 ng/g. Results should not be quantified, however, at levels below which no recovery samples have been analyzed.

For *cis*- and *trans*-CAAC, calculated statistics support an LOD of 0.07 and 0.06 ng/g, respectively.

2. Confirmation of Residue Identity

Confirmation of the presence of residues is described in Section I.3.b. For *cis*- and *trans*-CAAC, confirmation is by comparison of the chromatographic retention time as well as the peak area ratios resulting from mass selective ion monitoring. Positive confirmation of *cis*- and *trans*-CAAC is indicated when the confirmation ratio for the sample is in the range of $\pm 20\%$ of the average determined for the standards.

3. Assay Time

A typical analytical run would consist of a minimum of four standards encompassing the expected range of sample concentrations, a reagent blank, a control (a non-fortified sample), a minimum of two fortified controls (one of which must be at the LOQ), and ten samples. This typical analytical run could be prepared in approximately 10 hours, and the chromatographic analysis take place the same evening.

There are several acceptable "stopping points" in the method where sample preparation (Section I.) may be suspended without deleterious effects on the sample analysis. These are indicated below:

- a. Step I.1.o.
- b. Step I.1.q.(7). If the samples are to be stored overnight, the vials should be sealed with PTFE-lined caps.
- c. Step I.1.r.

4. Standardization of SPE Elution Profiles

Variation in the ion-exchange and silica gel SPE columns may influence the elution profile of *cis*- and *trans*-CAAC. If method performance degrades significantly with a change in the lot number of SPE columns used, an elution profile of the columns should be carried out. The following procedures can be used:

a. Ion-exchange SPE Profile

- (1) In an 8-mL vial, add 25 μ L of the 10 μ g/mL *cis*- and *trans*-CAAC spiking solution (Section G.2.a.) to 5 mL of DI water.
- (2) Place an ion-exchange SPE column on the vacuum manifold box.
- (3) Rinse the SPE column with approximately 5 mL of methanol.

Effective Date: July 26, 1995

- (4) Condition the SPE column with approximately 5 mL of DI water. (Do not allow the column bed to dry.)
- (5) Transfer the sample solution from Step M.4.a.(1) to the SPE column and, with the aid of vacuum, slowly pull the sample through the column.
- (6) Rinse the 8-mL vial with 2 mL of DI water, transfer the rinse to the SPE column and, with the aid of vacuum, slowly pull the rinse through the column.
- (7) Elute the *cis*- and *trans*-CAAC with 10.0 mL of 0.1 N hydrochloric acid solution, collecting 1-mL aliquots in 8-mL vials.
- (8) For each fraction collected, add 100 μ L 2.0 N HCl, approximately 2 g of sodium chloride and extract two times with 2.5 mL of 0.025% acetic acid in MTBE.
- (9) Combine both extracts in an 8-mL vial, add approximately 2 mL of isooctane and evaporate at ambient temperature under a gentle flow of nitrogen to approximately 0.5 mL.
- (10) Adjust the volume to 1.0 mL with isooctane by comparison to two 8-mL vials each containing a measured volume of 1.0 mL isooctane. Add 25 μ L of MTBSTFA and vortex the sample for 10 seconds.
- (11) Transfer to an autosampler vial and seal the vial with a cap and crimper.
- (12) Perform Steps I.1.y. and z.
- (13) Calculate the percentage recoveries for each analyte as described in Section I.3. Use actual standard concentrations in determination of the standard curve equation and fraction recoveries.

Evaluation of results:

If less than 90% of the analytes are recovered in the first five 1-mL fractions, the method recoveries may fall below acceptable values. In such a case it would be imperative to check the validity of reagents, particularly the 0.1 N hydrochloric acid solution used to elute the column. If similar results are obtained after checking reagents, a second lot of SPE columns should be evaluated.

b. Silica gel SPE Profile

- (1) In an 11-mL vial, add 25 μ L of the 10 μ g/mL *cis*- and *trans*-CAAC spiking solution (Section G.2.a.) to 5 mL of 0.1 N HCl.
- (2) Proceed with Steps I.1.n. through p.
- (3) Place a silica gel SPE column on the vacuum manifold box.
- (4) Rinse the SPE column with 5 mL of MTBE.
- (5) Condition the SPE column with 5 mL of isooctane. (Do not allow the column bed to dry.)
- (6) Transfer the MTBE layer (top layer) of the sample solution from Step M.4.b.(2) to the SPE column and, with the aid of vacuum, slowly pull the sample through the column. (Do not allow the column bed to dry.)
- (7) Repeat Steps I.1.o. and p. without the shaking procedure and transfer the

Effective Date: July 26, 1995

MTBE to the SPE column. With the aid of vacuum, slowly pull the sample through the column.

- (8) Elute the analytes with 10 mL of the 0.025% acetic acid in MTBE solution, collecting 1-mL aliquots in 8-mL vials.
- (9) Add 2 mL of isooctane to each of the elution fractions in Step M.4.b.(8).
- (10) Concentrate the fractions at ambient temperature under a gentle flow of nitrogen to approximately 0.5 mL.
- (11) Adjust the volume of the fractions to 1.0 mL with isooctane by comparison to two 8-mL vials each containing a measured volume of 1.0 mL isooctane. Add 25 μ L of MTBSTFA and vortex the samples for 10 seconds.
- (12) Transfer to an autosampler vial and seal the vial with a cap and crimper.
- (13) Perform Steps I.1.y. and z.
- (14) Calculate the percentage recoveries for each analyte as described in Section I.3. Use actual standard concentrations in determination of the standard curve equation and fraction recoveries.

Evaluation of results:

If less than 90% of the analytes are recovered in the ten 1-mL fractions, the method recoveries may fall below acceptable values. In such a case it would be imperative to check the validity of reagents, particularly the 0.025% acetic acid in MTBE solution used to elute the column. If similar results are obtained after checking reagents, a second lot of SPE columns should be evaluated. The presence of acetic acid in the MTBE is critical to eluting the analytes from the silica gel, increasing the eluting solution to 0.04% acetic acid in MTBE may be evaluated if the preceding suggestions fail to give acceptable results.

N. Notes

1. Equipment, glassware, materials, reagents, and chemicals considered to be equivalent to those specified may be substituted with the understanding that their performance must be confirmed by appropriate tests. Common laboratory supplies are assumed to be readily available and are, therefore, not listed.
2. The filters are used in the carrier gas supply lines to purify the helium entering the gas chromatograph.
3. Several tuning, or calibration, options are available for the Model 597X series of MSDs. The "Maximum Sensitivity Autotune" feature was found to consistently yield approximately 5-10 times the sensitivity compared to that of the "Standard Autotune".
4. Depending on the number of samples being prepared, one may elute the CAAC from each SPE column individually, using either gravity-feed or pressurized elution, or as a group, using the vacuum manifold box.

Effective Date: July 26, 1995

O: References

1. Certificate Of Analysis, Aldrich Chemical Co., January 20, 1994.
2. Certificate Of Analysis, Aldrich Chemical Co., January 20, 1994.
3. *HP-41C/41CV Standard Applications Handbook*, Hewlett-Packard Publication No. 00041-90402, 1982, pp 42-48.
4. Keith, L.H.; Crummett, W.B.; Deegan, J.; Libby, R.A.; Taylor, J.T.; Wentler, G., "Principles of Environmental Analysis", *Anal. Chem.*, 55, 2210-2218 (1983).

The information herein is presented in good faith, but no warranty, expressed or implied, is given nor is freedom from any patent owned by DowElanco or by others to be inferred. In the hands of qualified personnel, the procedures are expected to yield results of sufficient accuracy for their intended purposes, but recipients are cautioned to confirm the reliability of their techniques, equipment, and standards by appropriate tests. Anyone wishing to reproduce or publish the material in whole or in part should request written permission from DowElanco.

Effective Date: July 26, 1995

Table I. Recovery of *cis*-CAAC from Soil

Sample Number	Date of Analysis	<i>cis</i> -CAAC, ng/g		Percent Recovery
		Added	Found ^a	
14767201	08-Sept-1994	0.000	0.0223	--
14767202	08-Sept-1994	0.000	0.0177	--
14767203	08-Sept-1994	0.000	0.0194	--
14767201	15-Sept-1994	0.000	0.0284	--
14767202	15-Sept-1994	0.000	0.0251	--
14767203	15-Sept-1994	0.000	0.0231	--
14767202	21-Oct-1994	0.000	0.0232	--
14767203	21-Oct-1994	0.000	0.0174	--
14767202	21-Oct-1994	0.060	0.051	NA ^b
14767203	21-Oct-1994	0.060	0.052	NA
14767201	08-Sept-1994	0.200	0.1457	73
14767201	08-Sept-1994	0.200	0.1574	79
14767202	08-Sept-1994	0.200	0.1236	62
14767202	08-Sept-1994	0.200	0.1690	84
14767203	08-Sept-1994	0.200	0.1659	83
14767203	08-Sept-1994	0.200	0.1618	81
14767201	15-Sept-1994	0.200	0.1892	95
14767202	15-Sept-1994	0.200	0.1871	94
14767203	15-Sept-1994	0.200	0.1869	93
14767201	08-Sept-1994	1.000	0.7508	75
14767202	08-Sept-1994	1.000	0.7371	74
14767203	15-Sept-1994	2.000	1.610	81
14767202	15-Sept-1994	2.000	1.465	73
14767202	08-Sept-1994	5.000	4.112	82
14767203	08-Sept-1994	5.000	4.320	86
14767202	15-Sept-1994	10.00	7.862	79
14767203	15-Sept-1994	10.00	7.821	78
14767201	08-Sept-1994	20.00	15.41	77
14767203	08-Sept-1994	20.00	15.33	77
14767201	15-Sept-1994	100.00	81.27	81
14767203	15-Sept-1994	100.00	83.64	84
14767202	15-Sept-1994	500.0	398.9	80
14767203	15-Sept-1994	500.0	366.4	73
14767201	15-Sept-1994	2000.0	1486	74
14767202	15-Sept-1994	2000.0	1604	80
\bar{x} =				80
s =				7
n =				25

^a Fortified samples corrected for average of controls in set.

^b NA = not applicable. The residue was below the 0.20 ng/g LOQ.

Effective Date: July 26, 1995

Table II. Recovery of *trans*-CAAC from Soil

Sample Number	Date of Analysis	<i>trans</i> -CAAC, ng/g		Percent Recovery
		Added	Found ^a	
14767201	08-Sept-1994	0.000	0.0204	--
14767202	08-Sept-1994	0.000	0.0174	--
14767203	08-Sept-1994	0.000	0.0164	--
14767201	15-Sept-1994	0.000	0.0197	--
14767202	15-Sept-1994	0.000	0.0187	--
14767203	15-Sept-1994	0.000	0.0168	--
14767202	21-Oct-1994	0.000	0.0126	--
14767203	21-Oct-1994	0.000	0.0121	--
14767202	21-Oct-1994	0.060	0.051	NA ^b
14767203	21-Oct-1994	0.060	0.058	NA
14767201	08-Sept-1994	0.200	0.1641	82
14767201	08-Sept-1994	0.200	0.1738	87
14767202	08-Sept-1994	0.200	0.1291	65
14767202	08-Sept-1994	0.200	0.1776	89
14767203	08-Sept-1994	0.200	0.1785	89
14767203	08-Sept-1994	0.200	0.1750	88
14767201	15-Sept-1994	0.200	0.1861	93
14767202	15-Sept-1994	0.200	0.1806	90
14767203	15-Sept-1994	0.200	0.1961	98
14767201	08-Sept-1994	1.000	0.8135	81
14767202	08-Sept-1994	1.000	0.8074	81
14767203	15-Sept-1994	2.000	1.695	85
14767202	15-Sept-1994	2.000	1.493	75
14767202	08-Sept-1994	5.000	4.130	83
14767203	08-Sept-1994	5.000	4.357	87
14767202	15-Sept-1994	10.00	8.113	81
14767203	15-Sept-1994	10.00	8.239	82
14767201	08-Sept-1994	20.00	16.06	80
14767203	08-Sept-1994	20.00	16.15	81
14767201	15-Sept-1994	100.00	83.20	83
14767203	15-Sept-1994	100.00	86.10	86
14767202	15-Sept-1994	500.0	408.1	82
14767203	15-Sept-1994	500.0	393.2	79
14767201	15-Sept-1994	2000.0	1562.	78
14767202	15-Sept-1994	2000.0	1652.	83
\bar{x} =				84
s =				6
n =				25

^a Fortified samples corrected for average of controls in set.

^b NA = not applicable. The residue was below the 0.20 ng/g LOQ.

Effective Date: July 26, 1995

GRM 94.17

Table III. Calculated Limits of Detection and Quantitation for the Determination of *cis*-CAAC in Soil

Sample Number	Date of Analysis	<i>cis</i> -CAAC, ng/g	
		Added	Found
14767201	08-Sept-1994	0.200	0.1457
14767201	08-Sept-1994	0.200	0.1574
14767202	08-Sept-1994	0.200	0.1236
14767202	08-Sept-1994	0.200	0.1690
14767203	08-Sept-1994	0.200	0.1659
14767203	08-Sept-1994	0.200	0.1618
14767201	15-Sept-1994	0.200	0.1892
14767202	15-Sept-1994	0.200	0.1871
14767203	15-Sept-1994	0.200	0.1869
		\bar{x} =	0.1652
		s =	0.0216
		LOD ^a (3s) =	0.065
		LOQ ^b (10s) =	0.22

^a LOD = Limit of Detection.

^b LOQ = Limit of Quantitation.

Effective Date: July 26, 1995



GRM 94.17

Table IV. Calculated Limits of Detection and Quantitation for the Determination of *trans*-CAAC in Soil

Sample Number	Date of Analysis	<i>trans</i> -CAAC, ng/g	
		Added	Found
14767201	08-Sept-1994	0.200	0.1641
14767201	08-Sept-1994	0.200	0.1738
14767202	08-Sept-1994	0.200	0.1291
14767202	08-Sept-1994	0.200	0.1776
14767203	08-Sept-1994	0.200	0.1785
14767203	08-Sept-1994	0.200	0.1750
14767201	15-Sept-1994	0.200	0.1861
14767202	15-Sept-1994	0.200	0.1806
14767203	15-Sept-1994	0.200	0.1961
		\bar{x} =	0.1734
		s =	0.0188
		LOD ^a (3s) =	0.056
		LOQ ^b (10s) =	0.19



^a LOD = Limit of Detection.

^b LOQ = Limit of Quantitation.



cis-CAAC *t*-Butyldimethylsilyl Ester
Formula: C₉H₁₇ClO₂Si
Molecular Weight: 220


Figure 1. Mass Spectrum of *cis*-CAAC *t*-Butyldimethylsilyl Ester



***trans*-CAAC *t*-Butyldimethylsilyl Ester**
Formula: C₉H₁₇ClO₂Si
Molecular Weight: 220

Figure 2. Mass Spectrum of *trans*-CAAC *t*-Butyldimethylsilyl Ester

Calibration Curve




Equivalent <i>cis</i> -CAAC Conc. ng/g	<i>m/z</i> 163 Peak Area Position in Sequence	
	Start of Sequence	End of Sequence
0.10	283	323
0.20	621	652
1.00	3152	3206
2.50	8482	8449
10.0	34470	33903
25.0	86524	88377

Power Regression Equation: $X = \left[\frac{Y}{3243.3} \right]^{1/1.0244}$

Coefficient of Determination (r^2): 0.9997

Figure 3. Typical Calibration Curve for the Determination of *cis*-CAAC in Soil Samples

Effective Date: July 26, 1995




Equivalent <i>trans</i> -CAAC Conc. ng/g	<i>m/z</i> 163 Peak Area Position in Sequence	
	Start of Sequence	End of Sequence
0.10	325	340
0.20	654	689
1.00	3323	3348
2.50	8868	8802
10.0	35818	35134
25.0	88737	90484

Power Regression Equation: $X = \left[\frac{Y}{3424.0} \right]^{1/1.0147}$

Coefficient of Determination (r^2): 0.9999

Figure 4. Typical Calibration Curve for the Determination of *trans*-CAAC in Soil Samples

Effective Date: July 26, 1995



Data File : 0201002.D
ALS Bottle : 2
Date : 8 Sep 94 3:04 pm
Data Path : C:\CHEMPC\DATA\F090894A\
Instrument : GC/MSD - GC serial#3126A36485

Sample Name: 4 ng/ml cis/trans 3-Chloroacrylic acid
Sample Info:
Operator : SEF

cis-3-Chloroacrylic acid Retention Time: 10.79
PEAK AREA (M/Z 163) : 621
PEAK AREA (M/Z 165) : 224


cis-3-Chloroacrylic Acid:
RATIO OF M/Z 165/163: 0.3607

Equivalent cis-CAAC Concentration: 0.20 ng/g

Average Standard Confirmation Ratio: 0.3610

Figure 5. Typical Chromatogram of a 4.0 ng/mL Standard, Equivalent to 0.20 ng/g cis-CAAC in Soil

Effective Date: July 26, 1995



Data File : 1001010.D
ALS Bottle : 10
Date : 8 Sep 94 6:18 pm
Data Path : C:\CHEMPC\DATA\F090894A\
Instrument : GC/MSD - GC serial#3126A36485


Sample Name: control 14767202
Sample Info:
Operator : SEF

cis-3-Chloroacrylic acid Retention Time: 10.79
PEAK AREA (M/Z 163) : 43
PEAK AREA (M/Z 165) : 28
cis-3-Chloroacrylic Acid:
RATIO OF M/Z 165/163: 0.6512

Vial Weight Empty: 11.82 g	Normalized Peak Area (m/z 163): 52
Vial Weight with Sample: 12.24 g	Normalized Peak Area (m/z 165): 34
Final Sample Weight: 0.42 g	cis-CAAC Concentration: 0.018 ng/g
Final Sample Volume: 0.608 mL	Average Std. Confirmation Ratio: 0.3610

Figure 6. Typical Chromatogram of a Control Soil Sample for the Determination of cis-CAAC

Effective Date: July 26, 1995



Data File : 1701017.D
ALS Bottle : 17
Date : 8 Sep 94 9:06 pm
Data Path : C:\CHEMPC\DATA\F090894A\
Instrument : GC/MSD - GC serial#3126A36485

Sample Name: 0.2 ng/ml spiked control
Sample Info:
Operator : SEF


cis-3-Chloroacrylic acid Retention Time: 10.79
PEAK AREA (M/Z 163) : 415
PEAK AREA (M/Z 165) : 152

cis-3-Chloroacrylic Acid:
RATIO OF M/Z 165/163: 0.3663

Vial Weight Empty: 11.74 g	Normalized Peak Area (m/z 163): 577
Vial Weight with Sample: 12.22 g	Normalized Peak Area (m/z 165): 211
Final Sample Weight: 0.48 g	cis-CAAC Concentration: 0.1690 ng/g
Final Sample Volume: 0.695 mL	Average Std. Confirmation Ratio: 0.3610
	Recovery: 84%

Figure 7. Typical Chromatogram of a Control Soil Sample Fortified with 0.20 ng/g cis-CAAC

Effective Date: July 26, 1995



Data File : 0201002.D
ALS Bottle : 2
Date : 8 Sep 94 3:04 pm
Data Path : C:\CHEMPC\DATA\F090894A\
Instrument : GC/MSD - GC serial#3126A36485

Sample Name: 4 ng/ml cis/trans 3-Chloroacrylic acid
Sample Info:
Operator : SEF

trans-3-Chloroacrylic acid Retention Time: 10.01
PEAK AREA (M/Z 163) : 654
PEAK AREA (M/Z 165) : 236


trans-3-Chloroacrylic Acid:
RATIO OF M/Z 165/163: 0.3609

Equivalent *trans*-CAAC Concentration: 0.20 ng/g

Average Standard Confirmation Ratio: 0.3617

Figure 8. Typical Chromatogram of a 4.0 ng/mL Standard, Equivalent to 0.20 ng/g *trans*-CAAC in Soil

Effective Date: July 26, 1995



Data File : 0901009.D
ALS Bottle : 9
Date : 8 Sep 94 5:53 pm
Data Path : C:\CHEMPC\DATA\F090894A\
Instrument : GC/MSD - GC serial#3126A36485

Sample Name: control 14767201
Sample Info:
Operator : SEF


trans-3-Chloroacrylic acid Retention Time: 10.01
PEAK AREA (M/Z 163) : 52
PEAK AREA (M/Z 165) : 25

trans-3-Chloroacrylic Acid:
RATIO OF M/Z 165/163: 0.4808

Vial Weight Empty: 11.91 g	Normalized Peak Area (m/z 163): 66
Vial Weight with Sample: 12.35 g	Normalized Peak Area (m/z 165): 32
Final Sample Weight: 0.44 g	trans-CAAC Concentration: 0.020 ng/g
Final Sample Volume: 0.637 mL	Average Std. Confirmation Ratio: 0.3617

Figure 9. Typical Chromatogram of a Control Soil Sample for the Determination of trans-CAAC

Effective Date: July 26, 1995



Data File : 1301013.D
ALS Bottle : 13
Date : 8 Sep 94 7:30 pm
Data Path : C:\CHEMPC\DATA\F090894A\
Instrument : GC/MSD - GC serial#J126A36485

Sample Name: 0.2 ng/ml spiked control
Sample Info:
Operator : SEF

trans-3-Chloroacrylic acid Retention Time: 10.01
PEAK AREA (M/Z 163) : 460
PEAK AREA (M/Z 165) : 168

trans-3-Chloroacrylic Acid:
RATIO OF M/Z 165/163: 0.3652

Vial Weight Empty: 11.79 g	Normalized Peak Area (m/z 163): 613
Vial Weight with Sample: 12.25 g	Normalized Peak Area (m/z 165): 224
Final Sample Weight: 0.46 g	trans-CAAC Concentration: 0.1641 ng/g
Final Sample Volume: 0.666 mL	Average Std. Confirmation Ratio: 0.3617
	Recovery: 82%

Figure 10. Typical Chromatogram of a Control Soil Sample Fortified with 0.20 ng/g trans-CAAC

Appendix C
Characterization of Soils

DOWELANCO

Residue Research Sample Submission

Grower and Location WILKISON & Associates		Address 3506 EXCHANGE SUITE A		Phone #		Page of	
City NAPLES		State FL		Zip 33942		No. of Samples 1	
Total No. Boxes 1		Date Shipped 3-7-94					
Carrier		Waybill No.		Preparation of Samples for Shipment:		<input type="checkbox"/> Ambient <input type="checkbox"/> Refrigerated <input type="checkbox"/> Frozen (Dry Ice) <input type="checkbox"/> Frozen (No Dry Ice) <input type="checkbox"/> Dry <input type="checkbox"/> Other	
Part Submitted		Sample Type		Variety		Date of Planting	
Crop Stage		Soil Texture		<input type="checkbox"/> Fine <input type="checkbox"/> Medium <input type="checkbox"/> Coarse			
Protocol No. RES94071		Experiment No.		Account No. / CHEMTRAC		Date of Treatment	
I certify that the samples were collected as specified in the protocol.							
Signed _____				Date _____			
Remarks							

Internal Use Only	Sample No.	Trt No.	Rep No.	Sub Samples	Tube Diameter	Depth	Project Treatment Formulation	Rate	Date of Sampling
14762201	M465A						Soil Samples From test site. Sample #15 were assigned 6-13-94 by Dave Miller. GSP Soil Characterization done by A-L Labs. Document attached. See water Submission for copy of original submission. For any information prior to 6-13-94 See Notes G406024 in Lynn Kennedy possession S.E. Fisher 6-13-94		
14762202	M465B								
14762203	M465E								
Samples stored in chiller in 004A prior to tracking done with DSTAR S.E. Fisher 6-13-94									

147672 SGN		ANALYTICAL CHEMISTRY USE ONLY							
Date Received		By Whom		Inspected by		Date		Condition of Samples	
RDG 3/10/94									



A & L GREAT LAKES LABORATORIES, INC.


3505 Conestoga Drive • Fort Wayne, Indiana 46808-4413 • Phone 219-483-4759

GLP SOIL CHARACTERIZATION

Report Number: F94102-122
Date of Report: 05/20/94

Protocol Number:
ENV 93068

DOWELANCO
JACK R. MILLER
R & D BUILDING A-2/763
9410 ZIONSVILLE ROAD
INDIANAPOLIS IN 46268-1053



Sample ID:	A&L Check	M463	M464	M465-A	M465-E
Lab Number:	50660	50661	50662	50663	50664
pH	5.5	7.6	7.8	7.7	8.5
CEC (meq/100g)	8.46	8.35	12.74	2.20	0.29
O.M. (%)	1.86	0.71	1.80	1.15	0.11
WHC (%) @ 1/10 Bar	23.34	35.24	34.80	22.78	17.57
WHC (%) @ 1/3 Bar	20.80	16.46	22.77	2.86	2.23
WHC (%) @ 1 Bar	15.41	8.83	16.42	2.62	2.11
WHC (%) @ 15 Bar	6.27	4.93	8.92	2.46	1.94
Sand (%)	38.0	34.0	30.0	96.0	96.0
Silt (%)	37.2	53.2	41.2	1.2	1.2
Clay (%)	24.8	12.8	28.8	2.8	2.8
Soil Classification	Loam	Silt Loam	Clay Loam	Sand	Sand
Bulk Density (g/cc)	1.52	1.22	1.30	1.50	1.60

Verified By: Haylene J. Matthews Date: 5-20-94



A & L GREAT LAKES LABORATORIES, INC.


3505 Conestoga Drive • Fort Wayne, Indiana 46808-4413 • Phone 219-483-4759

GLP SOIL CHARACTERIZATION

Report Number: F94102-122
Date of Report: 05/20/94

Protocol Number:
ENV 93068

DOWELANCO
JACK R. MILLER
R & D BUILDING A-2/763
9410 ZIONSVILLE ROAD
INDIANAPOLIS IN 46268-1053



Sample ID:

Lab Number:

pH

CEC (meq/100g)

O.M. (%)

WHC (%) @ 1/10 Bar

WHC (%) @ 1/3 Bar

WHC (%) @ 1 Bar

WHC (%) @ 15 Bar

Sand (%)

Silt (%)

Clay (%)

Soil Classification

Bulk Density (g/cc)

	50665	50666	50667	50668
pH	6.6	7.4	7.8	
CEC (meq/100g)	6.01	5.15	6.80	
O.M. (%)	2.78	1.53	0.44	
WHC (%) @ 1/10 Bar	25.48	31.00	27.64	
WHC (%) @ 1/3 Bar	5.07	11.39	15.74	20.77
WHC (%) @ 1 Bar	3.87	7.23	9.92	16.03
WHC (%) @ 15 Bar	2.56	3.43	4.66	6.45
Sand (%)	94.0	60.0	56.0	
Silt (%)	1.2	31.2	29.2	
Clay (%)	4.8	8.8	14.8	
Soil Classification	Sand	Sandy Loam	Sandy Loam	
Bulk Density (g/cc)	1.39	1.35	1.50	

Verified By: Baylene J. Matthias

Date: 5-20-94