## Cover Sheet for

# **ENVIRONMENTAL CHEMISTRY METHOD**

Pestcide Name: Naptalam

**MRID** #: 400691-01

Matrix: Soil

Analysis: HPLC/ELCD

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method is not an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews all Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at <a href="mailto:flynt.elizabeth@epa.gov">flynt.elizabeth@epa.gov</a>.

UNIROYAL Chemical Crop Protection Research and Development Chemistry Section	Project No. 8627
TITLE: ANALYTICAL METHOD FOR DETERMINING ALA	MAP RESIDUE TV SOTT
(For SCIL DISSIPATION STUDY, Project	No. 8552)
This study meets the requirements for 40 CFR Pa	art 160
Submitter UNIROYAL chemical Crop Frote	ection
Sponsor UNIROYAL Chemical Crop Prote	ection
Stricy Director W Salonews	Date 1-13-87
Notebook and Page Number(s): AC-917 40-62	,
Raw Data may be obtained from the sponsor upon r	request.
Report Approved by:	
James C. Mart	Date: 1, 13 19 97
J.V. Mertz	(000)
Section Manager	Date: 1/11/57
H.M. Balba	<u> </u>
SIMPANY Larer II FALL	T. Stiffy 1

#### QUALITY ASSURANCE STATEMENT

•	Project Number	8627
	•	
This report has been reviewed by the Quality Assurance Unit. The final report reflect the raw data generated during the methods and procedures used to generate the		ction curately dy and the
This review was completed on 1/K-187	and 1/20/9	57
· · ·	•	
	•	,
Michaellege Signature, Q.A. Officer		
1/20/87		
Date		
	•	

Project No. 8627

#### Principle

ALANAP, also commonly known as the sodium salt of naptalam (2-[(1-naphtalenylamino)commonyl]benzoic acid, sodium salt),

CAS No. 132-66-1, is hydrolyzed in a boiling sodium hydroxide solution and 1-naphthylamine (1-NA) is recovered by steam distillation. The distillate is injected into the HPIC and analyzed for 1-NA using an electrochemical detector fitted with a glassy carbon electrode.

#### Method

#### 1. Reagents

Water: Distilled and HFL2-grade

NaOH: 50% solution

Titarous trichlcride: 20% solution, technical

Mossy zinc: Technical

Silicone Surfactant: Antifcam B\*, J.T. Baker Chemical Co.

(or ≅uivalent)

Phosphoric Acid: 85% solution, reagent-grade

1-naphthylamine, 99.4%, AC-917-58

ALANAP, 89.4%, AC-921-54

#### 2. Apparatus

One-liter, flat-bottomed boiling flasks

Distillation heris with water corriensers suitable for distillation of water under an atmospheric pressure and receivers of 250 mL capacity, such as graduated Erlermeyer flasks.

Magnetic stirrer, hot-plams

Magnetic stir bars

High pressure liquid chromatograph with an electrochemical detector

#### Method (Continued)

### HPIC Operating Conditions

Pump:

Model M-6000 Solvent Delivery System, Waters Associates

(or equivalent)

Injector: Model 7120, Rheodyne, 100 µL loop (or equivalent)

Pecorder: Strip chart type, Linear Instruments

Chart speed 20 cm/hour

Detector: Electrochemical detector model EC/230, IRM (or

equivalent)

Electrode: Glassy carbon working electrode set at +850 mV vs.

Ag/Ag Cl reference electrode and 5 to 10 µA (depending

on the sensitivity required).

Column:

PRP-1, CDS, Semi-preparative, 10  $\mu$ m, 305  $\times$  7  $\mu$ m,

Hamilton, protected by a 0.45 µm inlet filter.

Mcbile Phase: 30/70% (v/v) acetcnitrile/0.15 M phosphoric acid

filtered through a 0.45 µm membrare filter and

degassed.

Injection Volume: 100 uL

Flcw Rate:

1.5 mL/mirute

1-MA Retention Time: 10 minutes

### 4. Sample Preparation

Soil samples are removed from the freezer and allowed to thaw for one to two hours. The whole sample is placed into an aluminum foil par and allowed to dry at ambient temperature overnight. To insure uniformity, the whole dried sample is ground to a fine powder and thoroughly mixed.



#### Method (Continued)

### 5. Hydrolysis and Distillation

A sample of the dry homogenized soil, approximately 50 g, is accurately weighed on a calibrated balance into a 1-L, flat-bottomed flask containing a magnetic stir bar. Distilled water, 100 mL, is added and the flask is placed on the magnetic stirrer, hot-plate. While the mixture is stirring, the following are added:

Massy zunc: /

about 1 q

Antifcam B\*:

0.5 to 1.0 mL

50% NaCH Solution:

250 mL

20% Titancus Trichloride: 5 mL

The boiling flask is fitted with the distillation apparatus and the distillate is collected in a 250 mL Erlermeyer flask until 200 mL is distilled. It takes about two hours to complete the distillation.

#### 6. Determination

The analyte, 100 uL, is injected into the MPIC. The 1-NA peak height is measured and the concentration, ug/mL 1-NA, is obtained from a standard calibration plot, or preferably, using a linear regression analysis.

Calculation:

ppm ALANAP in the sample =  $\frac{k \times C}{W \times R}$ 

where.

 $k = 437.6 = \frac{\text{analyte volume, mL}}{\text{mol. wt., l-NA}} \times \text{mol. wt., ALANAP}$ 

C = 1-NA,  $\mu\sigma/mL$ 

W = Sample weight, g

R = 1-NA recovery = parts recovered

#### Method (Continued)

## 7. Standard Solution for a Calibration Plot


A stock solution is prepared by accurately weighing on a calibrated analytical balance about 50 to 60 mg 1-NA standard into a 100 mL volumetric flask. The flask is filled with acetonitrile to the mark and sonicated briefly. A series of dilutions with IC-water are made to obtain standard solutions of approximately 0.15, 0.10, 0.05 and 0.025 µg/mL concentration. One hundred µL of each solution is injected into the HPIC to obtain a calibration plot. 1-NA standard solutions are stable for at least a week when kept in the refrigerator.

# 8. Standard Solution for Recovery Determinations

A stock solution is prepared by accurately weighing on a calibrated analytical balance approximately 800 mg of ALANAP standard into a 100 mL volumetric flask. The flask is filled with water to the mark and sonicated briefly. The stock solution is sarial diluted with IC-water 1/10 twice to give a solution of about 8 µg/mL ALANAP.

# 9. Linearity of Response and Limit of Detection

A linear relationship between the peak height and 1-NA concentration was observed over the range from about 0.15 to 0.025  $\mu$ g/mL, which is the expected range of concentrations in the analyte from the soil analysis. The limit of detection (signal to background ratio of 3) was found to be 0.1 ppm ALANAP in a sample.



### 10. Recovery Determinations

- a. 1-Naphthylamine Recovery from Fortified Reagent Blank
  The reagent blank was fortified with 1, 3, 7 and 10 mL of a
  standard solution containing 8.23 µg/mL ALANAP and the blank
  was treated according to the method. The recoveries ranged from
  99% to 102%, except for the sample fortified at the lowest level
  (7 µg) which gave a 119% recovery (Table I).
- b. 1-Naphthylamine Recovery from Fortified Control Soil
  Samples of the untreated Ockley silt loam soil, approximately 50 g,
  were fortified with 1, 4 and 7 mL of a standard solution containing
  approximately 8.23 µg/mL ALANAP. The samples were then treated
  according to the method. The recoveries ranged from 97% to 105%,
  except for the sample fortified at the lowest level (0.15 ppm)
  which gave a 113% recovery (Table II).

Typical chromatograms of the analyte from a control and a fortified Ockley silt loam soil are shown in the Figure.



Table I. 1-Naphthylamine Recovery from the Reagent Blank Fortified with ALANA

* Recovery	AL ANA P		118.9	99.1	101.1	101.1	
ALANAP Paprivalent	F-1		رن. 1 در	21.03 52.03	74.40		
1-NA Recovered	(2)	0 0	20.0	0.12	0.17		
1-N	E E	15	45	102	150		
AI ANAP Added	119 (1)	7.36	22.07	51.52	73.58	٠	
AIANAP 8.23 µg/mī,	mľ	1	m <sub>.</sub>	7	. 10	ē	
( ) Sample ID	AC-917-54	•	2	m	<b>ਧਾਂ</b>		

1. ALANAP, µg/mL x mL added x standard purity (0.894)

Values obtained by linear regression analysis of the 1-NA standards.

1-NA Std. Plo ug/mL		0.177
3. $\mu g/mL$ (1-NA) x 200 mL, where 200 mL is the	analyte volume and 0.457 is the correction	ror mol. wt.

(AC-917-54) peak ht.	158 104 47 23	
Plot		
ug/mī	0.177 0.117 0.059 0.030	r = 0.9997

Table II. 1-Naphthylamine Recovery from Control Soil Fortified with ALANAP

	± (S.D.)				<b>.</b>		<del>.</del>	. 1.3	<u>!</u>	-	<del>-</del> ~.	4.6			
er.y	ıΣ		113.3				104.5				8.96				
*Recovery			. ייי	113.3	113.3		103 6		105.4		102.1	. 70	1.66	94.1	
ALANAP Denivalent	urid ,	(3)	0.17		0.17	-	0.58		0.59		0.99	96-0	•	96*0	
1-NA Recovered	שש ייני אל/יות איני אל/יות	15)	14 0.02		14 0.02		54 0.07	50 0		0100	0.12	98 0.11		100 0.11	
ALANAP Added 8.23 µg/mL	ppm II. (1)		1 0.15				T 0.56	4 0.56		7 0.97		7 1.02 9	7		
Sample ID size	б	1. (AC-917-53) En 1	1.00 (cc./10 ~ )	2. (AC-917-53) 50.2		3. (AC-917-61) 52.4		4. (AC-917-61) 52.2		5. (AC-917-61) 52.9		0.(//C-91/-61) 50.4	7. (AC-917-53) 50.4		

1. ALANAP, µg/mL x mL added x 0.894/sample wt, g, where 0.894 is the standard purity
2. Values obtained by linear regression analysis of the 1-NA standards.


ug/mL (1-NA) x 200 mL, (0.457) (sample wt)


where 200 mL is the analyte volume and 0.457 is the correction

for molecular weight.

1-NA Std. Plot (AC-917-54) 158 104 47 23 0.117 0.177

r = 0.9997





## Figure . Typical chromatograms

- The analyte from the untreated Ockley silt loam soil.
- 3. The analyte from the same soil fortified at 1 ppm with ALANAP.

PROPERTY OF

PANY