

## Cover Sheet for

# ENVIRONMENTAL CHEMISTRY METHOD

***Pesticide Name:*** Metolachlor

***MRID #:*** 413098-05

***Matrix:*** Soil

***Analysis:*** GC/MS

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at [flynt.elizabeth@epa.gov](mailto:flynt.elizabeth@epa.gov).

STUDY TITLE

Summary of Gas Chromatography  
Mass Spectrometry Confirmations  
of Metolachlor Field Dissipation Study,  
California- Site C  
Agriseach Study No. 1265C

STUDY REQUIREMENT

Field Dissipation Study For Metolachlor

AUTHOR

James E. Whetzel

ANALYTICAL PHASE COMPLETED

November, 1989

ANALYTICAL REPORT ISSUED

November, 1989

CONDUCTED BY

EMS Laboratories, Inc.  
4132 Pompano Rd.  
Charlotte, N. C. 28216

EMS PROTOCOL

04-1089

-1-

AGRISEARCH 1265-C - 264

EMS

## TABLE OF CONTENTS

|                                              | <u>Page Number</u> |
|----------------------------------------------|--------------------|
| TITLE PAGE                                   | 1                  |
| TABLE OF CONTENTS                            | 2                  |
| GOOD LABORATORY PRACTICE STATEMENT           | 3                  |
| QUALITY ASSURANCE STATEMENT                  | 4                  |
| STUDY PERSONNEL                              | 5                  |
| <br>                                         |                    |
| I. Introduction                              | 6                  |
| II. Experimental                             | 6                  |
| A.) Materials                                | 6                  |
| B.) Shipping, Receipt and Storage Of Samples | 6                  |
| C.) Analytical Method                        | 7                  |
| D.) Quality Control and Archiving of Data    | 8                  |
| III. Results and Discussion                  | 9                  |
| IV. Figures                                  | 11                 |
| V. Tables                                    | 17                 |

(Total Pages: 28)

-2-

AGRISEARCH 1265-C - 265

EMS LABORATORIES, INC. 4132 POMEROY ROAD, CHARLOTTE, NC 28216 (704) 393-1853

TOTAL NUMBER OF PAGES IS 294

PAGE



GOOD LABORATORY PRACTICE STATEMENT

To the best of my knowledge, the analytical data reported within was generated by EMS Laboratories, Inc. in compliance with generally accepted EPA Good Laboratory Practice Standards, 40 CFR Part 160, July 1, 1988.

Study Director: James E. Whetzel  
James E. Whetzel  
EMS Laboratories, Inc.

Date: 11/2/89

Study Initiation Date: October 1989  
Study Completion Date: November 1989

-3-

AGRISEARCH 1265-C - 266

EMS

TOTAL NUMBER OF PAGES IS 294

PAGE

266



QUALITY ASSURANCE STATEMENT

The analytical portion of this study was inspected by the EMS Quality Assurance Unit for compliance with GLP standards and EMS Standard Operating Procedures as follows:

Inspection Dates

10-18-89

11-10-89

Dates Submitted  
To Study Director And  
To The Management

10-18-89

11-10-89

Helmuth Janssen

Helmuth Janssen  
Quality Assurance Officer

11-10-89

Date

-4-

AGRISEARCH 1265-C - 267

EMS

TOTAL NUMBER OF PAGES IS 294

PAGE

267

STUDY PERSONNEL

|                      |                                                 |
|----------------------|-------------------------------------------------|
| James Whetzel        | - Study Director                                |
| Steven H. Guptill    | - Senior Scientist                              |
| Helmuth M.B. Janssen | - Quality Assurance Officer                     |
| Lynn Butler          | - Document Control Officer/<br>Sample Custodian |
| Michele E. Sakwa     | - Laboratory Director                           |

AGRISEARCH 1265-C - 268

-5-

EMS

## I. Introduction

The purpose of this report is to summarize the results of gas chromatography-mass spectrometry (GC/MS) confirmation of selected soil samples from metolachlor field dissipation study, California-Site C. Samples for GC/MS confirmation were provided by Agriseach, Inc., Frederick, Maryland after gas chromatography analyses using a flame ionization detector (GC/FID). Representative sample extracts showing positive detections of metolachlor or any of the metabolites of metolachlor: CGA-50720, CGA-40919, CGA-40172 or CGA-51202, were analyzed by GC/MS using single ion monitoring (SIM) approach. The screening level for GC/MS confirmation was 0.1ng/ul. This report describes the procedure used and summary of data for GC/MS confirmation of samples requested by Ciba-Geigy Corporation, Greensboro, North Carolina. Questions regarding the GC/MS confirmations mentioned in this report can be directed to James E. Whetzel, EMS Laboratories, phone number 704/393-1853.

## II. Experimental

### A). Materials

- 1) Autosampler vials, 100 ul
- 2) Hexane (pesticide grade)
- 3) Gas Chromatographic Standards which were prepared from an Agriseach, Inc. 25 ug/ml stock solution of metolachlor, metabolites: CGA-40172, CGA-40919, and methylated metabolites: CGA-51202, CGA-50720.

### B). Shipping

- 1) All extracts were shipped (in dry ice) over night in coolers from Agriseach, Inc. to EMS Laboratories.
- 2) Extracts were accompanied by a sample transmittal form (chain-of-custody) that had been filled out by the person relinquishing the extracts. A signed copy of this form was returned to Agriseach, Inc.
- 3) Extracts were logged into the LIMS (Laboratory Information Management System) and given unique identification numbers. Extracts were stored under refrigeration until removal for analysis. All retainer extracts were kept under refrigeration for 90 days after issuance of final report.

-6-

AGRISEARCH 1265-C - 169

EMS

c). Analytical Method

Analyses were performed using a gas chromatograph/mass spectrometer (GC/MS) in single ion monitoring (SIM) mode to specifically detect and confirm the presence of metolachlor and metabolites: CGA-50720, CGA-40172, CGA-40919 or CGA-51202.

Ions used for SIM analysis were chosen from mass spectra produced by scanning from 35 to 500 AMU each compound of interest. For each compound, the five most abundant ions greater than 100 AMU were chosen as the group to monitor. The most abundant ion of the group was chosen to be the quantitation ion. (Ions less than 100 AMU were disqualified to minimize interferences). Due to the close elution pattern and low sensitivity of CGA-40172 and CGA-51202, one group of four ions was selected to monitor both metabolites.

1) Gas Chromatographic Conditions

Instrument: Hewlett Packard 5890 Gas Chromatograph  
Column: Capillary J & W DB-5, 0.25mm ID, 30m length, 0.25um film  
Injector Temperature: 285°C  
Oven Temperature: 180°C  
Analysis Time: 19 minutes  
Carrier Flow: Helium at 1-2ml/min. at 10 psi

2) Detector Conditions

Detector: Hewlett Packard 5970 Mass Selective Detector in SIM mode.  
Detector Temp: 285°C  
Multiplier Voltage: 2200 volts  
Number of groups: 4

| Group 1 start/stop scan: | 5.00/7.50 |              |
|--------------------------|-----------|--------------|
|                          | M/Z       | Dwell (MSEC) |
| Quant ion CGA-50720      | 134.1     | 100          |
|                          | 146.0     | 100          |
|                          | 162.1     | 100          |
|                          | 163.1     | 100          |
|                          | 221.1     | 100          |



Group 2 start/stop scan: 7.50/10.00

| M/Z                       | Dwell (MSEC) |
|---------------------------|--------------|
| 146.0                     | 100          |
| 160.1                     | 100          |
| Quant ion CGA-40919 161.1 | 100          |
| 188.1                     | 100          |
| 233.1                     | 100          |

Group 3 start/stop scan: 10.00/13.00

| M/Z                       | Dwell (MSEC) |
|---------------------------|--------------|
| 146.1                     | 125          |
| Quant ion CGA-40172 162.1 | 125          |
| 220.1                     | 125          |
| Quant ion CGA-51202 248.1 | 125          |

Group 4 start/stop scan: 13.00/19.00

| M/Z                       | Dwell (MSEC) |
|---------------------------|--------------|
| 146.1                     | 100          |
| Quant ion Metachlor 162.1 | 100          |
| 211.1                     | 100          |
| 238.1                     | 100          |
| 240.1                     | 100          |

3. Instrument Calibration

Prior to the analysis of extracts, known quantities of metolachlor and the metabolites of interest were injected to construct an initial calibration curve for each compound. Calibration curves for concentrations between 0.10ng and 2.00ng were prepared by plotting the concentration of the standard versus the peak area of quantitation ion. Second degree equations of the curves were determined and used in all low level calculations. Standards containing metolachlor and each of the metabolites of interest were analyzed throughout the run to check the calibration.

4. Detection Limits

The method detection limit (MDL) was determined by the method listed in 40 CFR Part 136, App. B. The detection limit reported for the analysis is 0.1ng/ul.

-8-

AGRISEARCH 1265-C - 271

EMS

5. Calculation Method

The ng/ul values found for each extract were determined by the external standard method. The area of the quantitation ion of each compound found was substituted into the appropriate second degree equation as determined by the initial calibration, and ng/ul value was generated.

D. Quality Control and Archiving of Data

- 1) All analyses were performed in compliance with Good Laboratory Practices (40CFR, Part 160). All raw data was audited and approved by the Quality Assurance Officer.
- 2) All original raw data has been stored in the archives at EMS Laboratories, 4132 Pompano Rd., Charlotte, NC 28216. Copies of raw data have been provided to Ciba-Geigy Corporation.

III. Results and Discussion

A summary of the analytical results of the GC/MS confirmation and blind sample analysis is provided in Table 1 and 2 respectively. These tables provide information regarding sample code by Agriseach and EMS Laboratories, Agriseach project number, dates samples were received and analyzed and GC/MS results in units of ng/ul. The individual reports for the analyses of samples were provided to Ciba-Geigy Corporation.

The data in Table 1 and 2 show the use of the term "BDL". This expression means that the quantity is less than the detection limit.

The following Figures have been included:

- Figures 1-5: Typical calibration curves used for quantitation.
- Figures 6: Typical standard chromatogram.
- Figures 7-11: Reference spectra for Metolachlor and the metabolites of interest.
- Figure 12: Example chromatogram of sample with confirmed presence of CGA-51202.
- Figure 13: Example chromatogram of sample with all compounds of interest below the detection limit.

- Figure 14: Example of sample spectrum with confirmed CGA-51202.  
Figure 15: Ion chromatograms of sample with confirmed presence of CGA-51202.  
Figure 16: Ion chromatograms of sample showing absence of CGA-40172.

Positive confirmation was noted only if the following criteria was met:


- 1) Retention time of sample peak had to be comparable to the reference standard.
- 2) All ions of reference standard had to be present in the sample spectrum.
- 3) Ion chromatograms of each ion of interest had to be comparable to those of the reference standard.
- 4) After appropriate background subtraction, intensity ratios had to be reasonably comparable to those of the reference standard.

-10-

AGRISEARCH 1265-C - 273

EVS

FIGURE 1



00002

Compound # 1 Calib File: CALMB111EX

Compound: CGR-50720

Istd:

File: >B4506 >B4505 >B4504 >B4503

Conc: .10 .50 1.00 2.00

Rf: 38440. 125694 106463 103782

Average of 4 Rf's: 93595 (40.65 % Rsd) R= .0000000 Ry: .0000000

1st Degree Equation:  $y = -.005814 + .0000095(x)$

1st Degree Corr Coef: .9964559

2nd Degree Equation:  $y = .0538667 + .0000074(x) + .0000000(x^2)$

2nd Degree Corr Coef: .9985181

In the above equations:

$$y = \frac{\text{Conc Std}}{\text{Conc Istd}} \quad x = \frac{\text{Area Std}}{\text{Area Istd}}$$


Istd Conc for all calibration points is: 1.00

AGRISEARCH 1265-C - 274

FIGURE 2

Calib File: CALMB111EX Comp # 2  
 Calib Date: 891012 22133  
 Comp: CGA-40919

00003



Compound # 2 Calib File: CALMB111EX

Compound: CGA-40919  
 Istd:

Area: 34506 34505 34504 34503  
 Conc: 1.0 1.5 1.00 2.00  
 Area: 45600 35210 37450 11976

Average of 4 Area: 50059 120.68 2 Std 0.000000 3450000  
 1st Degree Equation:  $y = -1.23131 + .0000745(x)$   
 1st Degree Corr Coef: .987623  
 2nd Degree Equation:  $y = .0136653 + .0000171(x) + .0000000(x^2)$   
 2nd Degree Corr Coef: .996977

In the above equations:

Conc Std Area Std  
 $y = \frac{\text{Conc Std}}{\text{Conc Istd}} \times \frac{\text{Area Std}}{\text{Area Istd}}$

Istd Conc for all calibration points is: 1.00


AGRISEARCH 1265-C - 275

-12-

BEST AVAILABLE COPY

FIGURE 3

Calib File: CALMB111EX Comp = 3  
 Calib Date: 891012 22:33  
 Comp: CGA-40172



00004

Compound # 3 Calib File: CALMB111EX

Compound: CGA-40172

Istd:

File: >B4506 >B4505 >B4504 >B4503  
 Conc: .10 .50 1.00 2.00  
 Rf: 38140. 59622. 50239. 45110.

Average of 4 Rf: 48278. (18.73 % Rsd) Rx: .0001000 Ry: .0000010  
 1st Degree Equation:  $y = -.073799 + .0000224(x)$   
 1st Degree Corr Coef: .9945191  
 2nd Degree Equation:  $y = .0316397 + .0000146(x) + .0000007(x^2)$   
 2nd Degree Corr Coef: .9993496

In the above equations:

$$x = \frac{\text{Conc Std} \cdot \text{Area Istd}}{\text{Conc Istd} \cdot \text{Area Std}}$$


Istd Conc for all calibration points is: 1.00

AGRISEARCH 1265-C - 276

-13-

BEST AVAILABLE COPY

00025



Compound: CGA-51202

```

Coverage of 4 Rfs: 50335. (19.00 % Rsd) R= .9600000 Rv: .9900000
1st Degree Equation: y = -.105702 + .0000275(x)
1st Degree Corr Coef: .7213462
2nd Degree Equation: y = -.0042794 + .0000153(x) + .0000000(x^2)
2nd Degree Corr Coef: .9958446

```

$$y = \frac{\text{Conc Std}}{\text{Conc Istd}} \quad x = \frac{\text{Area Std}}{\text{Area Istd}}$$


AGRISEARCH 1265-C - 277

-14-

FIGURE 5

Calib File: CALMB11:EX Comp = 5  
 Calib Dates: 091012 22132  
 Comp: METOLACHLOR

00006



Compound # 5 Calib File: CALMB11:EX

Compound: METOLACHLOR  
 Istd:

| Area | B4506  | B4505  | B4504  | B4503  |
|------|--------|--------|--------|--------|
| Conc | .10    | .50    | 1.00   | 2.00   |
| Area | 107800 | 167832 | 123802 | 126467 |

Average of 4 Rfs: 132475 (19 59 2 Red) Rx: 1.000000 Ry: 1.000000  
 1st Degree Equation:  $y = -.051207 + .0000081(x)$   
 1st Degree Corr Coef: .9938732  
 2nd Degree Equation:  $y = -.000501 + .0000068(x) + .0000000(x^2)$   
 2nd Degree Corr Coef: .9948670

In the above equations:

|                                                |                                                |
|------------------------------------------------|------------------------------------------------|
| Conc Std                                       | Area Std                                       |
| $y = \frac{\text{Conc Std}}{\text{Conc Istd}}$ | $x = \frac{\text{Area Std}}{\text{Area Istd}}$ |

Istd Conc for all calibration points is: 1.00

AGRISEARCH 1265-C - 278


BEST AVAILABLE COPY



FIGURE 6

Typical chromatogram of standard solution of all  
five compounds at a concentration of 2.0ng/ul

00008



Data File: >B4503.D8  
Name: 2.0NG METB  
Misc: CALIBRATION

Quant Output File: >B4503.DT

SP # 1

Id File: MTBID.D8  
Title: METOLACHLOR-METABOLITES ID FILE  
Last Calibration: 891012 16:12

Operator ID: JIM  
Quant Time: 891012 22:13  
Injected at: 891012 19:01


AGRISEARCH 1265-C - 279

-16-

FIGURE 7

Reference spectra for CGA-50720

00010



Data File: B4503:08

Quant Output File: B4503:QT

Name: 2.0NG METB

Misc: CALIBRATION

Quant Time: 891012 22:13

Quant ID File: MTBID:08

Injected at: 891012 17:01

Last Calibration: 891012 16:12

Compound No: 1

Compound Name: CGA-50720

Scan Number: 146

Retention Time: 6.22 min.

Quant Ion: 162.1

Area: 207563

Concentration: 2.23 NG


q-value: 100

AGRISEARCH 1265-C - 280

-17-

FIGURE 8

Reference spectra for GCA-40919



00012

Data File: >B4503.D8  
Name: 2.0NG METB  
Misc: CALIBRATION  
Quant Time: 891012 22:13  
Injected at: 891012 19:01

Quant Output File: >B4503.QT

STL# 3


Quant ID File: METID.D8  
Last Calibration: 891012 16:12

Compound No: 2  
Compound Name: GCA-40919  
Scan Number: 457  
Retention Time: 9.87 min  
Quant Ion: 161.1  
Area: 83953M  
Concentration: 2.22 NG  
I-value: 100

AGRISEARCH 1265-C - 281

FIGURE 9

Reference spectra for CGA-40172



00014

Data File: B4503.D8  
Name: 2.0NG METB  
Misc: CALIBRATION  
Quant Time: 891012 22:13  
Injected at: 891012 19:01

Quant Output File: B4503.D8


Quant ID File: METBID.D8  
Last Calibration: 891012 14:12

Compound No: 3  
Compound Name: CGA-40172  
Scan Number: 755  
Retention Time: 11.67 min.  
Quant Ion: 162.1  
Area: 90221  
Concentration: 2.34 NG  
q-value: 100

AGRISEARCH 1265-C - 282

FIGURE 10

Reference spectra for CGA-51202



Data File: 84503:08  
Name: 2.0NG METB  
Misc: CALIBRATION  
Quant Time: 891012 22:13  
Injected at: 891012 19:01

Quant Output File: 84503:07


Quant ID File: METID:05  
Last Calibration: 891012 15:12

Compound No: 4  
Compound Name: CGA-51202  
Scan Number: 935  
Retention Time: 12.10 min.  
Quant Ion: 248.1  
Area: 87727  
Concentration: 2.23 NG  
S-value: 100

AGRISEARCH 1265-C - 283

FIGURE 11  
Reference spectra for Metolachlor

00018



Data File: >B4503:08  
Name: 2.0NG METB  
Misc: CALIBRATION  
Quant Time: 891012 22:13  
Injected at: 891012 19:01

Quant Output File: >B4503:08

Quant ID File: MTBID:08  
Last Calibration: 891012 16:12


Compound No: 5  
Compound Name: METOLACHLOR  
Scan Number: 1015  
Retention Time: 13.63 min.  
Quant Ion: 162.1  
Area: 252934  
Concentration: 2.44 NG  
Quality: 100

AGRISEARCH 1265-C - 284

-21-

FIGURE 12

Example chromatogram of sample extract (Agrisearch sample 2720-9R)  
with a confirmed detection of CGA-51202 above the detection limit.



Data File: 847001403  
Name: 817560  
File: AGRISARCH 2720-9R

Quant Output File: 847001403


File: 847001403  
Last Calibration: 971102

Operator: J. M. M. M. M.  
Quant Time: 971102  
Injected: 971102

AGRISARCH 1265-C - 285

FIGURE 13

Example chromatogram of a sample extract (Agriseach sample No. 7011-8R)  
which contained no compounds of interest above the detection limit.



00152

Data File: 84516:09

Quant Output File: 84516:12

Name: 817545

File: AGRISEACH 7011-8R

817545

ID File: MISID:08

Title: METOLACHLOR-METABOLITES ID FILE

Last Calibration: 891012 22:37

Operator ID: JIM

Quant Time: 891013 01:11


Injected at: 891013 00:47

AGRISEARCH 1265-C - 286




Example of sample spectrum (Agrisearch sample No. 2720-9R)  
which confirms the presence of CGA-51202.


REFERENCE STANDARD SPECTRUM



SAMPLE SPECTRUM (BACKGROUND SUBTRACTED)



SAMPLE SPECTRUM (UNALTERED)



Data File: >B4700:0A  
Name: B17560  
Misc: AGRISEARCH 2720-9R  
Quant Time: 09:10:07 10:10  
Injection: 09:10:07 11:10

Quant Output File: >B4700:0T  
BT0 1  
Quant ID File: M7610:02  
Last Calibration: B71102 09:07


Compound No: 4  
Compound Name: CGA-51202  
Scan Number: 044  
Retention Time: 12.33 min  
Quant Limit: 240.0  
Area: 20450M  
Concentration: 31 NG  
q-value: 100

**BEST AVAILABLE COPY**

AGRISEARCH 1265-C - 287

FIGURE 15

Ion chromatogram of Agriseach sample No. 2720-9R showing the presence of all four ions monitored for CGA-51202 at the CGA-51202 retention time of 12.33 minutes.



Date File: >B4700:DA  
 Name: 817560  
 MS# : AGRISEARCH 2720-9R  
 Quant Time: 891102 10:24  
 Injected at: 891102 10:00

Quant Output File: >B4700:QT

Quant ID File: MTB10:02  
 Last Calibration: 891102 09:03  
 BTL# 1

Compound No: 4  
 Compound Name: CGA-51202  
 Quant Ion: 248.1

| Hit | R.T.  | Scan# | Area   | Conc    | Units | g   |
|-----|-------|-------|--------|---------|-------|-----|
| 1   | 12.33 | 863   | 20450M | .31, NG |       | 100 |


BEST AVAILABLE COPY

AGRISEARCH 1265-C - 288

-25-

FIGURE 16

Ion chromatograms of Agriseach sample No. 7011.8R showing the absence of the four ions monitored for CGA-40172.



891017

Data File: >B4516:08  
Name: B17565  
Misc: AGRISEARCH 7011-8R  
Quant Time: 891013 01:11  
Injected at: 891013 00:47

Quant Output File: >B4516:08

BT.#15  
Quant ID File: MISID:08  
Last Calibration: 891012 22:17

Compound No: 3  
Compound Name: CGA-40172  
Quant Ion: 162.1

| Hit | R.T.  | Scan# | Area | Conc | Units | q   |
|-----|-------|-------|------|------|-------|-----|
| 1   | 11.59 | 775   | 89   | .03  | NG    | 100 |

BEST AVAILABLE COPY

AGRISEARCH 1265-C - 289