

2.0 INTRODUCTION

Independent laboratory validation of enforcement methods is required by the EPA (Draft Guidelines for the Independent Laboratory Validation of Environmental Chemistry Methods). The study was designed to demonstrate the utility, ruggedness, and efficiency of the subject methods and to identify any inherent weakness in the subject methods as written. The analytical methods DuPont-2549, "Analytical Method for the Determination of Hexazinone and Its Metabolites in Soil by GC/NPD Analysis", and DuPont-2292 Draft 2 (April 23, 1999), "Enforcement Analytical Method for the Determination of Hexazinone and Metabolites of Interest in Soil and Water Using Electrospray-LC/MS/MS", both applicable for the quantitation of hexazinone, IN-T3937, IN-A3928, IN-T3935, IN-G3453, IN-JS472, and IN-G3170 in soil were validated. Soil is a matrix that may be subject to

enforcement testing for hexazinone, IN-T3937, IN-A3928, IN-T3935, IN-G3453, IN-JS472, and IN-G3170.

Samples were acidified and concentrated by solid phase extraction (SPE) using disposable graphitized carbon cartridges. The analytes were eluted from the SPE column with acid-acetone solution. The eluate was evaporated to dryness and brought to the water phase for a second column cleanup using a C18 SPE column. The analytes were then eluted from the column with acid-methanol. The eluate was evaporated to about 1 mL, then brought to a volume of 2 mL with methanol. The final extract was split at this point and 1 mL was evaporated, redissolved in 1.0 mL HPLC water, sonicated, then filtered through a 0.2 micron Acro Disc filter. This extract was then analyzed by ESI-LC/MS/MS using positive mode multiple reaction monitoring (MRM) for hexazinone, IN-T3937, IN-A3928, IN-T3935, IN-G3453, IN-JS472, and IN-G3170. Quantitation was based on the integration of a single MRM transition response. The other 1 mL was evaporated and reconstituted into acetone, ethyl acetate, and toluene (total of 1 mL final volume) and analyzed by GC/NPD. Samples which were analyzed with only one instrument were not split and were brought to a final volume of 2 mL in the proper solvent(s).

3.0 MATERIALS AND METHODS

3.1 Test Substances

DPX-A3674

<u>DuPont Code Number:</u>	DPX-A3674
<u>Common Name:</u>	Hexazinone
<u>Reference Number:</u>	not available
<u>Chemical Name:</u>	S-triazine-2,4/1H,3H/-dione 3-cyclohexyl-1-methyl-6-dimethyl-amino
<u>CAS Registry No:</u>	51235-04-2
<u>Lot Number:</u>	Dash # 233
<u>Purity:</u>	99.9%
<u>Centre I.D. No.:</u>	99-08-59
<u>Storage Conditions:</u>	Refrigerated

IN-T3937

<u>DuPont Code Number:</u>	IN-T3937
<u>Common Name:</u>	Metabolite A of Hexazinone
<u>Reference Number:</u>	AG0084-154
<u>Chemical Name:</u>	S-triazine-2,4/1H,3H/-dione 6-dimethylamino-3-/4-hydroxycyclohexyl/-1-methyl-
<u>Lot Number:</u>	Dash # 3
<u>Purity:</u>	99.0%

Centre I.D. No.: 99-08-55
Storage Conditions: Frozen

IN-A3928

DuPont Code Number: IN-A3928
Common Name: Metabolite B of Hexazinone
Reference Number: E76965-35
Chemical Name: S-triazine-2,4/1H,3H/-dione 3-cyclohexyl-1-methyl-6-methylamino
CAS Number: 56611-54-2
Lot Number: Dash # 4
Purity: 96.3%
Centre I.D. No.: 99-08-56
Storage Conditions: Room Temperature

IN-T3935

DuPont Code Number: IN-T3935
Common Name: Metabolite C of Hexazinone
Reference Number: E5255-9
Chemical Name: S-triazine-2,4/1H-dione 3-/4-hydroxycyclohexyl/-1-methyl-6-methylamino-
CAS Number: 72585-88-7
Lot Number: Dash # 3
Purity: 95.6%
Centre I.D. No.: 99-08-57
Storage Conditions: Frozen

IN-G3453

DuPont Code Number: IN-G3453
Common Name: Metabolite A1 of Hexazinone
Reference Number: E70731-2
Chemical Name: S-triazine-2,4/1H,3H/-dione 3/trans-2hydroxycyclohexyl/-1-methyl-6-dimethylamino
Lot Number: Dash # 2
Purity: 99.4%
Centre I.D. No.: 99-08-58
Storage Conditions: Room Temperature

IN-JS472

DuPont Code Number: IN-JS472
Common Name: Metabolite 1 of Hexazinone
Reference Number: not available
Chemical Name: 1,3,5-triazine-2,4(1H,3H)-dione,6-(dimethylamino)-1-methyl-3-(4-oxycyclohexyl)
Lot Number: Dash # 2
Purity: 95.5%
Centre I.D. No.: 99-08-54
Storage Conditions: Room Temperature

IN-G3170

DuPont Code Number: IN-G3170
Common Name: Metabolite G3170 of Hexazinone
Reference Number: noi available
Chemical Name: S-triazine-2,4/1H,3H/-dione-1-methyl-6-dimethylamino-
Lot Number: Dash # 2
Purity: 92.0%
Centre I.D. No.: 99-08-53
Storage Conditions: Room Temperature

Characterization and certification records for the analytical standards will be archived by:

E.I. du Pont de Nemours and Company
Wilmington, DE 19898, U.S.A.

3.2 Test System

The analytical methods DuPont-2549, "Analytical Method for the Determination of Hexazinone and Its Metabolites in Soil by GC/NPD Analysis", and DuPont-2292 Draft 2 (April 23, 1999), "Enforcement Analytical Method for the Determination of Hexazinone and Metabolites of Interest in Soil and Water Using Electrospray-LC/MS/MS", both applicable for the quantitation of hexazinone, IN-T3937, IN-A3928, IN-T3935, IN-G3453, IN-JS472, and IN-G3170 in soil were validated.

Soil sample field ID 71189 (CAL 992713) was shipped to Centre Analytical Laboratories, Inc. from Morse Laboratories in Sacramento, CA on 3/29/99 and received on dry ice at Centre Analytical Laboratories, Inc. on 3/31/99, logged in, and given a unique identification number. It was then stored in a walk-in freezer at a temperature of $\leq -10^{\circ}\text{C}$ until extraction.

3.3 Equipment

All equipment used in this study followed the guidelines specified in the method. Equivalent changes in apparatus were made in the following instance:

- The method specifies an HP5890 gas chromatograph and HP6890 autosampler, instead a Varian 3500 Gas Chromatograph with a Varian 8200 autosampler was used.

3.4 Reagents

Reagents used during the validation included the following:

Acetone- J.T. Baker, HPLC grade

Envi-Carb SPE column, Supelco

C₁₈ SPE column, Varian

Dimethyldichlorosilane – Supelco

Ethyl Acetate - J.T. Baker, HPLC grade

Glacial Acetic Acid

Hexane- J.T. Baker, HPLC grade

Methanol - J.T. Baker, HPLC grade

Potassium phosphate – J.T. Baker

Sodium Chloride – J.T. Baker

Toluene - J.T. Baker, HPLC grade

All water was Type I distilled, deionized water (CAL)

3.5 Principles of the Analytical Methods

The analytes were extracted from the soil with an aqueous buffer-acetone solution. The resulting extract was passed through solid phase extraction (SPE) using disposable graphitized carbon cartridges. The analytes were eluted from the SPE column with acid-acetone solution. The eluate was evaporated to dryness and brought to the water phase for a second column cleanup using a C₁₈ SPE column. The analytes were then eluted from the column with acid-methanol. The eluate was evaporated to about 1 mL, then brought to a volume of 2 mL with methanol. The final volume was split at this point and 1 mL was evaporated, redissolved in 1.0 mL HPLC water, sonicated, then filtered through a 0.2 micron Acro Disc filter. This extract was then analyzed by ESI-LC/MS/MS using positive mode multiple reaction monitoring (MRM) for hexazinone, IN-T3937, IN-A3928, IN-T3935, IN-G3453, IN-JS472, and IN-G3170. Quantitation was based on the

integration of a single multiple reaction monitoring (MRM) transition response. The other 1 mL was evaporated and reconstituted into acetone, ethyl acetate, and toluene (total of 1 mL final volume) and analyzed by GC/NPD. Samples that were analyzed with only one instrument were not split and were brought to a final volume of 2 mL in the proper solvent(s).

3.6 Modifications, Interpretations and Critical Steps

These methods were run as written, with the modification described in the protocol Appendix. Protocol Deviation number 4 specifies a change made to the re-equilibration solvent from 100% acetonitrile (mobile phase B) to 100% aqueous 0.01M acetic acid (mobile phase A) due to an error in the original draft of the method. This protocol deviation did have an impact on the validity of the study. Without this deviation the instrument would not have been properly equilibrated between sample injections and the validation would not have passed. There was an additional change in the methods entailing a minor equivalent equipment substitution outlined in Section 3.3.

3.7 Instrumentation

DuPont-2550

1. Instrument: Varian 3500
Gas Chromatograph
2. Detector: Nitrogen Phosphorus Detector
3. Column: Rtx-35, 15m x 0.53mm x
0.5 μ m df
4. Integrator: HP Chemstation
5. Instrument
Conditions:
 - a. Gas Flow Rates: Hydrogen, 4 mL/min
Column, 5 mL/min.
Make-up, 15 mL/min.
 - b. Retention Time: ~ 4.8, 6.6, 6.9, 7.5, 7.8, 7.8,
8.1 min.
 - c. Run Time: ~ 9.5 min.
 - d. Injection Volume: 2 μ L
 - e. Injector
Temperature: 290°C

f. Detector
 Temperature: 285°C

g. Column
 Temperature: 140°C hold for 1 min., to 290°C at
 20°C /min, hold for 1 min.

DuPont-2292 Draft 2 (April 23, 1999)

Instrument: LC/MS/MS
 Micromass Quattro LC
 Electrospray Ion source
 Desolvation Temp.: 400°C
 Desolvation N₂ Flow Rate: 850 L/hr
 Source Temp.: 100°C
 Nebuliser N₂ Flow Rate: 100 L/hr

Computer: Digital 266i Personal Workstation

Software: Microsoft Windows NT
 Version 4 Build 1381: Service Pack 3
 Micromass Limited
 MassLynx 3.1 Build 004

HPLC Equipment: Hewlett Packard (HP) Series 1100
 HP Quat Pump
 HP Vacuum Degasser
 HP Autosampler
 HP Column Oven

HPLC Column: Zorbax Rx-C8, 4.6 mm x 25cm
 Column Temperature: 35° C
 Mobile Phase (A): Aqueous 0.01 M acetic acid
 Mobile Phase (B): Acetonitrile

<u>Time</u>	<u>% A</u>	<u>% B</u>
0.0	100	0
3.0	90	10
10.0	50	50
15.0	25	75
15.1	5	95
20.0	5	95
20.1	100	0
30.0	STOP	STOP

Total run time = 30 min.

Flow Rate: 1.0 mL/min (Split Flow 5:1)
 Injected Volume: 50 µL

Ions monitored :

Analyte	Parent ion	Daughter ion	Approximate Retention Time (min.)
hexazinone	253.0	171.0	13.1
IN-G3170	171.0	71.1	7.2
IN-T3935	255.0	156.8	8.7
IN-T3937	269.0	171.0	9.3
IN-JS472	267.0	171.0	9.9
IN-G3453	269.0	171.0	10.3
IN-A3928	239.3	156.9	12.2

A full-scan was performed to verify appropriate m/z responses and abundance ratios (See Appendix 1).

3.8 Calculations

The calculations were performed exactly as described in the analytical methods. For DuPont-2549, the ppb found and percent recoveries were determined from calculations using the peak area responses of each analyte in the sample. These data were plotted versus concentration (ng/mL) of the corresponding standard to obtain linear regression standard calibration curves. Standard curves were prepared each analysis day.

An example of the calculation is presented here using an actual sample of IN-G3170, Spk B, reported in Table 2:

$$1. \mu\text{g/mL analyte} = \frac{2643.627 - (-1395.245)}{78147.327} = 0.0517 \mu\text{g/mL}$$

$$2. \text{ppb analyte} = \frac{0.0517 \mu\text{g/mL} \times 1000 \text{ ng}/\mu\text{g} \times 1 \text{ mL final vol.} \times 2 \text{ GC aliquot fact.}}{5 \text{ g sample}} = 20.7 \text{ ppb}$$

$$3. \% \text{ Recovery} = \frac{20.7 \text{ ppb} - \left(\frac{(0 \text{ ppb control}(1) + 0 \text{ ppb control}(2))}{2} \right)}{20 \text{ ppb added}} = 103\%$$

For DuPont-2292 Draft 2 (April 23, 1999), the ppb found and percent recoveries were determined from calculations using the peak area responses of each analyte in the sample and the average response factor from the one preceding standard and the one following standard.

An example of the calculation is presented here using an actual sample of IN-G3170, Spk B, reported in Table 1:

1. The ppb Found calculated for Run 052799-010:

$$\begin{aligned} &\text{Response Factor of preceding Std.} \\ &21873/25.0 = 875 \text{ area/ng/mL} \end{aligned}$$

$$\begin{aligned} &\text{Response Factor of following Std.} \\ &39359/50.0 = 787 \text{ area/ng/mL} \end{aligned}$$

Average RF = 831

$$\text{Calc. Amt. (ng/mL)} = \frac{9109}{(831)} = 10.96 \text{ ng/mL}$$

$$\text{Amt (ppb)} = \frac{(10.96 \text{ ng/mL})(1 \text{ mL})(4 \text{ HPLC DF})(2 \text{ aliquot fract.})}{(5 \text{ g})}$$

$$= 17.5 \text{ ppb}$$

$$\begin{aligned} \% \text{ recovery} &= (\text{ppb found/ppb added}) \times 100 \\ &= (17.5 \text{ ppb}/20.0 \text{ ppb}) \times 100 = 88\% \end{aligned}$$