Residue Method for the Determination of NOA449280, and Metabolites SYN503780, CSCC163768, CSCD656831, and SCD642512 in Water.

Data Requirement:	EPA Guideline:	850.6100
	OECD Data Point:	IIA 4.5

Reports:

 Analytical Method: MRID 47842016. Braid, S. Richardson, M. and Eberhard Zietz. 2012. NOA449280 - Residue Method GRM030.06A for the Determination of NOA449280 and its metabolites SYN503780, CSCC163768, CSCD656832, CSCD642512 and CSAA806573 in Water. Analytical Method. Report No. GRM030.06A. Task no. T000995-08. Unpublished study prepared by Syngenta Ltd., Jealott's Hill International Research Centre, Bracknell, Berkshire, RG42 6EY, UK. Submitted by Syngenta Crop Protection, LLC, Greensboro, NC.

Independent Laboratory Validation: MRID 47841956. Zeitz, E. 2011. NOA449280 – Validation of an Analytical Method (DraftGRM030.06A) for the Determination of NOA449280 and its Metabolites SYN503780, CSCC163768, CSCD656832, CSCD642512 and CSAA806573 in Water. Method Validation. Report Number: IF-09/01378742. Study Number: IF-09/01378742. Task Number: TK0009384. Report Number: IF-09/01378742. Unpublished study prepared SGS INSTITUT FRESENIUS GmbH Im Maisel 14, D-65232 Taunusstein, Germany. Submitted by Syngenta Crop Protection, LLC, Greensboro, NC.

This amin's carry in submet, and bit 200. A set the inpact for the quantitative three requests in b(1), -4/0110 cancers increaling trans. C.A.S.W. 2020(1):403. A set the form of the case is the constant in the constant. The set product of the case is the increase a models and a set if moved is constant plan and where is to a Sense 5. S.M. casted, standard with restance. The statement is not welled to the variation and planequeries in the case of the definition of the statement is constant to the variation and planequeries in the later state with restance 1. The statement and form to submet in the statement in the statement is the state of the interval and the constant of the variation and plane to the plane of the state of the statement is restared by the statement of the statement of the statement is the state of the function in particular the statement statement of the statement of the statement is the state with the function in particular the statement statement of the statement. The statement of the stat

Table 1. Analytical Mediated Scattering

Statements:			
statements.	These studies were conducted in com	pliance with GLP practices	
Classification:	This study is classified Fully Reliable	e (EPA classification: Acceptable	e).
	However, although this analytical me	thod is valid for the determination	n of
	NOAA449280 and metabolites SYN5	03780, CSCC163768, CSCD656	5831
	SCD642512, and CSAA806573 in wa	ter, the ILV data submitted do not	ot
	qualify as an independent validation,	as they are the same data as in th	e
	initial methods study. Thus, this study		
	requirements, as it does not meet the		
	the data requirements with respect to		
	Additional studies (MRIDs 47841954		
	submitted that do address the parent a		
	bicyclopyrone, but they also lacked su	afficient ILV data.	
PC Code:	018986	umola des y recolection	
Primary		O al mat	_
	Paul Mastradone, Ph.D.	Signature: Tank I man	Leadon 01
Reviewer:	i aui mastrauone, i n.D.		toos
Reviewer:	Chemist (USEPA)	Date: June 30, 2014	tube
Reviewer: Secondary			Digitally signed by Sutton, Cheryl
			Digitally signed by Sutton, Cheryl DN: cn=Sutton, Cheryl, email=sutton, cheryl@epa.gov Date: 2014.11.1912-258.02.05'00'

Executive Summary

This analytical method, GRM030.06A, is designed for the quantitative determination of NOA449280 (parent bicyclopyrone; CAS# 352010-65-5) and metabolites SYN503780, CSCC163768, CSCD656831, and SCD642512 in water. The analytical chemistry method involves acidifying a well mixed water subsample and adding it to a Strata-X SPE cartridge, discarding the eluate. The analytes are sorbed to the column and subsequently quantitatively eluted with methanol, acetonitrile and formic acid, then evaporated to dryness. They are reconstituted in methanol/0.2% formic acid and are appropriately diluted for analysis by LC-MS/MS. Review of this analytical chemistry method indicated that it is acceptable. However, the Independent Laboratory Validation study is a repackaging of the original study data and is not acceptable as an ILV. Thus, the data requirement for two ILVs for the determination of the degradates in water has not been met.

Table 1. Analytical Method Summary

MRID		EDA		i				
Analyte(s) by Pesticide	Environmental Chemistry Method	Independent Laboratory Validation	EPA Lab Review	Matrix	Method Date	Registrant	Analysis	Limit of Quantitation (LOQ)
NOA449280 SYN503780 CSCC163768 CSCD656831 CSCD642512 CSAA806573	47842016	47841956	None	Water	1/31/2012	Syngenta	LC/MS- MS	0.01µg/L

I. **PRINCIPLE OF THE METHOD**

This method starts with a 10 ml sample of thoroughly mixed water. The 10 ml aliquot is then placed in a centrifuge tube and 200 ul of formic acid is added to each sample (pH must be <2, checked with suitable indicator paper) Use one Strata-X SPE cartridge per sample. Place in a vacuum manifold. Add 2ml methanol percolate or draw through column to level of top frit discarding the eluate. Add samples to column allowing flow thru rate of 1-2 ml per minute to level of top frit. analytes are retained on column. At no time should column be allowed to dry. Rinse tube with ultra pure water plus 2ml of 2% formic acid and add to column. Elute as before. Elution of analytes is accomplished adding 5 ml of 0.1% formic acid in methanol that is percolated or drawn through the column under low vacuum collecting the column eluate, Eluate is evaporated to dryness with clean air. Immediately upon dryness methanol (200 ul) is added and the sample ultrsonicated. Add 900 ul of 0.2% formic acid in ultra-pure water and mix well. Immediately prior to analysis dilute an aliquot of the sample twofold and add to auto-sampler vial. Analysis is by LC-MS/MS.

RECOVERY FINDINGS II.

Analyte	Fortification Level (ug/L)	and the second se	Recovery Range (%)	Mean Recovery (%)	Standard Deviation (%)	Relative Standard Deviation (%)
NOA449280	0.01	5	91-109	101	n/a	6.8
Primary transition	0.1	5	92-101	98	n/a	3.4
SYN503780	0.01	5	87-102	95	n/a	6.7
Primary transition	0.1	5	92-99	95	n/a	2.9
CSCC163768	0.01	5	87-98	93	n/a	4.6
Primary transition	0.1	5	96-98	97	n/a	0.7
NCSCD656832	0.01	5	73-107	91	n/a	13.8
Primary transition	0.1	5	75-83	84	n/a	4.7
CSCD642512	0.01	5	86-100	93	n/a	5.6
Primary transition	0.1	5	98-107	101	n/a	3.5
CSAA806573	0.1	5	82-94	90	n/a	5.2

Primary transition	0.01	5	90-96	94	n/a	2.4
NOA449280	0.01	5	89-99	93	n/a	4.3
Confirmatory transition	0.1	5	91-96	93	n/a	2.0
SYN503780	0.01	- 5	102-110	106	n/a	3.2
Confirmatory transition	0.1	5	88-102	97	n/a	6.8
CSCC163768	0.01	5	92-102	99	n/a	4.6
Confirmatory transition	0.1	5	95-101	98	n/a	2.3
CSCD656832	0.01	5	85-92	89	n/a	3.5
Confirmatory transition	0.1	5	78-83	81	n/a	2.4
CSCD642512	0.01	5	80-106	95	n/a	11.4
Confirmatory transition	0.1	5	92-109	97	n/a	6.9
CSCD806573	0.01	5	91-94	92	n/a	1.6
Confirmatory transition	0.1	5	90-96	94	n/a	2.2

Table 3. Initial Validation Method Recoveries for Analytes in Surface Water

Analyte	Fortification Level (ug/L)	Number of Tests	Recovery Range (%)	Mean Recovery (%)	Standard Deviation (%)	Relative Standard Deviation (%)
NOA449280	0.01	5	78-88	85	n/a	4.7
Primary transition	0.1	5	85-97	92	n/a	4.8
SYN503780	0.01	5	89-95	92	n/a	2.7
Primary transition	0.1	5	91-95	94	n/a	2.1
CSCC163768	0.01	5	91-95	94	n/a	1.6
Primary transition	0.1	5	92-98	95	aris marking	2.6
NCSCD656832	0.01	5	75-81	81	n/a	2.6
Primary transition	0.1	5	79-84	84	n/a	4.7
CSCD642512	0.01	5	87-109	93	n/a	10.0
Primary transition	0.1	5	88-94	91	n/a	3.1
CSAA806573 Primary transition	0.1	5	90-95	92	n/a	2.4
	0.01	5	88-95	92	n/a	3.2
NOA449280	0.01	5	92-110	103	n/a	7.5
Confirmatory transition	0.1	5	92-98	95	n/a	2.4
SYN503780	0.01	5	80-86	82	n/a	3.3
Confirmatory transition	0.1	5	86-92	88	n/a	3.0
CSCC163768	0.01	5	100-108	102	n/a	3.5
Confirmatory transition	0.1	5	90-102	96	n/a	5.0
CSCD656832	0.01	5	89-101	96	n/a	4.8
Confirmatory transition	0.1	5	75-93	84	n/a	7.6
CSCD642512	0.01	5	81-92	89	n/a	5.1
Confirmatory transition	0.1	5	83-90	87	n/a	3.3
CSAA806573	0.01	5	85-100	91	n/a	7.2
Confirmatory transition	0.1	5	92-95	93	n/a	1.5

Analyte	Fortification Level (ug/L)	CONTRACTOR OF A	Recovery Range (%)	Mean Recovery (%)	Standard Deviation (%)	Relative Standard Deviation (%)
NOA449280	0.01	5	82-90	86	n/a	3.9
Primary transition	0.1	5	86-93	91	n/a	3.6
SYN503780	0.01	5	94-109	102	n/a	6.1
Primary transition	0.1	5	96-98	97	n/a	0.8
CSCC163768	0.01	5	92-97	95	n/a	2.0
Primary transition	0.1	5	92-94	93	n/a	1.2
NCSCD656832	0.01	5	80-87	85	n/a	3.1
Primary transition	0.1	5	75-84	79	n/a	4.9
CSCD642512	0.01	5	86-101	95	n/a	6.1
Primary transition	0.1	5	86-92	89	n/a	3.2
CSAA806573 Primary transition	0.1	5	90-95	92	n/a	2.4
	0.01	5	88-95	92	n/a	3.2
NOA449280	0.01	5	98-108	103	n/a	5.1
Confirmatory transition	0.1	5	85-97	90	n/a	5.0
SYN503780	0.01	5	101-105	103	n/a	2.1
Confirmatory transition	0.1	5	91-99	94	n/a	3.5
CSCC163768	0.01	5	88-106	99	n/a	6.9
Confirmatory transition	0.1	5	88-97	92	n/a	3.9
CSCD656832	0.01	5	85-94	91	n/a	2.3
Confirmatory transition	0.1	5	72-77	75	n/a	3.2
CSCD642512	0.01	5	90-104	99	n/a	6.0
Confirmatory transition	0.1	5	82-105	94	n/a	9.7
CSAA806573	0.01	5	98-100	98	n/a	1.1
Confirmatory transition	0.1	5	86-93	90		3.0

Table 4. Initial Validation Method Recoveries for Analytes in Ground Water

III. METHOD CHARACTERISTICS

	NOA449280	SYN503780	CSCC163768	CSCD656832	CSCD642512	CSAA806573
Limit of Quantitation (LOQ)*	0.01 µg/L	0.01 µg/L	0.01 µg/L	0.01 µg/L	0.01 µg/L	0.01 µg/L
Limit of Detection (LOD) primary ion	0.0003 µg/L	0.0007 µg/L	0.0005 µg/L	0.001 μg/L	0.0004 µg/L	0.0008 µg/L
Limit of Detection (LOD) confirmatory ion	0.005 μg/L	0.005 μg/L	n/a	n/a	n/a	n/a
Linearity (calibration curve r ² and concentration range)	r ² = 0.99 0.05ng/ mL- #10ngm/L	r ² = 0.99 0.05ng/ mL- #10ngm/	r ² = 0.99 0.05ng/ mL- #10ngm/L			
Repeatable	Yes	Yes	Yes	Yes	Yes	Yes
Reproducible	Yes	Yes	Yes	Yes	Yes	. Yes
Specific	Yes	Yes	Yes	Yes	Yes	Yes

Table 5. Method Characteristics

* The limit of quantitation of the method is defined as the lowest analyte concentration in a sample at which the methodology has been validated and a mean recovery of 70-110% with a relative standard deviation of \leq 20% has been obtained.

IV. METHOD DEFICIENCIES AND REVIEWER'S COMMENTS

The data submitted as independent laboratory validation data are the same data reported for the analytical method validation and, as they were not determined by an independent laboratory, are classified **not acceptable**.

Attachment 1: Chemical Names and Structures:

.

Figure 1 NOA449280

Compound Code Number CAS Number **IUPAC** Name

NOA449280 : 352010-68-5 : 4-hydroxy-3-[2-(2-methoxy-ethoxymethyl)-6-(triflouromethyl)-pyridine-3-carbonyl]-bicyclo[3.2.1]oct-3-en-2-one

Molecular Formula Molecular Weight

: C19H20F3NO5 : 399.4

OH

Figure 2 SYN503780

Compound Code Number : SYN503780 CAS Number 380355-55-5 5 **IUPAC** Name 2-(2-methoxy-ethoxymethyl)-6-triflouromethyl-nicotinic acid : Molecular Formula : C11H12F3NO4 279.2 Molecular Weight 1 HO

Figure 3 CSCC163768

Compound Code Number	4	CSCC163768
CAS Number	1	Not in registry
IUPAC Name		6-(trifluoromethyl)pyridine-2,3-dicarboxylic acid
Molecular Formula	:	C8H4F3NO4
Molecular Weight		235.1

OH

0



CSCD656832

Compound Code Number
CAS Number
IUPAC Name
Molecular Formula
Molecular Weight

- : CSCD656832 : Not in registry
- : 6-(trifluoromethyl)pyridin-3-ol-2-carboxylic acid
- - : C7H4F3NO3 : 207.1





Figure 5 CSCD642512

Compound Code Number	1	CSCD642512
CAS Number	÷	Not in registry
IUPAC Name	:	Not known




Attachment 1: Chemical Names and Structures:

IUPAC Name: CAS Name:

CAS Number:

4-Hydroxy-3-[2-(2-methoxy-ethoxymethyl)-6-trifluoromethyl-pyridine-3carbonyl]-bicyclo[3.2.1]oct-3-en-2-one N/A 352010-65-5



IUPAC Name: CAS Name: CAS Number: 2-(2-Methoxy-ethoxymethyl)-6-trifluoromethyl-nicotinic acid N/A None listed

