
USEPA REGION II DATA VALIDATION SOP FOR SW-846 METHOD 8290 POLYCHLORINATED DIBENZODIOXINS (PCDDs) AND POLYCHLORINATED DIBENZOFURANS (PCDFs) BY HIGH-RESOLUTION GAS CHROMATOGRAPHY/ HIGH-RESOLUTION MASS SPECTROMETRY (HRGC/HRMS)



Page: 1 of 24 Date: December 2010 Revision 1.1

1.0 Introduction

- 1.1 The attached Standard Operating Procedure (SOP) is applicable to polychlorinated dibenzodioxin and polychlorinated dibenzofuran (PCDD/PCDF) data obtained using SW-846 Method 8290, Polychlorinated Dibenzodioxins (PCDDs) and Polychlorinated Dibenzofurans (PCDFs) by High-Resolution Gas Chromatography/High-Resolution Mass Spectrometry (HRGC/HRMS), Revision 0, November 1992. Its scope is to facilitate the data validation process of the data reported by the contracting laboratory and also to ensure that the data is being reviewed in a uniform manner.
- This SOP is based upon the quality control and quality assurance requirements specified in SW-1.2 846 Method 8290, Revision 0, November 1992. This SOP is based also upon additional QA/QC requirements prescribed in the Special Analytical Service (SAS) requests provided to the laboratory.

2.0 Responsibilities

- 2.1 The reviewer must be knowledgeable of the analytical method and its QC Criteria.
- 2.2 The reviewer must complete and/or file the following:
- Data Assessment Checklist The data reviewer must read each item carefully and must check yes 2.2.1 if there is compliance, no if there is non compliance and N/A if the question is not applicable to the data.
- Data Assessment Narrative The data reviewer must present professional judgement and must 2.2.2 express concerns and comments on the validity of the overall data package. The reviewer must explain the reasons for rejecting and/or qualifying the data.
- Rejection Summary Form The reviewer must submit the completed form using a ratio format. The 2.2.3 numerator indicates the number of dioxins/furans data rejected; the denominator indicates the number of dioxins/furans fractions containing rejected compounds.
- Organic Regional Data Assessment Summary The data reviewer is also required to submit the 2.2.4 completed Organic Regional Data Assessment Form.
- 2.2.5 Telephone Record Log - All phone conversations must be initiated by the technical project officer through SMO. If a phone call has been made, the reviewer must transcribe the conversation. After the data review has been completed, the white copy of the telephone log is mailed to the laboratory and the pink copy to SMO. The yellow copy is filed in the appropriate folder. A photocopy of the Telephone Record Log is attached to the Data Assessment Narrative.
- Forwarded Paperwork Upon completion of the review the following are to be forwarded to the 2.2.6 Regional Sample Control Center (RSCC):

USEPA REGION II DATA VALIDATION SOP FOR SW-846 METHOD 8290 POLYCHLORINATED DIBENZODIOXINS (PCDDs) AND POLYCHLORINATED DIBENZOFURANS (PCDFs) BY HIGH-RESOLUTION GAS CHROMATOGRAPHY/ HIGH-RESOLUTION MASS SPECTROMETRY (HRGC/HRMS)



PREPARED BY:		Date:
_	Russell Arnone, Chemist HWSS	
PEER REVIEWED	BY: Muhammad Sheikh, Chemist HW:	Date: SS
CONCURRED BY	: Michael Mercado, Acting Chief HV	Date: VSS
APPROVED BY: _	Robert Runyon, Chief HWSB	Date:
	Annual Review	
REVIEWED BY: _	Name	Date:
REVIEWED BY: _	Nama	Date:
	Name	

Page: 1 of 24 Date: December 2010 Revision 1.1

1.0 Introduction

1.1 The attached Standard Operating Procedure (SOP) is applicable to polychlorinated dibenzodioxin and polychlorinated dibenzofuran (PCDD/PCDF) data obtained using SW-846 Method 8290, Polychlorinated Dibenzodioxins (PCDDs) and Polychlorinated Dibenzofurans (PCDFs) by High-Resolution Gas Chromatography/High-Resolution Mass Spectrometry (HRGC/HRMS), Revision 0, November 1992. Its scope is to facilitate the data validation process of the data reported by the contracting laboratory and also to ensure that the data is being reviewed in a uniform manner.

1.2 This SOP is based upon the quality control and quality assurance requirements specified in SW-846 Method 8290, Revision 0, November 1992. This SOP is based also upon additional QA/QC requirements prescribed in the Special Analytical Service (SAS) requests provided to the laboratory.

2.0 **Responsibilities**

- 2.1 The reviewer must be knowledgeable of the analytical method and its QC Criteria.
- 2.2 The reviewer must complete and/or file the following:
- 2.2.1 Data Assessment Checklist - The data reviewer must read each item carefully and must check yes if there is compliance, no if there is non compliance and N/A if the question is not applicable to the data.
- 2.2.2 Data Assessment Narrative - The data reviewer must present professional judgement and must express concerns and comments on the validity of the overall data package. The reviewer must explain the reasons for rejecting and/or qualifying the data.
- 2.2.3 Rejection Summary Form - The reviewer must submit the completed form using a ratio format. The numerator indicates the number of dioxins/furans data rejected; the denominator indicates the number of dioxins/furans fractions containing rejected compounds.
- 2.2.4 Organic Regional Data Assessment Summary - The data reviewer is also required to submit the completed Organic Regional Data Assessment Form.
- 2.2.5 Telephone Record Log - All phone conversations must be initiated by the technical project officer through SMO. If a phone call has been made, the reviewer must transcribe the conversation. After the data review has been completed, the white copy of the telephone log is mailed to the laboratory and the pink copy to SMO. The yellow copy is filed in the appropriate folder. A photocopy of the Telephone Record Log is attached to the Data Assessment Narrative.
- 2.2.6 Forwarded Paperwork - Upon completion of the review the following are to be forwarded to the Regional Sample Control Center (RSCC):

HW-19 PCDDs/PCDFs using HRGC/HRMS Date: December 2010 Revision 1.1

Page: 2 of 24

- a. data package
- b. completed data assessment checklist and narrative (original)

The reviewer will forward one copy of the completed Data Assessment and one copy of the Organic Regional Data Assessment to the appropriate Regional TPO.

- 2.2.7 Filed Paperwork The following are to be submitted to the Monitoring Management Branch (MMB) files:
 - a. a photocopy of the Data Assessment Narrative
 - b. a photocopy of the Regional Data Assessment Summary
 - c. Telephone record Log (copy)
 - d. Rejection Summary Form
- 2.3 Rejection of Data All values determined to be unacceptable on the Organic Analysis Data Sheet (Form I) must be flagged with an "R". The qualifier R means that due to significant QA/QC problems the analysis is invalid and it provides no information as to whether the compound is present or not. Once the data are flagged with R any further review or consideration is unnecessary. The qualifier "J" is used to indicate that due to QA/QC problems the results are considered to be estimated.
 - The qualifier "NJ" indicates that there is presumptive evidence for the presence of the compound at an estimated value.
 - The data reviewer must explain in the data assessment narrative why the data was qualified. He or she must also indicate all items of contract non-compliance.
 - When 2,3,7,8- substituted TCDD, TCDF, PnCDD and PnCDF data are rejected (flagged "R") or qualified "J" the project officer must be notified promptly. If holding times have not been exceeded reanalysis of the affected samples may be requested.
- All qualifications and corrections to reviewed data must be made in red pencil.

USEPA Region II DV SOP for SW-846 Method 8290 PCDDs/PCDFs using HRGC/HRMS Page: 3 of 24 Date: December 2010 Revision 1.1

	AGE COMPLETENESS AND DELIVERABLES SITE:			
		YES	<u>NO</u>	N/A
I.0 <u>Data</u>	a Completeness and Deliverables			
1.1 1.2 1.3	Are the Traffic Report Forms present for all samples? Is the Narrative or Cover letter present? Are the Case Number and/or SAS numbers contained in the case narrative? Do the Traffic Reports or Lab Case Narrative indicate	[] [] []		-
1.4	problems with sample receipt, sample condition, analytical problems, or other comments affecting the quality of the data?	[_] _	
2.0	ACTION: Use professional judgement to evaluat the quality of the data. Reporting Requirements and Deliverables	e the effect of the noted	proble	ms or
2.1	All deliverables must be clearly labeled with the SMO number number. Missing or illegible or incorrectly labeled items must immediately be contacted and requested to submit the miss.	t be identified. The conti sing or incorrect items.	actor i	must
2.2	The following forms were taken from the CLP SOW, DFLM0 Request. Are these forms present?	1.1 and are specified in t	ne SA	5
	a. Sample Data Summary (Form I PCDD-1)	[] _		
	b. PCDD/PCDF Toxicity Equivalency Factor (Form I, PCDD-]		
	c. Second Column Confirmation Summary (Form I, PCDD-3) [] _		
	d. Total Homologue Concentration Summary (Form II PCDD)) [] _		
	e. PCDD/PCDF Spiked Sample Summary (Form III PCDD-1) [] _		
	f. PCDD/PCDF Duplicate Sample Summary (Form III PCDD	-2) [] _		
	g. PCDD/PCDF Method Blank Summary (Form IV-PCDD)	[] -		
	h. PCDD/PCDF Window Defining Mix Summary (Form V-PC	CDD-1) []		
	i. Chromatographic Resolution Summary (Form V PCDD-2)	[] .		
	j. PCDD/PCDF Analytical Sequence Summary (Form V PCD	DD-3) [] _		
	k. Initial Calibration (Form VI, PCDD-1, PCDD-2)	[]		
	I. Continuing Calibration (Form VII,PCDD-1, Form VII,PDD-2	2) []		

2.3	GC/MS Displays	YES NO N/A
	Are the following GC/MS displays present?	
	Standard and sample SIM chromatograms. SIM and TIC chromatograms must list date and time of analysis; the	[]
	file name; sample number; and instrument I.D. number b. Percent peak resolution valley c. GC column performance check raw data	
	 d. SIM mass chromatograms must display quantitation ion, confirmation polychlorinated diphenylether ion, where applicable. e. Integrated area and peak height must be listed for all peaks 2.5 times above background f. All peaks must show retention time at the maximum height 	ion, and [] []
2.4	Are the following Chain of Custody Records and in-house Laboratory Control Documents present?	
	 a. EPA Chain of Custody Records b. SMO Sample Shipment Records c. Sample log-in sheets d. GC/MS Standard and Sample Run Log in chronological order e. Sample Extraction Log 	
2.5	Was the sample data package paginated?	[]
	ACTION: If deliverables are missing call the lab for explanation/resubmitted provide missing deliverables, assess the effect on the validity of reviewers narrative.	
3.0	Holding Times	
3.1	Have any of the following holding times been exceeded?	
	For aqueous samples, 30 days from sample collection to extraction	[]
	For soil/sediment samples, 30 days from sample collection to extraction	[]
	For all samples 45 days from time of extraction to time of analysis	1 1

Page: 4 of 24 Date: December 2010 Revision 1.1

Date: December 2010 Revision 1.1 ACTION: If holding times are exceeded, flag all data as estimated ("J"). Holding time criteria do not apply to PE samples. All samples except fish and adipose samples must be stored in dark at 4°C. Fish and Note: adiposetissue must be stored at -20°C in the dark. YES NO N/A 4.0 **Instrument Performance** Mass Calibration - Mass calibration of the MS must be performed prior to analyzing calibration 4.1 solutions, blanks, samples, and QC samples. A static resolving power of at least 10,000 (10% valley definition) must be demonstrated at appropriate masses before any analysis is performed Static resolving power checks must be performed at the beginning and at the end of each 12 hour period of operation. Include in the narrative, minimum required resolving power of 10000 was obtained for perfluorokerosene (PFK) ion 380.9760. This is done by first measuring peak width at 5% of the maximum. This should not exceed 100 ppm, i.e., it should not exceed 0.038, for ion 380.9760. Resolving power, then is calculated using the formula. Resolving Power = m/m = 380.9760/0.038 = 10025. 4.1.1 Was mass calibration performed at the frequency given above? [__] ___ __ 4.1.2 Was the resolving power of PFK ion 380.9760 above 10000, when it was transmitted at the accelerating voltage corresponding to m/z ion 304.9824? 4.2 GC Column Performance Check Solution The GC Column Performance Check solution must contain the first and the last isomers of each homologue PCDD/PCDF, (the internal and recovery standards are optional). The solution also should contain a series of other TCDD isomers for the purpose of documenting the chromatographic resolution. 4.2.1 For analyses on a DB-5 (or equivalent) GC column, the chromatographic resolution is evaluated by the analysis of GC column performance check solution at the beginning of every 12 hour period. Was this performed accordingly? ACTION: If the GC column performance check solution was not analyzed at the required frequency, use professional judgement to determine the effect on the quality of the data.

4.2.2 Were all peaks labeled and identified on the Selected Ion Current Profiles (SICPs)?

Page: 5 of 24

USEPA Region II DV SOP for SW-846 Method 8290 PCDDs/PCDFs using HRGC/HRMS Date: December 2010 Revision 1.1

Page: 5 of 24

	<u>Y</u> 1	ES	NO	N/A
	For DB-5 or equivalent, the peak separation between the unlabeled 2378-TCl representing any other TCDD isomer shall be resolved with a valley of \leq 25 p criteria met?			•
	% Valley = $(x/y) x (100)$			
	Y = The peak height of 2,3,7,8-TCDD isomer			
	 X = The distance from the baseline to the bottom of the valley between peaks. ACTION: If the percent valley criteria are not met qualify all positive data J. D. detects. 		•	
4.2.4	Is the last eluting tetra chlorinated congener (1,2,8,9-TCDD) and the first elutichlorinated congener (1,3,4,6,8-PeCDF) separated properly, since they elute of each other?			econds ——
	ACTION: If one of the congener is missing, report that in the case narrative.			
	Initial 5-Point Calibration - The initial calibration standard solutions (HRCC1 analyzed prior to any sample analysis. They do not have to be analyzed daily continuing calibration standard met all criteria. However, initial calibration sho least once every week and/or whenever the continuing calibration standard docriteria. The calibration standards must be analyzed on the same instrument of GC/MS conditions that were used to analyze the GC column performance che the initial calibration performed at the frequency specified above?	, provuld be nessing	vided to a	he yzed at et all ame
5.1	The following MS/DS conditions must be used:	L	_,	
	Is mass calibration performed as per Section 4.1?	ſ	1	
	Is the total cycle time ≤ 1 second?	[-J]	
	Note: The total cycle time includes the sum of all the dwell times and voltage	rese	t times	 S.
	Were SIM data acquired for each of the ions listed in Table 6, including interfe			
5.1.5	analytical method)	[_]	
5.2	Were the following GC criteria met?			
5.2.1	The chromatographic resolution between the 2378-TCDD and the peaks reprunlabeled TCDD isomers must be resolved with a valley of \leq 25 percent.	esen [ting ar	ny other
5.2.2	In the HRCC3 solution, the chromatographic peak separation between 1,2,3, 1,2,3,6,7,8-HxCDD shall be resolved with a valley of < 50 percent.	4,7,8 []	-HxC[DD and

USEPA Region II DV SOP for SW-846 Method 8290 PCDDs/PCDFs using HRGC/HRMS Date: December 2010 Revision 1.1

YES NO 5.2.3 For all calibration solutions the retention times of the isomers must fall within the retention time windows established by the GC column performance check solution. In addition, the absolute retention times of recovery standards, ¹³C₁₂1234-TCDD and ¹³C₁₂-123789HxCDD shall not change by more than 10 seconds between the HRCC3 analysis and the analysis of any other standard. 5.2.4 The two SIM ions for each homolog must maximize simultaneously and within 3 seconds of the corresponding labeled isomer ions. 5.2.5 The relative ion abundance criteria for PCDDs/PCDFs listed in Table 8 (see analytical method) must be met. 5.2.6 The relative ion abundance criteria for the labeled internal and recovery standards listed in Table 8 must be met. 5.2.7 For all calibration solutions, including HRCC3, the signal to noise ratio (S/N) for the GC signal present in every SICP, including the ones for the labeled standards must be > 10. 5.2.8 The percent relative standard deviations (% RSD) for the the mean response factors (RRF) from the 17 unlabeled standards must not exceed + 20%, and those for the nine labeled reference compounds must not exceed + 30%. []

Page: 7 of 24

- ACTION: 1. If the 25% percent valley for TCDD and 50% valley for HxCDD requirement are not met, quality positive data J. Do not qualify non-detects. The tetra, pentas and hexas (dioxins and furans) are affected. Heptas and Octas are not affected.
 - 2. If the %RSD for each unlabeled isomer exceeds 20%, or the %RSD for each labeled isomer exceeds 30%, flag the associated sample positive results for that specific isomer as estimated ("J"). No effect on the non-detect data.
 - If the ion abundance ratio for an analyte is outside the limits, flag the results for that 3. analyte R (reject).
 - 4. If the ion abundance ratio for an internal or recovery standard falls outside the QC limits flag the associated positive hits with J. No effect on the non-detects.
 - 5. If the signal to noise ratio (S/N) is below control limits, use professional judgement to determine quality of the data.

YES NO N/A

Page: 8 of 24

6. If the selected monitoring ions specified in Table 6 were not used for data acquisition, the lab must be asked for an explanation. If an incorrect ion was used, reject all the associated data.

- 7. If mass calibration criteria as specified in Section 4.1 is not met, specify that in case narrative.
- 8. Non compliance of all other criteria specifiedabove should be evaluated using professional judgement.
- 5.2.9 Spot check response factor calculations and ion ratios. Ensure that the correct quantitation ions for the unlabeled PCDDs/PCDFs and internal standards were used. In addition, verify that the appropriate internal standard was used for each isomer.

To recalculate the response factor, use the equation:

RRFn =
$$(A_n^1 + A_n^2) \times Q_{is}$$

 $(A_{is}^1 + A_{is}^2) \times Q_n$

RRFis =
$$(A_{is}^{1} + A_{js}^{2}) \times Q_{rs}$$

 $(A_{rs}^{1} + A_{rs}^{2}) \times Q_{is}$

Where:

 A_n^1 and A_n^2 = integrated areas of the two quantitation ions of isomer of interest (Table 6).

 A_{is}^{1} and A_{is}^{2} = integrated areas of the two quantitation ions of the appropriate internal standard (Table 6).

 A_{rs}^{1} and A_{rs}^{2} = integrated areas of the two quantitation ions of the appropriate recovery standard (Table 6).

Q_n = quantity of the unlabeled PCDD/PCDF analyte injected (pg)

 Q_{is} = quantity of the appropriate internal standard injected (pg)

 Q_{rs} = quantity of the appropriate recovery standard injected (pg)

Page: 9 of 24 Date: December 2010 Revision 1.1

VEC	NO	N/A
169	NO	IN/A

6.0 Continuing Calibration (HRCC3)	TEO NO NA
The continuing calibration must be performed at the beginning of a 12 hours mass resolution and GC resolution performance checks. A continuing calibration run at the requirement of a 12 hour shift. Was the continuing calibration run at the requirement.	bration is also required at
6.1 Were the following MS/DS conditions used?	
6.1.1 The total cycle time was ≤ 1 second.	[]
6.1.2 SIM data were acquired for each of the ions listed in Table 6 including diphenylether interfering ions (see analytical method).	[_]
6.2 Were the following criteria met?	
6.2.1 For the continuing calibration solution the retention time of the isome retention time windows established by the GC column performance	
6.2.2 The absolute retention time of the recovery standards ¹³ C ₁₂ 1234-TC HxCDD shall not change by more than 10 seconds between the initia HRCC3 standard analyses.	DD and ¹³ C ₁₂ 123679 I HRCC3 and ending []
6.2.3 The two SIM ions for each homolog must maximize simultaneously seconds of the corresponding ions of the labeled isomers.	(<u>+</u> 2 sec) and within 3
6.2.4 For the HRCC3 standard solution, the signal to noise ratio (S/N) for ion shall be greater than 2.5.	the unlabeled PCDD/PCDF
6.2.5 For the internal standards and the recovery standards, the signal to r greater than 10.	noise ratio (S/N) shall be
6.2.6 The relative ion abundance criteria (Table 8 - analytical method) for met.	all PCDD/PCDF shall be
6.2.7 The relative ion abundance criteria for all internal and recovery stand method) must be met.	dards (Table 8 – analytical
6.2.8 The %Difference of RRF of each <u>unlabeled</u> analyte must be within <u>+</u> established during the initial calibration. The measured RRFs for estandards must be within <u>+</u> 30 percent of the mean RRF established	ach of the <u>labeled</u>
	[]

Spot check response factor calculations and ion ratios. Verify that the appropriate quantitation ions for the unlabeled PCDD/PCDFs and internal standards were used.

	YES	<u>NO</u>	N/A
6.2.9 Was the same internal standard used to calculate RRF for each PCDD/PC initial calibration?	DF hom	olog ir	the —
6.2.10 Was the chromatographic peak separation on DB-5 (or equivalent) column 2378-TCDD and the peaks representing any other unlabeled TCDD) isomovalley of ≤ 25 percent?			
6.2.11 Was the chromatographic peak separation between the 123478-HxCDD ar HxCDD in the HRCC3 solution resolved with a valley of <50 percent?	nd the 1 []	23678· 	

Page: 10 of 24

Date: December 2010 Revision 1.1

- ACTION: 1. If any of the requirements listed in sections 6.1.1, 6.1.2, 6.2.1, 6.2.2, and 6.2.9 are not met, use professional judgement to determine the validity of the data.
 - 2. If any requirements listed in sections 6.2.3, 6.2.4, 6.2.5, 6.2.6, and 6.2.7 are not met reject all data (flag R) directly affected by each specific problem.
 - 3. When the %D of the RRF is in between 30% and 50%, all the data for the outlier congeners are flagged J. Data with %D above 50% are rejected (R).
 - 4. If the continuing calibration standard was not analyzed at the required frequency, reject all the data. Contact TPO to initiate reanalysis.
 - 5. If the 25 percent valley (6.2.10) and 50 percent valley (6.2.11) criteria are not met, qualify all positive data with J. Do not qualify non-detects. Note: The tetras, pentas and hexas (dioxins and furans) are affected. Heptas and octas are not affected. If the percent valley is >75 percent and 2378-TCDD is non-detect but 1234-TCDD or an adjacent TCDD isomer is present, the data is questionable. The sample must be reanalyzed. Contact TPO. If the valley criteria for HxCDD are not met, but the valley criteria for TCDD are met or vice-versa, use professional judgement to determine which data must be qualified.
 - 6. If the HRCC3 standard performed at the end of the 12 hour shift did not meet criteria specified in Sections 6.2.1, 6.2.4, 6.2.5, 6.2.6, and 6.2.7, examine the samples which were analyzed prior to this standard and use professional judgement to determine if data qualification is necessary.
 - 7. For all other criteria, use professional judgement.

VEC	NO	NI/A
1 5	NU	IN/A

Page: 11 of 24

6.2.12 To recalculate RRFs for the unlabeled target analytes, and the RRFs for the nine labeled internal standards, use the following equations:

RRFn =
$$(An^1 + An^2) \times Qis$$

(Ais¹ + Ais²) x Qn

RRFis =
$$(Ais^1 + Ais^2) \times Qrs$$

(Ars¹ + Ars²) x Qis

An¹, An², Ais¹, Ais², Ars¹, Ars², Qn, Qis and Qrs are defined in Section 5.2.9.

To calculate percent difference use the following equation:

Where:

RRFi = Relative response factor established during initial calibration

RRFc = Relative response factor established during continuing calibration

7.0 <u>Sample Data</u>

- 7.I Were the following MS/DS conditions used?
- 7.1.1 The total cycle time was \leq 1 second.

___] ___

7.1.2 SIM data were acquired for each of the ions listed in Table 6 (see analytical method) including diphenylether interfering ions.

___] ____

- 7.2 Were the following identification criteria met?
- 7.2.1 For the 2378 substituted isomers found present and for which an isotopically labeled internal or recovery standard is present in the sample extract, the absolute retention time at the maximum peak height of the analyte must be within -1 to 3 seconds of the retention time of the corresponding labeled standard.

[__] ___

7.2.2 For the 2378 substituted isomer reported present, and for which a labeled standard does not exist, the relative retention time (RRT) of the analyte must be within ±.005 RRT units of the RRT established by the continuing calibration standard (HRCC3).

YES NO N/A

Page: 12 of 24

	substituted compounds (tetra throndow established by the GC colunomologue.		
labeled standar	ns listed in Table 6 (analytical meds must be present in the SICP. recovery standards must maximize	The two SIM ions for	or the analyte, the internal
•	l ion current for each characterist imes background noise and mus	_	•
•	l ion current for the internal and rebackground noise.	ecovery standard cha	aracteristic ions must be at
7.2.7 The relative ion present must be	n abundance criteria (Table 8 – a e met.	ınalytical method) for	all PCDDs/PCDFs found
	n abundance criteria for the interrallytical method).	nal and recovery stan	ndards must be met []
	ion of a GC peak as a PCDF can e same time in the corresponding on met?		
	concentration must be within the cring the concentration within the		
concent	alytical method clearly states that trations higher than 10 times the e, high resolution GC/low resolut	upper MCLs should I	
	Reject (flag R) all positive data for Sections 7.2.1, 7.2.2, 7.2.3, and 7	•	do not meet criteria listed in

3. If the requirements listed in section 7.2.6 are not met but all other requirements are met qualify the positive data of the corresponding analytes with "J".

all positive data of the specific analyte with J.

 If the analytes reported positive do not meet ion abundance criteria, section 7.2.7, reject (R) all positive data for these analytes. Change the positive values to EMPC (estimated maximum possible concentration).

2. If the criteria listed in section 7.2.5 are not met but all other criteria are met, qualify

YES NO N/A

Page: 13 of 24

5. If the internal standards and recovery standards do not meet ion abundance criteria (Table 8 - analytical method) but they meet all other criteria flag all corresponding data with "J".

- 6. If PCDF is detected but an interfering PCDPE is also detected (see Section 7.2.9) reject the PCDF data (R). The reported value of PCDF is changed to EMPC.
- 7. If the lab did not monitor for PCDPEs, qualify all positive furan data J.
- 7.2.11 Spot check calculations for positive data and verify that the same internal standards used to calculate RRFs were used to calculate concentration and EMPC. Ensure that the proper PCDDs/PCDFs and internal standards were used.

To recalculate the concentration of individual PCDD/PCDF isomers in the sample use the following equation:

ALL MATRICES OTHER THAN WATER

Cn (pg/g) =
$$\frac{\text{Qis x } (\text{An}^1 + \text{An}^2)}{\text{W x } (\text{Ais}^1 + \text{Ais}^2) \text{ x RRFn}}$$

WATER

Cn (ng/L) =
$$\frac{\text{Qis x } (\text{An}^1 + \text{An}^2)}{\text{V x } (\text{Ais}^1 + \text{Ais}^2) \text{ x RRFn}}$$

Where:

An¹ and An² = integrated ion abundances (peak areas) of the quantitation ions of the isomer of interest (Table 6).

Ais¹ and Ais² = integrated ion abundances (peak areas) of the quantitation ions of the appropriate internal standard (Table 6).

W= Weight (g) of sample extracted V= Volume (ml) of sample extracted

Qis= Quantity (pg) of the appropriate internal standard added to the sample prior to extraction

RRFn= Calculated relative response factor from continuing calibration (see Section 7.7 of the analytical method).

Note: See CLP/SOW DFLMO1.1, Section 15.3 for calculations when any internal standard in a diluted sample is less than 10% of the internal standard area in the continuing calibration standard.

YES NO N/A

Page: 14 of 24

7.3 <u>Estimated Detection Limits (EDL)</u>

- 7.3.1 Was an EDL calculated for each 2,3,7,8-substituted isomer that was not identified regardless of whether other non-2378 substituted isomers were present?

 [___] ____
- 7.3.2 Use the equation below to check EDL calculations:

ALL MATRICES OTHER THAN WATER

EDL (pg/g) =
$$\underline{2.5 \times \text{Qis} \times (\text{Hx}^1 + \text{Hx}^2) \times \text{D}}$$

W x (His¹ + His²) x RRFn

WATER

EDL (ng/L) =
$$2.5 \times \text{Qis } \times (\text{Hx}^1 + \text{Hx}^2) \times \text{D}$$

V x (His¹ + His²) x RRFn

Where:

 Hx^1 and Hx^2 = peak heights of the noise for both quantitation ions of the 2,3,7,8-substituted isomer of interest.

 His^{1} and His^{2} = peak heights of both the quantitation ions of the appropriate internal standards.

D = dilution factor (see Paragraph 10.4.3 of the SOW).

Qis, RRFn, W and V are defined in Section 7.2.11.

NOTE: The validator should check the EDL data to verify that peak heights and not areas were used for this calculation. If the area algorithm was used, the validator should contact the laboratory for recalculation. The TPO must be notified.

7.4 Estimated Maximum Possible Concentration (EMPC)

7.4.1 Was an EMPC calculated for 2378-substituted isomers that had S/N ratio for the quantitation and confirmation ions greater than 2.5, but did not meet all the identification criteria?

[__] ___ __

7.4.2 Use the equation below to check EMPC calculations:

ALL MATRICES OTHER THAN WATER

EMPC (ug/L) =
$$(Ax^1 + Ax^2) x Qis x D$$

$$\overline{(Ais^1 + Ais^2) x RRFn x W}$$

WATER

EMPC (ng/L) =
$$(Ax^1 + Ax^2) \times Qis \times D$$

$$(Ais^1 + Ais^2) \times RRFn \times V$$

Date: December 2010 Revision 1.1

<u>YES</u>

Page: 15 of 24

NO N/A

Where:

 Ax^1 and Ax^2 = areas of both quantitation ions.

Ais¹, Ais², Qis, RRF, W, and V are defined in Section 7.2.11. D is dilution factor defined in Section 10.4.3 of the CLP/SOW.

- Action: 1. If EDL or EMPC of an analyte which was not reported as present is missing, contact the laboratory for correction.
 - 2. If the spot check calculations yielded EDLs or EMPCs different from those reported in Form I,contact the laboratory for an explanation.
 - 3. If EDLs or EMPCs for the most toxic analytes (TEF≥ 0.05) are above CRQLs contact TPO for sample reanalysis.

7.5 Method Blanks

7.5.1	Has a method blank per matrix been extracted and analyzed with each be	oatch of 20 samples?
	If samples of some matrix were analyzed in different events (i.e. different one blank for each matrix been extracted and analyzed for each event?	nt shifts or days) has
7.5.3	Acceptable method blanks must not contain any signal of 2378-TCDD, of equivalent to a concentration of > 20 ppt for soils or 0.2 ppt for water sar	
met?	equivalent to a concentration of > 20 ppt for soils of 0.2 ppt for water sai	
	For other 2378- substituted PCDD/PCDF isomers of each homologue, the concentration in the method blank is less than 1/10 of the upper MCL specified or the area must be less than 2% of the area of the nearest intercriteria met?	ecified in Table 1 of the
7.5.5	For the peak which does not meet identification criteria as PCDD/PCDF area must be less than 5% of the area of the nearest Internal Standard.\	
ACTI	ON: 1. If the proper number of method blanks were not analyzed, notify	the contractor.If they are

unavailable, reject all positive sample data. However, the reviewer may also use professional judgement to accept or reject positive sample data if no blank was run.

Date: December 2010 Revision 1.1

<u>YES</u>

Page: 16 of 24

NO N/A

2. If the method blank is contaminated with 2378-TCDD, 2378-TCDF, 12378PeCDD,12378PeCDF or 23478 PeCDF at a concentration higher than the upper MCL listed in Table 1 of the method, reject all contaminant compound positive data for the associated samples (flag R) and contact the technical project officer to initiate reanalysis if it is deemed necessary.

3. If the method blank is contaminated with any of the above isomers at a concentration of less than the upper MCL specified in the method or of any other 2378-substituted isomer at any concentration and the concentration in the sample is less than five times the concentration in the blank, transfer the sample results to the EMPC/EDL column and cross-out the value in the concentration column. If the concentration in the sample is higher than five times the concentration in the blank, do not take any action.

7.6 Rinsate Blank

7.6.1 One rinsate blank must be collected for each batch of 20 soil samples of is more frequent. Was rinsate blanks collected at the above frequency?	r one per day whichever
7.6.2 Do any rinsate blanks show the presence of 2378-TCDD,2378-TCDF, amounts > .5 ug/L or any other analyte at levels > $1\mu g/L$?	nd 12378PeCDD at
ACTION If any rinsate blank was found to be contaminated with any of the PCDDs/PCDFs notify the technical project officer to discuss what proper action must be taken.	
7.7 <u>Field Blanks</u>	
7.7.1 The field blanks are PEM samples (blind blanks) supplied by EPA from I frequency of one field blank per 20 samples or one per samples collected week, which ever comes first. A typical "field blank" will consist of uncont blanks are used to monitor possible cross contamination of samples in the laboratory.	d over a period of one aminated soil.The field
Were the following conditions met?	
7.7.2 Acceptable field blanks must not contain any signal of 2378-TCDD, 2378 and 12378-PeCDF equivalent to a concentration of > 20 ppt.	8-TCDF, 12378-PeCDD []
7.7.3 For other 2378 substituted PCDD/PCDF isomers of each homologue the a in the field blank is less than the upper MCLs listed in the method.	allowable concentration

YES NO N/A

Date: December 2010 Revision 1.1

Page: 17 of 24

ACTION: When the field blank is found to be contaminated with target compounds, apply the same action as described for the method blank (section 7.5).

NOTE: Contact EPA EMSL/LV to verify that the PEM blank (field blank) did not contain any PCDD/PCDF isomers and ask their assistance in the evaluation of the PE field blank.

8.1	Were the samples spiked with all the internal standards as specified in the method?	[]
8.2	Were internal standard recoveries within the required (40 - 135%) limits?	[]
8.3	If not, were samples reanalyzed?	[]

ACTION:

- 1. If the internal standard recovery was below 25 percent, reject (R) all associated non detect data (EMPC/EDL) and flag with "J" all positive data.
- 2. If the internal standard recovery is above the upper limit (135 percent) flag all associated data (positive and non-detect data) with "J".
- 3. If the internal standard recovery is less than 10%, qualify all associated data R (Reject). when highly toxic isomers (TEF≥ 0.05) are affected, notify TPO to initiate reanalysis.

Recalculate the percent recovery for each internal standard in the sample extract, Ris, using the formula:

Ris =
$$(Ais^1 + Ais^2 \times Qrs \times 100\%)$$

(Ars¹ + Ars² x RRFis x Qis

Ais¹, Ais², Ars¹, Ars₂, Qis, Qrs and RRFis are defined, previously.

9.0 Recovery Standards

There are no contractual criteria for the Recovery Standard area. However, because it is very critical in determining instrument sensitivity, the <u>Recovery Standard</u> area must be checked for every sample.

9.1	Are the recovery standard areas for every sample and black	ank within the upper and lower limits of
	each associated continuing calibration? Area upper limit=	+100% of recovery standard area .Area
	lower limit= -50% of recovery standard area.	[]

validity.

			<u>YES</u>	<u>NO</u>	<u>N/A</u>
9.2		retention time of each recovery standard within onds of the associated daily calibration standard?	Į.] _	
AC	TION: 1	 If the recovery standard area is outside the upper or lower and non-detect data (EMPC/EDL) with "J" regardless wh recoveries met specifications or not. 			
		extremely low area counts (<25%) are reported flag all assousable (R) and the positive data J.	ociated non-de	tect da	ta as
	ca m	the retention time of the recovery standard differs by more talibration use professional judgement to determine the effectore than 10 seconds may cause certain analytes to elute out indow established by the GC column performance check so	t on the result utside the rete	s. A tim	ne shift of
10.0	PEM Ir	nterference Fortified Blanks			
10.1	LV, is of sample week p	nown blank usually an interference fortified soil/sediment sandesignated by the sampling team for the laboratory for spiking is one per group of 20 environmental samples or one per speriod, whichever occurs first. The sample is spiked by the less of the matrix spiking solution and then extracted and analy	ng.The freque samples colled laboratory with	ncy of toted oven	his QC er one ppropriate
10.2		fortified PEM blank analyzed at the frequency bed above?	[
10.3		ne percent recovery of 2378-TCDD and other 2378-substitutercent control limits?	ed compound [s within	the 50 to
ACT	ION:	1. If the recovery of a 2,3,7,8-substituted isomer falls outs limit, flag all positive and non-detect data of the same and homolog series with J. However, if the recovery is below 2 detects R.Notify the TPO. Reanalysis may be initiated.	related isome	rs in the	e same

Page: 18 of 24

Date: December 2010 Revision 1.1

NOTE: This blank, as prescribed above in Section 10.1, however, is not given in the analytical method.

2. If no fortified PEM blank was analyzed, use professional judgement to assess data

USEPA Region II DV SOP for SW-846 Method 8290 PCDDs/PCDFs using HRGC/HRMS Date: December 2010 Revision 1.1

												YES		<u>NO</u>	N/A
11.0	<u>Matrix</u>	x Spik	e (Field	Sample	<u>e)</u>										
11.1	Was a	a matr	x spike	analyze	d at the	e frequ	uency	of one	e per S	SDG sa	amples	per ma [trix? _]		
11.2	Was th			covery o	f 2378	-TCDD) and	other 2	2378-	substi	tuted F	PCDDs/F []	PCD —	Fs wi	ithin 50 t
AC	CTION:	qualiof the	y of the spiked	data. T	he 50- only. T	-150% I The ma	limits atrix s	of the	matri ata of	x spik the Pl	e data		use	ed to f	sess the lag data e
12.0	Enviro	<u>onme</u>	ntal Du	olicate S	Sample	<u>es</u>									
	For ev less, th above f	here r	nust be							•		ne week, amples [
	id result mers ar				•	_		25% re	elative	differ	ence fo	or 2,3,7,8 	3 su _	bstitu —–	ited
AC	CTION:		•					•				ata. If n recoveri		ts are)
13.0	<u>Perfor</u>	rman	e Evalı	ıation S	Sample	<u>es</u>									
amo PE s	Include ounts of samples o sample	unlab s are p	eled 237 rovided	'8-TCD[by the F	O or a r Region	mixture n, and m	e of 23 must b	378-T0 be ana	CDD a llyzed	and oth at the	er PĊl freque	DD/PCD ncy of o	F is	omer set pe	s. The er batch
13.2	. Were t	the ar	alytical	results v	within t	the EPA	A 99%	% acce	eptand	e crite	ria?				
AC	CTION:			amples ine for P			dated a	as if th	ey we	ere env	vironme	ental sar	nple	es. Th	ere is no

Page: 19 of 24

2. PE samples containing only 2378-TCDD When 2378-TCDD was not qualitatively identified, or if the reported concentration is outside the 99% acceptance window all positive and negative (EMPC/EDL) data for all associated samples are rejected.

USEPA Region II DV SOP for SW-846 Method 8290 PCDDs/PCDFs using HRGC/HRMS Page: 20 of 24 Date: December 2010 Revision 1.1

<u>YES</u>	NO	N/A

- 3. PE samples containing a mixture of PCDD/PCDF isomers
 - When the reported concentration of any analyte is outside the EPA 99% confidence interval, all positive and negative (EMPC/EDL) data of the 2378 substituted isomers within the same homologue for all associated samples are rejected.
- 4. When PCDD/PCDF data are rejected because of PE results, the EPA technical project officer must be notified. Reanalysis may be initiated.
- 5. For PE blind blanks see Section 7.7 (Field Blanks).

14.0 <u>Second Column Confirmation</u>

14.1	Was a second column confirmation performed?	[]		
	Was the sample extract reanalyzed on a 30 m DB-225, fused silica capillated. TCDF using the GC/MS conditions given in Section 7.9.7.1.2 of the analy	-	•	2,3,7,8
	NOTE: The concentration of 2,3,7,8 TCDF obtained from the primary cobe used for qualification, due to better QC data associated with the primary that the confirmation and quantitation of 2,3,7,8-TCDD may be accomplised to be accomplisation.	ary colu	mn. É Al	so note
	ACTION: If confirmation is missing, use professional judgement, or conf	act TPC) for as	sistance.
14.3	Did the second column meet the calibration and linearity specification in above?	Section []	s 5.0 a ——	nd 6.0
14.4	Was the % D of the quantitation results of the two columns less than 50%	? []		

15.0 <u>Sample Reanalysis</u>

15.1 The Region II TPO will evaluate the need for reanalyzing the samples with qualified data based on site-specific Regional Data Quality Objectives. The rerun may be billable or non billable as specified in the SOW. SMO should be notified of all reruns.

YES NO N/A

Page: 21

Date: December 2010 Revision 1.1

of 24

15.2 Due to a variety of situations that may occur during sample analysis the laboratory is required to reanalyze or reextract and reanalyze certain samples. If a reanalysis was required but was not performed, contact TPO to initiate reanalysis. List below all reextractions and reanalyses and identify the PCDD/PCDF sample data summaries (Form I) which must be used by the data user (when more than one is submitted).

16.0 Isomer Specificity and Toxicity Equivalency Factor (TEF) -

When calculating the 2378-TCDD Toxicity Equivalency of a sample only those 2378 substituted isomers that were positively identified in the sample must be included in the calculations. The sum of the TEF adjusted concentration is used to determine when a second column confirmation is required to achieve isomer specificity.

16.1	Did the lab include EMPC or EDL values in the toxicity equivalency calculation [_	_		
	Were all samples, whose toxicity equivalency exceeded the required values a confirmation column to establish isomer specificity?		•	zed on

ACTION: 1. If the toxicity equivalency calculations were not performed properly notify TPO.

2. If the toxicity equivalency exceeded the required limits (0.7 ppb for soil/sediment, 7ppt for aqueous and 7ppb for chemical waste samples), and the lab failed to reanalyze the samples on a specific secondary column, notify TPO.

	PCDFs/PCDDs Data Assessment							
CASE NO	LABORATORYSAMPLE NO							
DATA ASSESSMENT:								
"J" (estimated). Rejected data do	except those values which have been qualifices not imply the analyte is not present. It me id and it provides no information as to whether	ans that due to significant						
All action is detailed below and on	the attached sheets.							
Reviewer's Signature:		Date://20						
Verified By:	Date://20							

Page: 22 of 24

Date: December 2010 Revision 1.1

USEPA Region II DV SOP for SW-846 Method 8290 PCDDs/PCDFs using HRGC/HRMS Date: December 2010 Revision 1.1

Page: 23 of 24

Case#	_	
Site:		
Lab:		

Overall Assessment

USEPA Region II DV SOP for SW-846 Method 8290 PCDDs/PCDFs using HRGC/HRMS Date: December 2010 Revision 1.1

Page: 24 of 24

Case#	_	
Site:		
Lab:	_	

Contract Problems/Non-Compliance