STANDARD OPERATING PROCEDURE FOR THE VALIDATION OF ORGANIC DATA ACQUIRED USING METHOD 524.2(Revision 4.1, 1995) MEASUREMENT OF PURGEABLE ORGANIC COMPOUNDS IN WATER BY GAS CHROMATOGRAPHY/MASS SPECTROMETRY (GC/MS) CAPILLARY COLUMN



Prepared by: RusseM Arnone, Chemist	Date:	01/12/1
Peer Reviewed by: Muhammad Sheikh, Chemist Hazardous Waste Support Section	Date:	01/12/11
Concurred by: Michael Mercado, Arting Chief	Date	01/12/1
Approved by: Robert Runyon, Chief Hazardous Waste Support Branch	Date:	01/13/11
Annual Review Reviewed by:	Date:	
Name Reviewed by:	Date:	
Name		

TABLE OF CONTENTS

TITLE PAG	E	1
TABLE OF	CONTENTS	2
INTRODUCT	ION	3
I.PACKAGE	COMPLETENESS AND DELIVERABLES	4
1.0	Data Completeness and Deliverables	4
2.0	Cover Letter, SDG Narrative	4
II.VOLATI	LE ANALYSES	4
1.0	Traffic Reports and Laboratory Narrative	4
2.0	Holding Times	
3.0	Surrogate Recovery (CLP Form II Equivalent)	6
4.0	Laboratory Fortified Blanks (CLP Form III Equivalent)	7
5.0	Laboratory Fortified Sample Matrix (LFM)	8
6.0	Laboratory Reagent Blanks (LRB)	9
7.0	Contamination	0
8.0		
9.0	GC/MS Instrument Performance Check (CLP Form V Equivalent)1	. 2
10.0		
11.0	Tentatively Identified Compounds (TIC) (CLP Form I/TIC Equivalent) 1	
	Compound Quantitation and Reported Detection Limits	
	Standards Data (GC/MS)	
14.0	GC/MS Initial Calibration (CLP Form VI Equivalent)	. 8
15.0	GC/MS Calibration Verification (CLP Form VII Equivalent)	. 9
	Internal Standards (CLP Form VIII Equivalent)	
	Field Duplicates	
DEFINITIO	NS	22
	nyms:	
	Oualified Definitions:	

STANDARD OPERATING PROCEDURE FOR THE VALIDATION OF ORGANIC DATA ACQUIRED USING METHOD 524.2(Revision 4.1, 1995) MEASUREMENT OF PURGEABLE ORGANIC COMPOUNDS IN WATER BY GAS CHROMATOGRAPHY/MASS SPECTROMETRY (GC/MS) CAPILLARY COLUMN



Prepared by:			
	Date:		
	Russell Arnone, Chemist		
	Hazardous Waste Support Section		
Peer Reviewed	by:	Date:	
	Muhammad Sheikh, Chemist		
	Hazardous Waste Support Section		
			
Concurred by:_		Date	
	Michael Mercado, Acting Chief		
	Hazardous Waste Support Section		
Approved by: _		Date:	
	Robert Runyon, Chief		
	Hazardous Waste Support Branch		
	Annual Review		
Reviewed by: _		_Date:	
	Name		
Reviewed by: _		_ Date:	
	Name		

TABLE OF CONTENTS

TITLE PAG	距	1
TABLE OF	CONTENTS	2
INTRODUC	TION	3
T PACKAGI	COMPLETENESS AND DELIVERABLES	4
1 0	<u>Data Completeness and Deliverables</u>	4
2.0	Cover Letter, SDG Narrative	4
2.0	<u> </u>	• • •
II.VOLAT	ILE ANALYSES	4
1.0	Traffic Reports and Laboratory Narrative	
2.0	Holding Times	
3.0	Surrogate Recovery (CLP Form II Equivalent)	6
4.0	Laboratory Fortified Blanks (CLP Form III Equivalent)	
5.0	Laboratory Fortified Sample Matrix (LFM)	
6.0	Laboratory Reagent Blanks (LRB)	9
7.0		
8.0		
9.0	GC/MS Instrument Performance Check (CLP Form V Equivalent)	.12
10.0	Target Analytes (CLP Form I Equivalent)	
11.0	Tentatively Identified Compounds (TIC) (CLP Form I/TIC Equivalent	<u>)</u> 15
12.0	Ocompound Quantitation and Reported Detection Limits	.17
	Standards Data (GC/MS)	
14.0	O GC/MS Initial Calibration (CLP Form VI Equivalent)	.18
15.0	GC/MS Calibration Verification (CLP Form VII Equivalent)	.19
16.0	Internal Standards (CLP Form VIII Equivalent)	.20
17.0	Field Duplicates	.21
חביבדאדיידי/	ONS	2.2
	onyms:	
	onyms:	

INTRODUCTION

Scope and Applicability

This SOP offers detailed guidance in evaluating laboratory data generated according to the USEPA Method 524.2. The validation methods and actions discussed in this document are based on the requirements set forth in USEPA Method 524.2 and "USEPA Contract Laboratory Program National Functional Guidelines for Organic Data Review", October 1999 (EPA - 540/R-99-008). This document covers technical as well as method specific problems; however situations may arise where data limitations must be assessed based on the reviewer's own professional judgement.

Summary

To ensure a thorough evaluation of each result in a data case, the reviewer must complete the checklist within this SOP, answering specific questions while performing the prescribed "ACTIONS" in each section. Qualifiers (or flags) are applied to questionable or unusable results as instructed. The data qualifiers discussed in this document are defined on page 23.

The reviewer must prepare a detailed data assessment to be submitted along with the complete SOP checklist. The Data Assessment must list all data qualifications, reasons for qualifications, instances of missing data, and contract noncompliance.

I. PACKAGE COMPLETENESS AND DELIVERABLES
--

CASE NUMBER:	LAB:
SITE NAME:	

YES NO NA

- 1.0 Data Completeness and Deliverables
 - 1.1 Has all data been submitted in CLP deliverable

US EPA Region II Date: November 2010 Method 524.2 (Rev.4.1, 1995) SOP HW-29, Rev. 2 format or CLP Forms Equivalent? ____ If not, note the effect on review of the data in the Data Assessment narrative. Cover Letter, SDG Narrative 2.0 2.1 Is a laboratory narrative, signed release, or Ц____ cover letter present? 2.2 Are case number and SDG number(s) contained in the [] ____ narrative or cover letter? II. **VOLATILE ANALYSES** Traffic Reports and Laboratory Narrative 1.0 1.1 Are the Traffic Reports, Chain of Custodies, or signed releases from the field samplers present for all samples? <u>[]</u> ____ If no, contact the laboratory/sampling team ACTION: for replacement of missing or illegible copies. 1.2 Is a sampling trip report present (if required)? 1.3 Sample Conditions/Problems 1.3.1 Do the Traffic Reports, Chain of Custodies, or Lab Narrative indicate any problems with sample receipt, condition of samples, analytical problems or special notations affecting the quality of the data? If all the VOA vials for a sample have air ACTION: bubbles or the VOA vial analyzed had air bubbles, flag all positive results "J" and all non-detects "R". ACTION: If samples were not iced or if the ice was melted upon receipt at the laboratory and the temperature of the cooler was elevated (>10°C), flag all positive results "J" and all non-detects "UJ". YES NO NA Holding Times 2.0 Have any volatile holding times, determined from date of collection to date of analysis, been ___ [] exceeded?

4

The holding time for aqueous samples is 14 days.

Date: November 2010 Method 524.2 (Rev.4.1, 1995) SOP HW-29, Rev. 2

If unpreserved, aqueous samples maintained at NOTE: 4°C for aromatic hydrocarbons analysis must be analyzed within 7 days. If preserved with acid to a pH <2 and stored at 4°C, then aqueous samples must be analyzed within 14 days from whether or not samples were preserved.

If holding times are exceeded, flag all ACTION: positive results as estimated ("J") and sample quantitation limits as estimated ("UJ"), and document in the narrative that holding times were exceeded.

> If analyses were done more than 14 days beyond holding time, either on the first analysis or upon re-analysis, the reviewer must use professional judgement to determine the reliability of the data and the effects of additional storage on the sample results. At a minimum, all results should be qualified "J", but the reviewer may determine that non-detect data are unusable ("R"). If holding times are exceeded by more than 28 days, all non-detect data are unusable (R).

3.0 <u>Surrogate Recovery (CLP Form II Equivalent)</u>

3.1		the volatile surrogate recoveries been ed on Surrogate Recovery forms?			
3.2		o, are <u>all the samples listed</u> on the opriate Surrogate Recovery forms?	<u>[]</u>		
ACTI	ON:	If large errors exist, deliverables are unavailable or information is missing, document the effect(s)in Data Assessments and contact the laboratory/project officer/appropriate official for an explanation/ resubmittal, make any necessary corrections and document effect in the Data Assessment.	d YES	NO	NA
3.3	Were	outliers marked correctly with an asterisk?	11		
ACTI	ON:	Circle all outliers with a red pencil.			
3.4 Were one or more volatile surrogate recoveries outside required limits for any sample or method blank (Surrogate recovery is 70-130% for aqueous samples)			<u>[]</u>		

	NOTE:	crit	can use their developed in house acceptance ceria, (See Method 8000B Sect.8.7) if none, use 70-130%.			
		If y	ves, were samples reanalyzed?	11		
		Were	e method blanks reanalyzed?	<u>[]</u>		
	ACTIO	N:	If all surrogate recoveries are > 10% but 1 or more compounds do not meet method specifications:			
	;	1. 2. 3.	Flag all positive results as estimated ("J") Flag all non-detects as estimated detection limits ("UJ") when recoveries are less than the lower acceptance limit. If recoveries are greater than the upper acceptance limit, do not qualify non-detects.			
			If any surrogate has a recovery of < 10%:			
		1. 2.	Positive results are qualified with ("J"). Non-detects for that should be qualified as unusable ("R").			
	NOTE:		Professional judgement should be used to qualify data that have method blank surrogate recoveries out of specification in both original and reanalyses. Check the internal standard areas.	:		
			there any transcription/calculation errors een raw data and reported data?		11	
	ACTIO	N:	If large errors exist, take action as specified in section 3.2 above.	YES	NO	NA
4.0	Labora	atory	Fortified Blanks (CLP Form III Equivalent)			
		(LFB)	the volatile Laboratory Fortified Blanks recoveries been listed on the laboratory ting form?	11		
	NOTE:		If the data has not been reported, then contact the laboratory/project officer to obtain the information necessary to evaluate the spike recoveries in the MS, MSD, and LFB. The required data which should have been provided by the lab include the analytes and concentrations used for			

spiking, background concentrations of the spiked analytes (i.e., concentrations in unspiked sample), methods and equations used to calculate the QC acceptance criteria for the spiked analytes, percent recovery data for all spiked analytes.

The data reviewer must verify that all reported Equations and percent recoveries are correct before proceeding to the next section.

NOTE: The LFB spike is spiked with the same analytes at the same concentrations as a calibration standard (Method 524.2-16, Sect.9.3) if different, make note in Data Assessment.

4.2	Were Laboratory Fortified Blanks analyzed at	
	the required frequency (1 LFB per 20 samples)?	_[]

ACTION: If any LFB data are missing, take the action specified in section 3.2 above.

4.3 How many LFB volatile spike recoveries are outside QC Limits?

Water ____ out of ____

ACTION: Circle all outliers with a red pencil.

4.4 Were one or more of the volatile LFB recoveries 70-130% recovery as per Method 524.2-17, Sect.9.6 outside

YES NO NA

ACTION: 1. If the recovery is > upper in-house limit (or 130%), only positive values for the Affected analytes of the compound(s) are flagged "J".

2. If the recovery is < lower in-house limit (or 70%), flag positive values for the Affected analytes of the compound(s) "J" and non detects "J".

NOTE: All analytes in associated sample results are qualified for the following criteria:

1. If 25% of the LFB recoveries were < lower in-house limit (or 70%) qualify all positive results "J" and all non-detects "R".

 If two or more LFB recoveries were < 10% qualify all positive results "J" and all non-detects "R".

5.0 Laboratory Fortified Sample Matrix (LFM)

NOTE: Analysis of a laboratory fortified sample matrix (LFM) is required ONLY if the criteria in section 9.4 are not met. "The integrated areas of the quantitation ions of the internal standards and surrogate in all samples, continuing calibration checks and blanks should remain reasonably constant over time". An abrupt change may indicate a matrix effect and a laboratory fortified duplicate sample must be analyzed to test for matrix effect.

5.1 Have the volatile Laboratory Fortified Sample Matrix (LFM) recoveries been listed on the laboratory reporting form?

[] _____

NOTE: The required data which should have been provided by the lab include the analytes and concentrations used for spiking, background concentrations of the spiked analytes (i.e., concentrations in unspiked sample), methods and equations used to calculate the QC acceptance criteria for the spiked analytes, percent recovery data for all spiked analytes.

The data reviewer must verify that all reported Equations and percent recoveries are correct before proceeding to the next section.

NOTE: The laboratory should use one matrix spike and a duplicate analysis of an unspiked field sample if target analytes are expected in the sample. If the sample is not expected to contain target analytes, a Laboratory Fortified Duplicate Sample (LFM) should be analyzed (Method 524.2-17, Sect.9.4)

ACTION: No action is taken on LFM data alone. However using professional judgement, the validator may use the LFM results in conjunction with other QC criteria and qualify data for that matrix following the guidelines addressed in Sections 4.3 to 4.4.

6.0 <u>Laboratory Reagent Blanks (LRB)</u>

6.1 Is the LRB Summary form present?

	<i>c</i> 0					
	6.2	Has samp	uency of Analysis: a Laboratory reagent blank been reported for les of similar matrix, or concentration level for each extraction batch?	, []		
	6.3	Has used	a LRB been analyzed for each GC/MS system ?			
	ACTI(ON:	If any LRB data are missing, take action as specified in section 3.2. If not available use professional judgement to determine if the associated sample data should be qualified.			
	6.4	chro	matography: review the blank raw data - matograms (RICs), quant reports or data system touts and spectra.	m		
		stab	he chromatographic performance (baseline ility) for each instrument acceptable for the tiles?			
	ACTI(ON:	Use professional judgement to determine the effect on the data.			
7.0	Cont	<u>amina</u>	<u>tion</u>	YES	NO	NA
			there field reagent blanks (FRB) associated y sample?	<u></u>		
	ACTI	ON:	If no, note in Data Assessment that there is no associated field reagent blank. For analy with high concentrations, use professional judgement on qualification of these values and make note in Data Assessment. Duplicate FRB's must be handled along with each sample set, which is composed of the samples collected from the same general site at approximately the same time.	tes		
	7.2	blan	ny Laboratory reagent blank/Field reagent ks have positive results for target analytes or TICs?			
		conc	applied as described below, the contaminant entration in these blanks are multiplied by sample dilution factor.			

ACTION: Prepare a list of the samples associated

Date: November 2010

Method 524.2 (Rev.4.1, 1995) SOP HW-29, Rev. 2

with each of the contaminated blanks. (May attach a separate sheet.)

NOTE:

All field reagent blank results associated with a particular group of samples (may exceed one per case) must be used to qualify data. Blanks may not be qualified because of contamination in another blank. Field reagent blanks/ Laboratory reagent blanks must be qualified for outlying surrogates, poor

spectra, instrument performance or calibration

QC problems.

ACTION:

Follow the directions in the table below to qualify sample results due to contamination. Use the largest value from all the associated

blanks.

US EPA Region II

Date: November 2010 Method 524.2 (Rev.4.1, 1995) SOP HW-29, Rev. 2

			YES NO NA
	Sample conc > CRQL but < 10x blank value	Sample conc < CRQL & <10x blank value	Sample conc > CRQL & >10x blank
Methylene Chloride Acetone Toluene 2-Butanone	Flag sample result with a "U"	Report CRQL & qualify "U"	No qualification is needed
	Sample conc > CRQL but < 5x blank	Sample conc < CRQL & is < 5x blank value	Sample conc > CRQL value & > 5x blank
Other contam- inants	Flag sample result with a "U"	Report CRQL & qualify "U"	No qualification is needed

NOTE: The reporting of TIC compounds may or may not be required.

For TIC compounds, if the concentration ACTION: in the sample is less than five times the concentration in the most contaminated associated blank, flag the sample data "R" unusable.

8.0 GC/MS Apparatus and Materials

Did the lab use the proper gas chromatographic column(s) for analysis of volatiles by Method 524.2? Check raw data, instrument logs or contact the lab to determine what type of column(s) was (were) used.

For the analysis of volatiles, the method requires the use of 60 m. x 0.75 mm capillary column, coated with VOCOL(Supelco) or equivalent column. (Method 524.2-9, Sect. 6.3.2)

ACTION: If the specified column, or equivalent, was not used, document the effects in the Data Assessment. Use professional judgement to determine the acceptability of the data.

> YES NO NA

9.0 <u>GC/MS Instrument Performance Check (CLP Form V Equivalent)</u>

US EPA Region II

Date: November 2010 Method 524.2 (Rev.4.1, 1995) SOP HW-29, Rev. 2 9.1 Are the GC/MS Instrument Performance Check forms present for Bromofluorobenzene (BFB), and do these forms list the associated samples with date/time analyzed? [] 9.2 Are the enhanced bar graph spectrum and mass/charge (m/z) listing for the BFB provided for each twelve hour shift? [] 9.3 Has an instrument performance check solution (BFB) been analyzed for every twelve hours of sample analysis per instrument? (Method 524.2-18, Sect. 10.1) [] ACTION: List date, time, instrument ID, and sample analyses for which no associated GC/MS tuning data are available. DATE TIME INSTRUMENT SAMPLE NUMBERS If the laboratory/project officer/appropriate ACTION: all data generated outside an acceptable twelve hour calibration interval.

If mass assignment is in error, flag all ACTION: associated sample data as unusable, ("R").

- 9.4 Have the ion abundances been normalized to m/z 95? [_] ____
- 9.5 Have the ion abundance criteria been met for each [] ____ instrument used?

List all data which do not meet ion abundance ACTION: criteria (attach a separate sheet).

ACTION: If ion abundance criteria are not met, take action as specified in section 3.2.

9.6 Are there any transcription/calculation errors

between mass lists and reported values? (Check at least two values but if errors are found, check more.)

9.7 Have the appropriate number of significant Figures (two)been reported?

NO NA YES

Date: November 2010

SOP HW-29, Rev. 2

US EPA Region II Method 524.2 (Rev.4.1, 1995)

If large errors exist, take action as ACTION: specified in section 3.2. 9.8 Are the spectra of the mass calibration compound acceptable? Use professional judgement to determine ACTION: whether associated data should be accepted, qualified, or rejected. 10.0 Target Analytes (CLP Form I Equivalent) 10.1 Are the Organic Analysis reporting forms present with required header information on each page, for each of the following: a. Samples and/or fractions as appropriate Laboratory Fortified Sample Matrix b. Blanks C. d. Laboratory Fortified Blank [] ___ __ 10.2 Are the Reconstructed Ion Chromatograms, mass spectra for the identified compounds, and the data system printouts (Quant Reports) included in the sample package for each of the following? Samples and/or fractions as appropriate a. b. Laboratory Fortified Sample Matrix (Mass spectra not required) Blanks c. d. Laboratory Fortified Blanks [] ACTION: If any data are missing, take action specified in 3.2 above. YES NO NA 10.3 Is chromatographic performance acceptable with respect to: Baseline stability? Resolution? [] ___ [Peak shape?

		Date	e: 1	Novem	nber	2010
		SOP	HW-	-29,	Rev.	2

Full	Full-scale graph (attenuation)?						
Othe	er:						
ACTION:	Use professional judgement to determine the acceptability of the data.						
	the lab-generated standard mass spectra of ntified volatile compounds present for each ole?						
ACTION:	If any mass spectra are missing, take action specified in 3.2 above. If the lab does not generate their own standard spectra, make a note in the Data Assessment. If spectra are missing, reject all positive data.						
RRT	the RRT of each reported compound within 0.06 units of the standard RRT in the continuing bration?						
at a most	6 Are all ions present in the standard mass spectrum at a relative intensity greater than 10% (of the most abundantion) also present in the sample mass spectrum?						
ions cor	the relative intensities of the characteristics in the sample agree within ± 30% of the responding relative intensities in the erence spectrum?						
ACTION:	Use professional judgement to determine acceptability of data. If it is determined that incorrect identifications were made, all such data should be rejected ("R"), flagged ("N") - Presumptive evidence of the presence of the compound) or changed to non detected ("U") at the calculated		NO	NA			
	detection limit. In order to be positively identified, the data must comply with the criteria listed in 9.6, 9.7, and 9.8.						
ACTION:	When sample carry-over is a possibility, Professional judgement should be used to determine if instrument cross-contamination has affected any Positive compound identifica-	ation					

11.0 Tentatively Identified Compounds (TIC) (CLP Form I/TIC Equivalent)

NOTE: Use this section only if TIC are required.

11.1	form numb	all Tentatively Identified Compound reporting as present; and do listed TIC's include scan per or retention time, estimated concentration a qualifier?	
NOTE	:	Add "N" qualifier to all TIC's which have CAS number, if missing.	
11.2	comp incl	the mass spectra for the tentatively identified ounds and associated "best match" spectra uded in the sample package for each of the owing:	
	a.	Samples and/or fractions as appropriate []	
	b.	Blanks []	
ACTIO	ON:	If any TIC data are missing, take action specified in 3.2 above.	
ACTIO	ON:	Add "JN" qualifier only to analytes identified by a CAS #.	
NOTE		TIC's are present in the associated blanks take tion as specified in section 7.2 above.	
11.3		any priority pollutants listed as TIC compounds e., an BNA compound listed as a VOA TIC)? []	
ACTIO	ON:	If yes, document in the data assessment that non VOA Compounds are present in the sample(s).	
11.4	with	all ions present in the reference mass spectrum a a relative intensity greater than 10% (of the abundant ion) also present in the sample mass YES NO	NA
	spec	trum? []	
11.5		PIC and "best match" standard relative ion ensities agree within \pm 20%?	
ACTIO	ON:	Use professional judgement to determine acceptability of TIC identifications. If it is determined that an incorrect identification was made, change the identification to "unknown" or to some less specific identification (example: "C3 substituted benzene") as appropriate. Also, when a compound is not found in any blank, but is a suspected artifact of a common laboratory contaminant, the result should be qualified as unusable, "R". (Common	

lab contaminants: CO₂(M/E 44), Siloxanes (M/E 73), Hexane, Aldol Condensation Products, Solvent Preservatives, and related byproducts).

12.0

Comp	ound	Quantitation and Reported Detection Limits		
12.1	orgar least inter initi	chere any transcription/calculation errors in nic analysis reporting form results? Check at two positive values. Verify that the correct rnal standard, quantitation ion, and average all RRF/CF were used to calculate organic rsis reporting form result. Were any errors d?		
NOTE:	insubetwas irawincl	ictural isomers with similar mass spectra, but afficient GC resolution (i.e. percent valley ween the two peaks > 25%) should be reported someric pairs. The reviewer should check the data to ensure that all such isomers were suded in the quantitation (i.e., add the areas the two coeluting peaks to calculate the total centration).		
12.2		the method CRQL's adjusted to reflect sample ions?		
ACTIC	N:	If errors are large, take action as specified in section 3.2 above.		
ACTIC	DN:	When a sample is analyzed at more than one dilution, the lowest detection limits are used (unless a QC exceedance dictates the use of YES the higher detection limit from the diluted sample data). Replace concentrations that exceed the calibration range in the original analysis by crossing out the "E" and it's associated value on the original reporting form (if present) and substituting the data from the analysis of the diluted sample. Specify which organic analysis reporting form is to be used, then draw a red "X" across the entire page of all reporting formsnthat should not be used, including any in the summary package.	NO	NA

13.0 Standards Data (GC/MS)

13.1 Are the Reconstructed Ion Chromatograms, and data System printouts (Quant Reports) present for initial and continuing calibration?

	ACTION:	If any calibration standard data are missing take action specified in section 3.2 above	J,		
14.0	GC/MS Ir	nitial Calibration (CLP Form VI Equivalent)			
	14.1 Are and	ent []			
	ACTION:	If any calibration forms or standard raw data are missing, take action specified in section 3.2 above.			
	14.2 Are	all average RRFs > 0.050?	11		
	ACTION:	Circle all outliers with red pencil.			
	ACTION:	For any target analyte with average RRF < 0.05, qualify all positive results for that analyte "J" and all non-detect results for that analyte "R".			
	ranç devi	response factors stable over the concentration of the calibration. The % relative standard fation (%RSD) $\leq 20.0\%$ as per Method 524.2-20, t. 10.2.6.1.			
	ACTION:	Circle all outliers with a red pencil.			
	ACTION:	If the % RSD is > 20.0%, qualify positive results for that analyte "J" and non-detects using professional judgement. When RSD > 90 qualify all positive results for that analyte "J" and all non-detect results for that analyte "R".)응,	NO	NA
	CC	nalytes previously qualified "U" due to blank ontamination are still considered as "hits" nen qualifying for calibration criteria.			
	14.4 Was	the % RSD determined using RRF or CF?			
	line	no, what method was used to determine the earity of the initial calibration? Document effects to the case in the Data Assessment.			
	The	there any transcription/calculation errors in reporting of RRF or % RSD? (Check at least values but if errors are found, check more.)	<u> </u>		

US EPA Region II Date: November 2010 Method 524.2 (Rev.4.1, 1995) SOP HW-29, Rev. 2

ACTION: Circle errors with a red pencil.

ACTION: If errors are large, take action as

specified in section 3.2 above.

15.0	GC/MS	Calibration	Verification	(CLP	Form	VII	Equivalent)
------	-------	-------------	--------------	------	------	-----	-------------

- 15.1 Are the Calibration Verification reporting forms present and complete for all compounds of interest?[]
- 15.2 Has a calibration verification standard been analyzed for every twelve hours of sample analysis per instrument? [] ____
- The mean response factors calculated during Initial calibration are used for sample quantitation (Method 524.2-26, Sect.12.1.1).
- If any forms are missing or no calibration ACTION: verification standard has been analyzed twelve hours prior to sample analysis, take action as specified in section 3.2 above. If calibration verification data are not available, flag all associated sample data as unusable ("R").
- 15.3 Was the % D determined from the calibration verification YES NO NAdetermined using RRF and by CF? []

If no, what method was used to determine the calibration verification? Document any effects to the case in the Data Assessment.

Do any volatile compounds have a % D (difference or drift) between the initial and continuing RRF or CF which exceeds 30% (Method 524.2-21, Sect. 10.3.5).

___ [_] ___

15.4

ACTION: Circle all outliers with a red pencil.

ACTION: Qualify both positive results and non-detects for the outlier compound(s) as estimated, "J". When %D is above 90%, qualify all positive results for that analyte "J" and all non-detect results for that analyte "R".

	PA Region od 524.2 (Dat				have a RRF < 0.05? [with a red pencil. fy all positive results and all non-detect lyte "R". con/calculation errors in en initial and continuing t two values but if re) [red pencil. take action as specified VIII Equivalent) areas on the internal of every sample and blank r limits (-50% to + 100%) t calibration and cresponding continuing d 524.2-21, Sect.10.3.4)?	Date: November SOP HW-29, Rev.		
	15.5 Do a	any volat:	ile con	npounds h	ave a RRF	< 0.05?	-			
	ACTION:	Circle a	all out	cliers wi	th a red	pencil.				
	ACTION:	for that	analy	/te "J" a	nd all no		ults			
	the RRF's	reporting	g of %I Check a) between it least [.]	initial two values	and cont				
	ACTION:	Circle 6	errors	with a r	ed pencil	•				
	ACTION:			large, t 2 above.	ake actio	n as spe	cified			
16.0	Internal	Standar	ds (CLI	P Form VI	II Equiva	lent)				
	stan with for (-30 cali The	ndard repondent the upeach initial each init	orting oper ar tial mi 0%) of check mits fo	forms of nd lower dd point the corr (Method or intern	every sate limits (-calibrati esponding 524.2-21, al standa	mple and 50% to + on and continu Sect.10	blank 100%) ing .3.4)?	YES	NO	NA
		the next p		100	oo 0110 d. '		_1_	1		
	ACTION:		, take				ction			
	ACTION:	List ead	ch outl	Lying int	ernal sta	ndard be	low.			
	Sample ID) IS	#	Area Lo	wer Limit		Upper	Limi	it	

(Attach additional sheets if necessary.)

ACTION: 1. If the internal standard area count is outside the upper or lower limit, flag with "J" all positive

results quantitated with this internal standard.

- 2. Do not qualify non-detects when the associated IS Area is above the upper limit (+ 100%).
- 3. If the IS area is below the lower limit (- 50% for initial calibration and -30% for the corresponding continuing calibration), qualify all associated non-detects "UJ".
- 4. If extremely low area counts are reported (< 25%) or if performance exhibits a major abrupt drop off, flag all associated non-detects as unusable "R" and positive results as estimated "J".
- 16.2 Are the retention times of all internal standards within 3 standard deviations of the mean retention compounds in the associated initial mid-point calibration standards, Method 524.2-25, Sect.11.6)?[] _____

ACTION: Professional judgement should be used to qualify data if the retention times differ by

YES NO NA more than 3 standard deviations.

17.0 Field Duplicates

17.1	Were	any	field	duplicates	submitted	for	volatile		
	analysis?							<u>[]</u>	

ACTION: Compare the reported results for field duplicates and calculate the relative percent difference.

ACTION: Any gross variation between field duplicate results must be addressed in the Data Assessment. However, if large differences exist, take action specified in section 3.2 above.

Date: November 2010 Method 524.2 (Rev.4.1, 1995) SOP HW-29, Rev. 2

DEFINITIONS

Acronyms:

bromofluorobenzene

calibration factor (without internal standards)

contract laboratory program

BNA - base neutral acid

CCC - calibration check compound

CF - calibration factor (without

CLP - contract laboratory program

CRQL - contract required quantitat

% D - percent difference or percent contract required quantitation limit percent difference or percent drift GC/MSgas chromatography/mass spectroscopy

IS internal standard

1 liter

LFB laboratory fortified blank LRB laboratory reagent blank LFM -FRB laboratory fortified matrix

field reagent blank

Kq kilograms

meter m

mm – millimeter

m/z mass to charge ratio

QC quality control

RIC reconstructed ion chromatogram relative percent difference RPD -

RRF relative response factor (requires internal standard)

RRT relative retention time

RSD relative standard deviation

retention time RT -

sample delivery group

standard operating procedure

SDG -SOP -SPCC -TIC -TCLP system performance check compound tentatively identified compound

toxicity characteristic leach procedure

micrograms

VOA - volatile organic acid

DEFINITIONS

Data Qualified Definitions:

- U -The analyte was analyzed for, but was not detected above the reported sample quantitation limit.
- J -The analyte was positively identified; the associated numerical value is the approximate concentration of the analyte in the sample.
- N -The analysis indicates the presence of an analyte for which there is presumptive evidence to make a "tentative identification".
- NJ -The analysis indicates the presence of an analyte that has been "tentatively identified" and the associated numerical value represents its approximate concentration.
- UJ -The analyte was not detected above the reported sample quantitation limit. However, the reported quantitation limit is approximate and may or may not represent the actual limit of quantitation necessary to accurately and precisely measure the analyte in the sample.
- R -The sample results are rejected due to serious deficiencies in the ability to analyze the sample and meet quality control criteria. The presence or absence of the analyte cannot be verified.