MATERIALS AND METHODS

Chemicals

Authentic reference standards of Sumithion and Sumioxon were supplied by Sumitomo Chemical Company, LTD, at the following specified purities: Sumithion (Lot No. 40926) -96.6%: Sumioxon was supplied at unspecified purity; Chlorfenvinphos (Lot No. S 502) was obtained from the EPA Analytical Reference Standard Repository (Research Triangle Park, NC 27711) at 91.6% purity, for use as an internal standard. Toluene (HPLC grade) was obtained from Fisher Scientific.

Samples and Sample Storage

Soil samples were received frozen on dry ice from the two field study sites:

- Watson Consulting Services & Enterprises, Inc. Field Research Station 21 W042 Finley Road Downers Grove, IL 60515
- Wayne L. Curry Weed Systems, Inc. Route 3 Box 3115 Melrose, FL 32666

Dissipation study reports from the Illinois and Florida studies are attached as Appendix I and II (Soil and Turf dissipation) and Appendix III (Turf dissipation), respectively. Upon receipt, samples were identified and stored frozen (-20°C) until the time of analysis.

Fortification of Soil Samples with Surithion, Sumioxon, and Chiorfenvinghos

Control soil samples from both study sites were fortified with authentic reference standards of Sumithion, Sumioxon, and Chlorfenvinphos at various concentrations to facilitate the determination of standard curves covering the ranges of each compound's concentration in the actual test samples. The concentration of each standard added to control samples are given in Table 1. Authentic standards were initially dissolved in toluene and volumetrically added to control soil samples.

789000 - PTFL FIDE 1 GF

Extraction of Soil Samples

Soil samples were thoroughly mixed and 5.0 or 10.0 gm (± 0.01 gm) subsamples were removed and placed into 20 ml glass scintillation vials. To each sample was added 200µl of internal standard working solution (50 µg Chlorfenvinphos/ml) followed by 10 ml toluene and extracted on a wrist-action shaker for 40 minutes. The solvent was then removed and centrifuged at 200 x g for 10 minutes to remove particulates. Aliquots (2 to 5 µl) of the supernatant were then injected into the GLC.

Gas-Liquid Chromatography

All sample extracts were analyzed using a Varian 3400 gas-liquid chromatograph equipped with a Model 800 Autosampler, nitrogen-phosphorous selective detector and a Hewlett-Packard 3392 A Integrator. Analyses utilized a 10 meter HP-17 crosslinked 50% Ph Me Silicone column (0.54 mm i.d.; 2.0 µm film thickness). GLC analysis conditions were as follows.

Column Flow - 3.0 ml/min N₂

Injector temperature - 250°C

Detector temperature - 300°C

Temperature Program

initial temp	200~C	held 0 min		
to 🌣	250°C	rate 50/min		
to	275°C	rate 250/min		

Analysis or Chromatographic Data

Standard Curves. Standard curves were determined concurrently with each extracted series of samples (i.e., samples included in a single chromatographic analysis set, consisting of 15 to 45 samples). A representative standard curve is presented in Figure 1.

Test Samples. Concentrations of Sumithion and Sumioxon were calculated using regression analysis comparing concentrations of standards to the relative response of compound/internal standard (chlorfenvinphos).

730999

127,012 150

Table 1. Concentrations of Authentic Reference Standards in Fortified Control Soil Sample Used For Representative Standard Curve.

Matrix		Sunithion		Sumic	xon	•		
SOIL		0,		0			•	
: }	* .	20	n	20	•	•		
∴ . :		40 💬 🧢		40				
		75	,	75	* . '	,		k a
\$		125		125	•	•		1 1
	_	250		250				• •.
	•	375		375	a.	e.	•	,
•	* 1 E	1000		1000		,		

Soil samples weighed 5.00 gm \pm 1.1 or 10.0 gm \pm 0.1.