

1.0 INTRODUCTION

Two analytical methods were used for analysis of daminozide and its residues in soil; one for determining daminozide following derivatization to cyclic daminozide and the other for determining UDMH (unsymmetrical dimethyl hydrazine) residues. A GC equipped with a nitrogen-phosphorus detector was used for determining cyclic daminozide following derivatization of daminozide with acetic anhydride to form cyclic daminozide. The analytical method for UDMH involves distilling UDMH from soil followed by derivatization and analysis by GC/MS. These methods were developed in 1985-1986, prior to the implementation of the GLPs. The analytical method was taken from the soil field dissipation study and therefore does not strictly follow the EU or EPA requirements for analytical methods.

The objective of this report is to provide analytical method(s) for determining daminozide and its residues in soil and address the method specificity, linearity, precision (repeatability) and accuracy/recovery

2.0 MATERIALS AND METHODS

2.1 CHEMICALS

Structure of Daminozide, cyclic Daminozide and UDMH are provided in Figure 1.

2.1.1 Test Substance

The test substance is Daminozide (succinic acid-2,2-dimethylhydrazide).

2.1.2 Solvents and Reagents

- Daminozide analytical standard
- Acetic anhydride
- Sodium bicarbonate
- Hydrochloric acid
- Dichloromethane
- Anhydrous sodium sulfate

2.2 SAFETY AND HEALTH

This method should be performed by trained chemical personnel. Hazards associated with the chemicals used in this analytical method are shown in the MSDS sheets in Appendix 1.

2.3 SUPPLIES

- Gas Chromatograph with Mass Selective Detector
- Module heater
- Circulating bath, Refrigerated
- Grinder
- Vortex junior mixer
- Magnetic stirrer and Hot plate with ceramic Top
- Lab-Jack
- Condenser, Liebig 14/20 JTS
- Distillation flasks
- Adapter, Enlarging
- Adapter, 75 ° Bend
- Adapter, 105 ° Bend
- Receiving Tube
- Syringe
- Vials
- Volumetric pipets
- Volumetric flasks
- Graduated Cylinders
- Stoppers-Glass
- Magnetic stirring bars
- Rubber tubing
- Sodium hydroxide, 50% Aqueous
- Titanium Trichloride, 20% Aqueous
- Calcium Chloride
- Daminozide (Succinic acid 2,2-dimethylhydrazide)
- Salicylaldehyde
- Sodium chloride
- Methylene Chloride

2.4 TEST MATRIX

2.4.1 Soil

This method should be applicable to most soil types. In Uniroyal Chemical Inc., Project 8647, soil obtained from Uniroyal Crop Protection field station, Bethany, CT was used. The soil was characterized as Fine sandy loam soil (Spodosol to Alfisol/gray-brown Podzolic soil). The soil characterization is as follows:

pH	Sand	Silt	Clay	Organic Matter
6.1	54%	34%	10%	3.3%

2.4.2 Sampling

Samples were taken from three random locations in the plot. Soil around the immediate sampling site was pushed away using a shovel and an 8x8x6 inch sample was taken for the 0-6 inch level. A 6-12 inch sample was then taken using a shovel from the same hole. Each soil sample was mixed thoroughly on a plastic sheet.

2.4.3 Soil Fortification

The following soil samples were spiked with daminozide:

<u>Soil Depth</u>	<u>Amount of Soil</u>	<u>Fortification level (ppm)</u>	<u>Residue Analyzed</u>	<u>Analytical Method</u>
0-6 inch	20 g	1 ppm	Cyclic daminozide	GC-NPD
6-12 inch	20 g	1 ppm	Cyclic daminozide	GC-NPD
0-6 inch	50 g	0.4 ppm	UDMH Residues	GC/Mass Spec.
6-12 inch	50 g	0.4 ppm	UDMH Residues	GC/Mass Spec.
0-6 inch	10 g	2 ppm	UDMH Residues	GC/Mass Spec.
6-12 inch	10 g	2 ppm	UDMH Residues	GC/Mass Spec.

2.4.4 Soil Extraction and Analysis

The laboratory spiked soil samples along with the field treated soil samples were extracted as described in section 2.6 and analyzed as described in section 2.5.

2.5 EQUIPMENT

2.5.1 Gas Chromatograph-Nitrogen Phosphorous Detector

The following GC conditions were used for analyzing Cyclic Daminozide

Instrument: Gas Chromatograph
Detector: Nitrogen-Phosphorus Detector (NPD)
Column: 1.8 m x 2 mm (i.d.) coiled glass column packed with 5% SP-100 (Supelco Inc.), or equivalent on Chromosorb W 100/120.

Gas Flow rate:
Hydrogen: 3 mL/min
Air: 100 mL/min
Helium (Carrier gas): 40 mL/min

Injector Temperature: 230°C

Detector Temperature 300°C

Oven Temperature:
Initial Temperature: 160°C for 8 min.
Final Temperature: 180°C for 7 min.

2.5.2 Gas Chromatograph – Mass Spectrometry

The following GC-MS conditions were used for analyzing UDMH Residues

Instrument: Gas Chromatograph
Column: 30 m x 0.32 mm DB-5
Carrier gas: Helium @ 4 psi head pressure
Injector Temperature: 225°C
Detector Temperature 300°C

Uniroyal Chemical Company, Inc.
Uniroyal Project No.: 99001
Analytical Method No.: AC-6004

Oven Temperature:
Initial Temperature: 145°C (no initial hold time)
Final Temperature: 190°C (no final hold time)
Rate: 4.0°C/min.

Injection: Split mode, split ratio = 7.5:1

Injection Volume: 2µL

Mass Spectrometer Conditions:

Scan:	Multiple Interval Detection (MID)
Mode:	EI (GC)
Manifold Temp:	92°C
Ion Source Temp:	150°C
High Vacuum:	3.2×10^{-6} Torr
Separator Oven Temp.:	208°C
Transfer line Temp.:	213°C
Filament Current: -	0.3 mA
Electron energy:	70.0 eV.
Electron Multiplier Voltage:	-1025 V
Calibrated in EI mode using FC-43 reference standard	

Instrument Operation

- The instrument must be tuned daily.
- Prepare data acquisition parameters to perform selected ion monitoring of m/e 44 (NC₂H₆) and m/e 164 (molecular ion of salicylaldehyde dimethylhydrazone). Quantification is performed on the m/e 164 ion chromatogram and the m/e 44 ion is used as a "confirmation" ion (to make sure the ratio of the m/e 44 to 164 is the same for the standards and the samples).
- Rapidly inject 2.0 µL of the methylene chloride extract into the chromatographic column and immediately activate the data system.
- After the 10 minute analysis time, measure the area under the peak for the hydrazone (RT ~6 min). Compare that peak area with the standard curve generated that day.

The GC/MS with a switchable EI/CI source was operated under the INCOS (ver: 4.07) operating system in the EI mode using capillary GC conditions. The instrument was tuned daily. The mass spectra were obtained using MID (Multiple Interval Detection) scanning program for mass 164 (molecular ion of salicylaldehyde dimethylhydrazone).

The raw ion chromatogram (MID RIC) for each sample was printed for the scan range of 800 to 1100. Typical scan number for the apex of desired peak was 980. The data was then enhanced using the Finnigan MSDS program "NHA" from scan 800 to 1100, covering masses 162-166 m/z, using no smoothing. Background was subtracted using the programs default parameters. The mass list for the desired peak was printed. The intensity of the 164 m/z ion from the mass list was used for quantitation purposes (corresponds to peak height).

2.6 ANALYTICAL METHODOLOGY

2.6.1 Preparation of Calibration Standards for Standard Curve

1. Weigh 0.1015 g of daminozide into a 100 mL volumetric flask. Dilute to the mark with distilled water. Stir for 30 minutes using a magnetic stir bar to ensure complete dissolution. This 1015 $\mu\text{g/mL}$ stock solution may be kept for two days in a 4°C refrigerator.
2. Pipet 1.0 mL of the 1015 $\mu\text{g/mL}$ stock solution into a 100 mL volumetric flask and dilute to the mark with distilled water. This 10.15 $\mu\text{g/mL}$ solution must be prepared daily.
3. Pipet 5.0 mL of the 1015 $\mu\text{g/mL}$ solution into a 100 mL volumetric flask and dilute to the mark with distilled water. This 50.45 $\mu\text{g/mL}$ solution must be prepared daily.
4. Distill four daminozide standards ranging from 0 to 100 μg from 160 mL of 50% NaOH containing 5 mL of 20% TiCl_3 solution and 20 g of CaCl_2 . Follow steps 2 to 7 of the distillation section. 2.6.2.

Note: When TiCl_3 is added to the flask, an exothermic reaction occurs. Wait for the bubbling to subside before spiking the daminozide standard into the flask.

Note: The amounts selected for the standard curve should bracket the amount of daminozide that you expect in 10 to 50 g of sample.

$$\text{e.g.: } \frac{10 \mu\text{g daminozide}}{10 \text{ g sample}} = 1 \text{ ppm}$$

5. Prepare a standard curve with "peak area" on the Y-axis and μg of daminozide" spiked into the flask on the X-axis. The standard curve should be linear.

6. The standard should be distilled every day and the daily slopes recorded. Once the analyst feels comfortable that the data analysis and the distillation techniques can be reproduced in that laboratory, then single standards can be distilled daily to check the response.

2.6.2 Distillation Procedure

1. Add the following reagents to a 100 mL distillation flask containing a teflon stir bar:

10 to 50 \pm 0.01 g Homogeneous sample
20 g Calcium chloride
160 mL NaOH (50% aqueous solution)
5 mL* Titanium Trichloride (20% solution) – reducing agent

*Add TiCl_3 just before assembling distillation apparatus.

2. Place the flask on a “cold” hot plate.

Note: A “cold” hot plate is used to prevent a loss of UDMH that may be hydrolyzed upon the first application of heat. If repetitive distillations are being performed, a second “cold” hot plate should be used in order to give the “hot” hot plate a chance to cool off. Also, ceramic hot plates cool off faster than metal hot plates.

3. Add 20 μL of neat salicylaldehyde to the empty receiving tube prior to distillation.

4. Assemble the distillation apparatus using the following items:
 - Flask
 - Reducing adapter
 - Elbow – 75° bend
 - Liebig condenser – cooled to 0°C with a circulating 50/50 mixture of methanol/ H_2O .
 - Elbow with vent port - 105° bend.
 - Receiving tube – 25 mL graduated tube with 14/20 joint.

Use liberal amount of stopcock grease on all of the joints except the joint to the receiving vessel.

Recommended using clips at the ground glass joints to prevent separation and leakage during distillation.

5. Carefully heat the contents of the flask while stirring rapidly (low heat and stirring will help control the foaming).
6. After 10 mL of distillate is collected, turn off the hot plate. Remove the receiving tube, stopper and place in the 50°C module heater for 90 minutes to allow the derivatization to go to completion.
7. After the 90-minute derivatization step, the receiving tube is allowed to cool. The contents of the receiving tube are quantitatively transferred, using 2 mL of water, to a 20-mL glass scintillation vial containing 3.5 g of NaCl. The salt saturated aqueous phase is then extracted with 4.0 mL methylene chloride by placing the vial on a vortex mixer for 1 minute. The contents were mixed on a vortex mixer for 30 seconds. The phases were allowed to separate and 2.0 µL aliquots of the of the methylene chloride phase (lower layer) were injected into the GC.

The data was acquired using the scan program noted above. The raw ion chromatogram (MID RIC) for each sample was printed for the scan range of 800 to 1100 (typical scan number for the apex of desired peak was 980). The data was then enhanced using the Finigan MSDS program "ENHA" from 800 to 1100, covering masses 162-166 m/z, using no smoothing. Background was subtracted using the program default parameters. The mass list for the desired peak was printed. The intensity of 164 m/z ion from the mass list was used for quantitation purposes (corresponds to peak height).

2.6.3 Preparation of Derivatized Standard (Cyclic Daminozide) for Fortification

A stock solution was prepared by weighing 0.1015 g daminozide and dissolving it in 100 mL water. Aliquots of the stock solution were diluted 1:100 and 5:100 in water to obtain 10.15 µg/mL and 50.75 µg/mL, respectively.

Two milliliters of aliquots of each of the above two dilutions were used to spike 20 g of an aqueous blank (standard curve) or control soil samples, 0-6" and 6-12" depths for fortified sample recovery.

For the standard curve, 15 mL of the 0.01N HCl was added to the standard. All of the sample was derivatized with acetic anhydride. For the soil samples, treated and spiked, 30 mL of 0.01 N HCl was added to the soil. Following a 30 minute extraction, 15 ml (one-half of the extractant volume) was derivatized and assayed. A factor of 2 was used to calculate the recovery and test sample ppm levels measured against the standard curve.

2.7 SAMPLE BRACKETING

Calibration curves using daminozide and UDMH standards were generated using a standard bracketing technique. Four concentration (0, 10.15, 20.30 and 101.5 µg) of daminozide and five different concentrations (0, 5, 10, 30 and 50 µg) of UDMH were used for generating the standard curves. A linear regression analysis was done and the correlation coefficient 'R²' was determined. A second calibration curve was generated at day 7 (168 hours) following the first one and was used for quantification of the 7 day treated samples. The results are provided in Table II/ and Figure 5.

2.8 POTENTIAL INTERFERENCES

This method could have interferences from some soils and was observed with soil samples from different locations. In such cases one should modify the method. In the field dissipation study where the treated soil samples were extracted and analyzed concurrently with the spiked soil samples, the capillary column used was changed from the OV-1701 liquid phase to an equivalent SE-52 phase.

2.9 TIME REQUIRED FOR ANALYSIS

Following the above procedure, in one day, a reference sample (to check reproducibility) along with test samples and analytical standards for generating a calibration curve can be done.

2.10 DATA CALCULATIONS

The peak area for the sample was compared to the standard curve. The amount of daminozide (measured as cyclic daminozide) or UDMH residues in a sample was calculated using regression analysis as follows:

$$y = mx + b$$

where:

y = peak area
m = slope from standard curve
b = y-intercept
x = μg of analyte

The concentration of the analyte in a given sample is calculated as follows:

$$\frac{\mu\text{g Daminozide}}{\text{g Sample}} \times \frac{100}{\% \text{ Recovery}} = \text{ppm Daminozide}$$

Based on duplicate (in most cases) injections of standards and samples, the peak area (or peak height) of the analyte of interest, was determined.

Percent Recovery (% R) is calculated as follows:

$$\% R = [(C_F - C_U)/C_A] \times 100$$

where:

C_F = Concentration of analyte measured in fortified sample.
 C_U = Concentration of analyte measured in unfortified sample.
 C_A = Concentration of analyte added to fortified sample.

2.11 METHOD TESTING

2.11.1 Specificity

Two different analytical methods were used for the separation of daminozide and UDMH. The GC-NPD method used for analysis of daminozide is specific in separating it from any interfering peaks. The GC-MS method used for UDMH is specific in not only separating it from interfering peaks but the molecular ion of UDMH has been confirmed by Mass Selective detection.

2.11.2 Linearity

The linearity of the test procedure is the ability to obtain test results proportional to the concentration of analyte within the sample. Linearity data were obtained for daminozide using a four-point calibration curve to get a correlation coefficient. For UDMH, linearity data were obtained

using a five-point calibration curve to get a correlation coefficient. Each concentration containing the analyte was injected in duplicate and the average value was used to determine the linearity of the method. Out of the four to five concentrations used, the concentration of the analyte of interest was in the middle with one to two concentrations above and below the fortification level selected to ensure the linearity of the method.

2.11.3 Precision

Precision or Repeatability is the closeness of agreement between independent test results obtained with the same method on identical test material in the same laboratory by the same operator using the same equipment within short intervals of time.

The precision of the method was determined by making duplicate injections of the 0-6 inch and 6-12 inch control soil spiked with daminozide at the 1ppb level on two different days.

2.11.4 Accuracy/Recovery

The accuracy/recovery of the method is the degree to which the observed results correspond to the true value of the analyte in the sample.

Two control untreated soil (0-6 inch and 6-12 inch) samples were each fortified with 1 ppm daminozide .

Each of the fortified soil samples were extracted and analyzed in duplicate as described in Sections 2.5 and 2.6. The recovery of the method was determined based on the fortified soil samples.

The samples were analyzed in the following order:

- Calibration standards for standard curve
- Control soil (0-6 inch).
- Control soil (0-6 inch) spiked with 1 ppm daminozide.
- 0 Time point treated soil (0-6 inch)
- Control soil (6-12 inch)
- Control soil (0-6 inch) spiked with 1 ppm daminozide.
- 0 Time point treated soil (6-12 inch)
- 3-Hour Time point treated soil (0-6 inch)
- 6-Hour Time point treated soil (0-6 inch), etc.