

ENVIRONMENTAL CHEMISTRY METHOD FOR THE DETERMINATION OF DPX-KN128 RESIDUES IN SOIL USING GC-MSD

Gary L. Westberg & Daniel R. Vincent

REASON FOR REVISION

This method was revised to remove the analysis of IN-KG433. Although this analyte was observed in laboratory soil metabolism studies, it has not been observed in field dissipation studies. Therefore, it is not necessary to include it as an analyte in soil analyses for either future soil studies or for environmental testing.

1.0 ABSTRACT

DPX-MP062 and DPX-JW062 Experimental Insecticides are formulated with an active ingredient that is a mixture of two enantiomers: DPX-KN128 and IN-KN127; DPX-KN128 is insecticidally active, while IN-KN127 is not. DPX-MP062 is comprised of approximately 3 parts of DPX-KN128 and approximately 1 part of IN-KN127. DPX-JW062 is a racemic mixture (50% of each) of the two enantiomers.

An environmental chemistry method has been established to determine residues of the insecticide DPX-KN128/IN-KN127 (parent) and its known metabolite DPX-JT333 in soil; enantiomers are not resolved in the analysis for either analyte. Soil samples are twice-extracted in acetonitrile (ACN) and water (9:1, v:v), with mixing energy supplied by a wrist-action shaker. Two phases are developed in the extract by adding solid sodium chloride. An aliquot of the ACN layer is removed from the salt-separated extract in order to analyze for parent and DPX-JT333. After evaporating the ACN, the sample is subjected to liquid:liquid partitioning and silica solid-phase extraction (SPE). Final extract is analyzed by single ion monitoring (SIM) GC-MSD. Parent and DPX-JT333 molecular ions are monitored in the analysis. Standards are prepared in control matrix extract to overcome observed matrix enhancement. The method has an apparent Limit of Quantitation (LOQ) of 0.01 ppm of parent equivalents from either analyte.

2.0 INTRODUCTION

DPX-MP062 and DPX-JW062 Experimental Insecticides are being developed for control of various lepidopterous and other insect pests in cotton, sweet corn, pome fruits, vegetables, grapes, and root crops. The intended crop protection product will be used on a wide variety of crops, and the use will likely be multiple applications per season.


This environmental chemistry method is intended for the determination of DPX-KN128/IN-KN127 (parent) and its principal metabolite DPX-JT333 in soil samples. The principal use of the method has been to collect residue data for submission to regulatory authorities. Besides generating residue data, the method is intended to be used for regulatory surveillance.

DPX-KN128 is insecticidal, while its enantiomer, IN-KN127, is not. Since this chemical has a chiral center, the technical material exists as either a racemic mixture of both DPX-KN128 and IN-KN127 (DPX-JW062) or as an enhanced mixture containing approximately 3 parts of DPX-KN128 and 1 part of IN-KN127 (DPX-MP062). Known, significant, soil analytes include parent and DPX-JT333.

The analytical method does not resolve the enantiomers of either analyte. Concentrations of each analyte are measured, calculated, and reported as equivalents of either DPX-JW062 or DPX-MP062. In this method, parent refers to the starting reagent, whether it is DPX-JW062 or DPX-MP062; weights and concentrations are not adjusted for enantiomeric content. Parent and DPX-JT333 are detected as the particular molecular ion.

Structures, names and limited chemical data for the two analytes follow.

DPX-KN128




CAS Name

(S)-methyl 7-chloro-2,5-dihydro-2-[[[(methoxycarbonyl)[4-(trifluoromethoxy)phenyl]amino]carbonyl]indeno[1,2-e][1,3,4]oxadiazine-4a(3H)-carboxylate

CAS Registry Number 173584-44-6
Molecular Weight 527.84
Molecular Composition C₂₂ H₁₇ Cl F₃ N₃ O₇

DPX-JT333



CAS Name methyl 7-chloro-2,5-dihydro-2-[[[4-(trifluoromethoxy)phenyl]amino]carbonyl]indeno[1,2-e][1,3,4]oxadiazine-4a(3H)-carboxylate
CAS Registry Number 144171-39-1
Molecular Weight 469.81
Molecular Composition C₂₀ H₁₅ Cl F₃ N₃ O₅

3.0 MATERIALS

Equivalent equipment and materials may be substituted unless otherwise specified; note any specifications in the following descriptions before making substitutions. Substitutions should only be made if equivalency/suitability has been verified with acceptable control and fortification recovery data.

3.1 Equipment

Assorted laboratory glassware including: graduated cylinders, short stem glass funnels, pipettes, volumetric flasks, evaporating flasks, microliter syringes are needed for this method.

Gas Chromatograph: Hewlett-Packard 5890E gas chromatograph equipped with a Hewlett-Packard 5972 Mass Selective detector, an HP7673 Autosampler, and an HP G1034C MS ChemStation

GC Column: 10 m × 0.20 mm i.d. fused silica column crossbonded with a 0.33- μ m film of DB-1 (J&W Scientific, Folsom, CA)

Balances: Analytical balance capable of weighing to 0.1 mg for weighing analytical standards

Top-loading balance capable of weighing to 0.01 g, for all other weighing

Centrifuge:	IEC Model HN-SII (Damon IEC Division, Needham Hts., MA)
Centrifuge Bottles:	250 mL, Nalgene® polyethylene bottles, with screw caps
Centrifuge Tubes:	50 mL, Pyrex® conical tubes, silanized, with Teflon®-lined screw caps
Evaporators:	N-Evap Laboratory Sample Evaporator Model 115 attached to a nitrogen source (Organomation Associates, South Berlin, MA)
Extraction Apparatus:	Wrist action shaker (Burrell Corporation, Pittsburgh, PA)
GC autosampler vials:	2.0 mL, glass, with 400- μ L glass inserts (Alltech, Deerfield, IL)
Mixer:	Vortex Genie 2 (VWR Scientific, Bridgeport, NJ)
Separatory Funnels:	125 mL, Kimax, silanized, with Teflon® stopcocks
Solid-Phase Extraction Apparatus:	Vac Elut Model SPS 24 (Varian Analytical Instruments, Sunnyvale, CA)
Solid-Phase Extraction Columns:	Silica Bond Elut® Extraction Column, Part #1225-6026. 20 cc/5 g silica sorbent (Varian Analytical Instruments, Sunnyvale, CA) <u>Do not substitute.</u>
Syringe Needles:	4-inch, 13-gauge, stainless steel, Luer-Lok fitting (Popper & Sons, Inc., New Hyde Park, NY)
Syringes:	5 mL, glass, B&D Multifit, (VWR Scientific, Bridgeport, NJ)
Test Tubes:	13 \times 100 mm, borosilicate, silanized, with Teflon®-lined screw caps. 16 \times 150 mm, borosilicate, silanized
Assorted laboratory glassware	

3.2 *Reagents and Standards*

Acetonitrile:	Pesticide residue quality
Dimethyldichlorosilane:	Supelco, Catalog No. 3-3009 (Supelco, Inc., Bellefonte, PA)

Ethyl Acetate:	Pesticide residue quality
Hexane:	Pesticide residue quality
Isopropanol:	Pesticide residue quality
Toluene:	Pesticide residue quality
Sodium Chloride:	Reagent grade
Water:	HPLC grade
Silica wool:	Fused, catalog #20790 (Restek Corp., Bellefonte, PA)
DPX-KN128	Analytical standard grade DPX-KN128 (as either DPX-JW062 or DPX-MP062, available from DuPont Agricultural Products, Global Technology Division (E. I. du Pont de Nemours and Company, Wilmington, DE). The product used in this work was DPX-MP062-38, which had a purity of 97.1%.
DPX-JT333	Analytical standard grade DPX-JT333, available from DuPont Agricultural Products, Global Technology Division (E. I. du Pont de Nemours and Company, Wilmington, DE). The product used in this work was DPX-JT333-17, which had a purity of 98.3%.

3.3 *Safety and Health*

Each analyst must be acquainted with the potential hazards of the reagents, products and solvents used in this method before commencing laboratory work. All appropriate material safety data sheets should be read and followed, and proper personal protective equipment should be used.

4.0 METHODS

4.1 *Principle of the Analytical Method*

The soil sample is extracted twice with two volumes (2 mL/g of soil) of acetonitrile/water (9/1, v/v). The water and ACN in the extract are separated into two phases by adding solid sodium chloride. An aliquot of acetonitrile extract is then processed through additional analyte-specific cleanup procedures, including a hexane partition, a water/hexane partition, and a solid phase extraction (SPE) purification incorporating silica. Analytes are detected and quantified by gas chromatography using mass selective detection (MSD). The limit of quantitation for both analytes is 0.01 ppm.

4.2 *Analytical Procedure*

4.2.1 Glassware Cleaning and Preparation Procedures

The effectiveness of any cleaning procedure used should be demonstrated by preparation and analysis of reagent blanks. In general, all reusable glassware and plastic-ware should be washed in hot tap water with laboratory grade, non-phosphate detergent, rinsed several times with tap water, rinsed several times with de-ionized water, rinsed once with acetone, and allowed to dry before use. Care should be taken to avoid working with high levels of the analyte being monitored in the same laboratory where samples are being extracted and analyzed.

To avoid and/or reduce incidental adsorption of analytes to glassware, it is recommended that all glassware, except syringes, that contact sample extract be silanized. See Appendix 1 for detailed instructions on silanizing glassware.

4.2.2 Preparation of Reagent Solutions

2% Isopropanol/98% Hexane (v/v)

Add 20 mL isopropanol to 980-mL hexane. Mix.

5% Isopropanol/95% Hexane (v/v)

Add 50 mL isopropanol to 950-mL hexane. Mix.

5% Ethyl Acetate/95% Toluene (v/v)

Add 25 mL ethyl acetate to 475-mL toluene. Mix.

50% Ethyl Acetate/50% Hexane (v/v)

Add 500 mL ethyl acetate to 500-mL hexane. Mix.

10% Water/90% Acetonitrile (v/v)

Add 100 mL water to 900-mL acetonitrile. Mix.

4.2.3 Stock Standard Preparation and Stability

Weigh 12.5 mg (corrected for purity) of analytical standard and quantitatively transfer it to a 25-mL volumetric flask; separate solutions should be prepared for each analyte. Dissolve the standard in and bring to volume with ethyl acetate. The final concentration of the individual stock standard solutions is 500 µg/mL. Store the stock standard solutions at 1-8°C when not in use. Maximum storage life for the solution is 16 months for parent and 6 months for DPX-JT333.

4.2.4 Fortification Standard Preparation and Storage

Depending on the experimental design, prepare either separate or combined parent or DPX-JT333 standard solutions. Solutions at the following concentrations should be prepared every month.

- 100 µg/mL: Transfer 5.0 mL of 500-µg/mL stock standard solution(s) to a 25-mL volumetric flask. Bring to volume with ethyl acetate. Mix well. Store tightly capped at 1-8°C.
- 10 µg/mL: Transfer 2.5 mL of 100-µg/mL standard solution to a 25-mL volumetric flask. Bring to volume with ethyl acetate. Mix well. Store tightly capped at 1-8°C.
- 1 µg/mL: Transfer 2.5 mL of 10-µg/mL standard solution to a 25-mL volumetric flask. Bring to volume with ethyl acetate. Mix well. Store tightly capped at 1-8°C.

4.2.5 Chromatographic Standard Preparation and Storage

Prepare calibration standards in ethyl acetate. The following concentrations are needed:

- 0.5 µg/mL: Transfer 1.25 mL of the 10-µg/mL standard solution of each analyte to a 25-mL volumetric flask. Bring to volume with ethyl acetate and mix well. Store tightly capped at 1-8°C.
- 0.05 µg/mL: Transfer 125 µL of the 10-µg/mL standard solution of each analyte to a 25-mL volumetric flask. Bring to volume with ethyl acetate and mix well. Store tightly capped at 1-8°C.

These solutions serve as the source for preparation of all actual GC standard solutions.

Prepare the following concentrations of either separate or combined parent and DPX-JT333 standard solutions, as applicable, in control soil extract:

- 0.1 µg/mL: Transfer 40 µL of 0.5-µg/mL analyte standard solution to a 400-µL glass insert contained within a glass autosampler vial. Evaporate solvent under a stream of nitrogen at room temperature. Add 200 µL of control sample extract from Step 4.3.4, cap the vial, and mix by vortex.
- 0.05 µg/mL: Transfer 20 µL of 0.5-µg/mL analyte standard solution to a 400-µL glass insert contained within a glass autosampler vial. Evaporate solvent under a stream of nitrogen at room temperature. Add 200 µL of control sample extract from Step 4.3.4, cap the vial, and mix by vortex.
- 0.025 µg/mL: Transfer 10 µL of 0.5-µg/mL analyte standard solution to a 400-µL glass insert contained within a glass autosampler vial. Evaporate solvent under a stream of nitrogen at room temperature. Add 200 µL of control sample extract from Step 4.3.4, cap the vial, and mix by vortex.
- 0.006 µg/mL: Transfer 24 µL of 0.05-µg/mL analyte standard solution to a 400-µL glass insert contained within a glass autosampler vial. Evaporate solvent under a stream of nitrogen at room temperature. Add 200 µL of control sample extract from Step 4.3.4, cap the vial, and mix by vortex.

Store all GC standard solutions between 1°C and 8°C when not in use.

4.2.6 Source (and Characterization) of Samples

Since the method requires the use of control sample extract for preparing calibration standards, soil samples should be collected from both treated and untreated locations.

Samples used to validate this method were collected from two DuPont Agricultural Products sites in the United States (Bradenton, Florida, and Madera, California) during a dissipation study (AMR 3402-95) and from test site(s) in Canada, also during a dissipation study (CAN-96-092). Samples were shipped frozen to Morse Laboratories, Inc., in Sacramento, California, where analysis was performed. Soil from the following soil segments were tested: 0 to 15 cm, and 30 to 45 cm.

4.2.7 Storage and Preparation of Samples

Soil samples are stored frozen (-10°C or less) until analysis. Before chemical analysis, each sample should be homogenized in a mixer, and stones (i.e., greater than approximately 5-mm diameter) and debris (e.g., wood, plant material, etc.) should be removed.

4.2.8 Sample Fortification Procedure

Once the sample is weighed into the extraction vessel (bottle), an appropriate volume (≤ 1.0 mL) of fortification standard solution should be added to the soil surface. Allow the solvent to evaporate for approximately 5 min before proceeding.

4.2.9 Analyte Extraction Procedure

1. Weigh 10.0 g of soil sample into a 250-mL polyethylene centrifuge bottle. **Fortify appropriate samples at this time** (see Section 4.2.8).
2. Add 20.0 mL of extraction solution (10% water/90% acetonitrile) and cap the bottle.
3. Shake the sample on a wrist action shaker for 30 min.
4. Centrifuge for 15 min. at approximately 2400 rpm.
5. Carefully decant the supernatant into a 125-mL, silanized, separatory funnel.
6. Repeat Steps 2, 3, and 4, making sure the soil pellet formed at Step 4 is broken up prior to extraction.
7. Carefully decant the supernatant into the 125-mL separatory funnel containing the supernatant from the first extraction.
8. Add 3.0 g of sodium chloride to the combined extract in the separatory funnel and gently shake for approximately 30 sec. This step forces a separation of the water from the ACN. Allow the layers to separate (approximately 30 min.). Using a glass stirring rod, gently mix the upper acetonitrile layer to insure homogeneity.

9. Transfer a 3.6-mL aliquot of the acetonitrile layer to a silanized, 50-mL centrifuge tube. Label appropriately.

THE ANALYSIS CAN BE STOPPED AT THIS POINT. CAPPED EXTRACTS SHOULD BE STORED FROZEN (-8°C OR LOWER).

4.2.10 *Analyte Purification Procedure*

Check or calibrate the SPE cartridges prior to use in order to ensure optimum method performance. In general, check one cartridge per lot number. This assessment should be conducted well in advance of needing the columns for sample analysis. Analyte recovery greater than 90% indicates that a box of cartridges is suitable for use. Assessment analyses are conducted on a reagent blank (no soil). See Appendix 2 for detailed instructions on assessment of SPE cartridge.

4.2.10.1 *Hexane Washing of ACN extract aliquots*

1. Add 3.6 mL of hexane to each centrifuge tube containing the 3.6 mL of acetonitrile extract, cap and shake the tube for approximately 1 min. Remove and discard the upper hexane layer using a long-stemmed Pasteur pipette.
2. Repeat Step 1.
3. Evaporate the ACN remaining in each tube to near dryness using an N-Evap at 50°C. Gently blow dry with nitrogen at room temperature. Proceed to the next section (Section 4.2.10.2).

4.2.10.2 *Parent & DPX-JT333 Analysis*

Note: SPE cartridges are operated by attaching them to the Vac-Elut apparatus, setting the vacuum so that the effluent forms a steady flow of distinct droplets.

1. Add 2.0 mL of water to the tube containing the evaporated ACN extract (Step 3 of Section 4.2.10.1).
2. Vortex mix each sample (in 2.0-mL water) for approximately 30 sec.
3. Add 10 mL of hexane and mix thoroughly (hand shake) for approximately 30 sec.
4. Draw off the upper hexane layer with a long-stemmed Pasteur pipette and place in another silanized, 50-mL centrifuge tube.
5. Repeat Steps 3 and 4 placing the hexane layer into the 50-mL centrifuge tube containing the hexane layer from the first extraction.
6. Concentrate the combined hexane extracts to approximately 4.0 mL using the N-Evap (water bath set at 40°C).
7. Condition a 5-gram Silica Bond Elut® cartridge by passing 40 mL of hexane through the cartridge. Do not let the cartridge go to dryness after conditioning. (Stop elution when hexane reaches top of frit.) Discard conditioning solvent.

8. After conditioning, pass the sample extract from step 6 through the silica cartridge. Wash the centrifuge tube with two, 1-mL rinses of hexane, passing these through the cartridge as well. Do not let the cartridge go to dryness. Discard eluate.
9. Wash the cartridge with 10 mL of hexane. Discard this wash.
10. Wash the cartridge with 30 mL of 2% isopropanol/98% hexane. Discard this wash.
11. Add 35 mL of 5% isopropanol/95% hexane elution solvent. Discard the first 15 mL of this solvent.
12. Collect the remaining eluate in a silanized 16 × 150-mm test tube and reduce the volume to approximately 0.5 mL on the N-Evap (water bath set at 40°C).
13. Quantitatively transfer the sample to a silanized 13 × 100-mm test tube with small amounts of ethyl acetate and continue the evaporation using the N-Evap at 40°C to approximately 0.2 mL.
14. Evaporate the organic solvent with a stream of nitrogen at room temperature, and re-dissolve the sample residue in 1.0 mL of ethyl acetate. This solution represents 1.0 g of soil sample. Submit to GC analysis.

4.2.11 Derivatization Procedure

No derivatization is necessary in this method.

4.3 Instrumentation

4.3.1 Description

The method was developed and validated on a Hewlett-Packard model HP5890E gas chromatograph equipped with a HP5972 mass selective detector and a HP7673 autosampler, and an HP G1034C MS ChemStation. Instrument conditions follow.

4.3.2 Instrument Operating Conditions

Column:	10-m × 0.20-mm i.d. fused silica column crossbonded with a 0.33- μ m film of DB-1
Inlet Liner:	2 mm i.d. gooseneck splitless liner lightly packed with fused silica wool
Injection Volume:	4 μ L
Carrier Gas:	helium
Column Head Pressure:	12.8 psi
Purge Flow Timing:	on at 1.00 min.
Temperatures:	Injector: 300 °C GC/MSD transfer line: 295°C Column: Initial: 240°C Rate: 20°C/min.

Final: 295°C, hold for 2.75 min

Tuning: Prior to analysis, the instrument is tuned manually for ion m/e 502.

Ions Monitored: Parent m/e 527 Enter this ion 3 times into Group 1 of the SIM acquisition table.

DPX-JT333 m/e 469 Enter this ion 2 times into Group 2 of the SIM acquisition table. This arrangement will improve sensitivity for the instrument.

Dwell Time: 50 msec

Retention Times: approximately 3.0 minutes for Parent
approximately 3.5 minutes for DPX-JT333

4.3.3 Calibration Procedures

Prepare a four-point standard curve by injecting analyte standard solutions prepared in control-sample extract. Preparation of additional control-sample extract in each analytical set will be necessary to prepare GC standard. Guidance on the amount of control-sample extract is provided in the following section. All injections (analyses) of standard solutions are used to construct the standard curve.

4.3.4 Sample Analysis

This method relies on the preparation of calibration samples in control soil extract to offset matrix enhancement. Matrix enhancement is a phenomenon often observed in GC-MS analytical methods in which the signal of analytes in matrix extract is considerably higher than in pure solvents. The outcome of matrix enhancement is that recovery of analyte in fortified samples is consistently very high (e.g., 150% or higher). Several attempts (dilution, clean-up, injection additives, etc.) were made to eliminate enhancement, but were unsuccessful. Therefore, control matrix extract is required to "control" the enhancement phenomenon.

Inject a curve check standard after analyzing every 4 or 5 samples.

If the final-sample extract requires dilution for the response to be within the range of standard curve responses, then the final-sample extract should be diluted with control-sample extract and re-injected. The table below describes a typical analytical set and the corresponding amount of control-sample extract required. Experience with parent and DPX-JT333 analysis in soil has shown that final extracts of fortification samples typically require one dilution, while those of field samples may require up to three dilutions. It is unlikely that surveillance sample extracts would require dilution.

<u>Sample Set Description</u>	<u>Volume of control sample extract required</u>
GC standard solutions: (4 solutions × 200 µL/solution)200 µL/solution)	800 µL
1 control sample	0 µL
1 fortification requiring no dilution	0 µL
1 fortification requiring dilution	200 µL
8 samples, 6 of which are likely to require dilution (6 samples × 3 dilutions/sample × 200-µL/dilution)	<u>3600 µL</u>
Total	4600 µL

Starting at Step 9 of Section 4.2.9, process five additional 3.6-mL aliquots of control extract through the remainder of the procedure. This will yield five 1-mL portions (a total of 5 mL) of control sample extract. Combine all extracts generated into a single solution prior to use.

4.4 Calculations

4.4.1 Methods

Calculation of residue concentrations from instrumental data were conducted using a validated software application. Standard least-squares regression was used to calculate a best-fit line based on known concentrations (standards solutions) and the resultant peak area from the analysis of those solutions. The line was then used to determine concentrations of the analyte found during sample analysis.

The equation used for the calibration line was

$$y = mx + b$$

where: y = net peak area response (response in fortification *less* response in control),

x = µg analyte/mL,

m = slope of the regression line, and

b = y-intercept.

Therefore, the µg analyte/mL value for each sample was calculated according to the following adaptation of the regression equation:

$$x = (y - b) / m$$

The concentration (ppm) of each analyte in a residue sample was calculated according to the following equation.

$$\text{ppm of analyte} = \frac{\mu\text{g analyte} / \text{mL} \times \text{mL solvent} \times \text{mL of final volume} \times \text{GC dilution factor}}{\text{g of sample} \times \text{mL of aliquot}}$$

where, $\mu\text{g analyte/mL}$ was derived in the regression line from the area of the analyte peak,
 mL solvent was the volume of ACN portion of the combined extraction solvent,
 mL of aliquot was the volume of ACN extract taken through cleanup steps,
 g of sample was the weight of the portion of sample extracted,
 $\text{mL of final volume}$ was the volume of final extract submitted to analysis, and
 $\text{GC dilution factor}$ was the magnitude of dilution required to bracket the response of the sample within the standard curve responses. When the sample requires no dilution, the GC dilution factor equals 1.

The percent recovery for fortified control samples was calculated as follows:

$$\% \text{Recovery} = \frac{(\text{ppm of analyte measured in fortified sample} - \text{ppm measured in control sample}) \times 100\%}{\text{ppm of analyte added to sample}}$$

4.4.2 Examples

The following example calculations were for a DPX-KN128/IN-KN127 fortification at 1.0 ppm (Spike 43 of Set #16; Figure 7) extracted 17 December 1996 and analyzed 19 December 1996. Please refer to Table 2 for the results of the ppm-of-analyte and percent-recovery calculations.

$$\begin{aligned} \mu\text{g analyte/mL} &= ((113,833 \text{ area counts} - 0 \text{ area counts}) - (-3164 \text{ area counts})) \\ &\quad \div (2,537,000 \text{ area counts} \times \text{mL}/\mu\text{g}) \\ &= 116,997 \text{ area counts} \div (2,537,000 \text{ area counts} \times \text{mL}/\mu\text{g}) \\ &= 0.0461 \mu\text{g/mL} \end{aligned}$$

$$\begin{aligned} \text{ppm of analyte} &= (0.0461 \mu\text{g/mL} \times 36 \text{ mL} \times 1.0 \text{ mL} \times 20) \div (10.0 \text{ g} \times 3.6 \text{ mL}) \\ &= (33.19 \mu\text{g} \times \text{mL}) \div (36 \text{ g} \times \text{mL}) \end{aligned}$$

= 0.922 $\mu\text{g/g}$

= 0.92 ppm

% Recovery = $((0.92 \text{ ppm} - 0 \text{ ppm}) \times 100\%) \div 1.0 \text{ ppm}$
= 92%

5.3 *Modifications or Special Precautions*

It is imperative that a clean separation of phases be produced in the raw extract solvent (ACN:water). If this is not achieved, subsequent steps in the clean-up and analysis may be seriously compromised.

5.4 *Method Ruggedness*

5.4.1 *Stability*

The method appears to perform consistently over a long duration. The method has been successfully used for data generation purposes in a field dissipation program conducted in the U.S. since 1995 (DuPont Study Nos. AMR 3294-95 and AMR 3402-95) and is being used in a dissipation study in Canada (DuPont Study No. CAN-96-902).

5.4.2 *Specificity/Potential Interference*

5.4.2.a *Interference from glassware and reagents, matrices*

The method is a single-ion monitoring (SIM) method employing mass-spectral detection. Therefore, it has considerable inherent specificity. However, this method, as are other methods employing such detection, is subject to matrix enhancement phenomenon; hence, the use of standards in soil extract. A closely related soil (adjacent to treated field) could be used as a control in this assay. Since the test substance does seem to remain located in the upper 1 cm (approx.) of the soil surface, a control could be collected by carefully removing that upper surface and taking soil in the 5- to 15-cm horizon.

5.4.2.b *Interference from other pesticides*

This method was not tested for interference by other pesticides.

5.4.3 *Confirmatory Method*

This method has inherent confirmation of structure due to the detection by mass spectrometry. Additional ions or an entire mass spectrum could be monitored for structural identification. No confirmatory analysis by another technique should be necessary.

5.4.4 *Second Lab Tryout*

This method has not been tested in a second lab due to resource availability and timing issues.

APPENDIX 1 SILYLATION OF GLASSWARE

Silylation of Glassware

Purpose:

Silylation is a process used to chemically treat glassware in order to prevent or minimize binding of analyte residues to the glass surface.

Safety Considerations:

This procedure must be done in a fume hood. The person performing this procedure must wear heavy latex gloves. The silylating reagent, dimethyldichlorosilane, must NOT come in contact with water, otherwise chlorine gas and hydrogen chloride gas, both severely toxic, will be produced.

Procedure:

1. Prepare 100 mL of a 5% (v/v) solution of dimethyldichlorosilane (DMDCS) in hexane.

Add hexane (95 mL) to a glass stoppered glass container (approximately 200 mL volume). Slowly add 5 mL of DMDCS. Stopper the container and invert to mix.

Larger volumes may be prepared using the proportions discussed above. However, it is prudent to prepare amounts that will be entirely used in order to avoid storage and disposal of excess solution.

2. Pour a small amount of the DMDCS solution into the glassware to be treated. Rotate the glassware to thoroughly coat the inside surfaces. Pour excess solution into the next piece of glassware to be treated.
3. Allow the treated glassware to dry (approximately 20 minutes). Rinse with de-ionized water, then acetone. Again allow to dry.
4. Glassware is now ready for use.

- Notes:
- Any glassware that is cleaned with a brush after it has been silylanized, must be re-treated.
 - Store pure DMDCS at room temperature.
 - 5% solutions of DMDCS in hexane are stable for 5 days when stored well-stoppered at room temperature.

APPENDIX 2 QUALITY CONTROL PROCEDURE FOR SILICA SPE CARTRIDGES

Quality Control for SPE Cartridges

Silica Cartridges

DPX-KN128/IN-KN127 & DPX-JT333

Add 100 μL of 10 $\mu\text{g}/\text{mL}$ combined DPX-JW062 and DPX-JT333 standard (in ethyl acetate) to a 50-mL silanized centrifuge tube. Gently evaporate to dryness with nitrogen, and add 4 mL of hexane. Follow Steps 7 through 14 of Section 4.2.10.2 in procedure. Final concentration should be 1.0 $\mu\text{g}/\text{mL}$ for each analyte.

IN-KG433

Add 100 μL of 10 $\mu\text{g}/\text{mL}$ IN-KG433 standard (in ethyl acetate) to a 50-mL silanized centrifuge tube. Gently evaporate to dryness with nitrogen and add 4 mL of 5% ethyl acetate/95% toluene. Follow Steps 2 through 8 of Section 4.2.10.3 in procedure. Final concentration is 2.0 $\mu\text{g}/\text{mL}$.