

Analytical Methods TO-14a and TO-15: What are the Differences?

Determining the presence of volatile organic compounds in air can be complex given the many choices available. Many options exist for sample collection (sorbent tubes, bags, filters, and canisters), as well as, a variety of analytical techniques. To assist in providing technical assistance with the options for the organic compounds EPA has available a compendia of methods in the *Compendium of Methods for the Determination of Toxic Organic Compounds in Ambient Air -- Second Edition* (EPA/625/R-96/010b, January 1999)¹. Specifically this guidance document includes two methods using the specially-treated canister: Method TO-14a and Method TO-15. Method TO-15 is a new method available in the second edition of the compendium and Method TO-14a is a revision to Method TO-14 available from the 1987 compendium. Although Methods TO-14a and TO-15 are similar, there are differences that may impact the outcome of the desired results.

For sampling, Method TO-15 and Method TO-14a are identical; therefore, the analyte list and detection limit become the deciding factors when selecting the method of choice.

Method Specifications	TO-14a	TO-15
Non-polar VOCs (e.g., toluene, benzene)	√	√
Polar VOCs (e.g., methanol, alcohols, ketones)		√
GC/MS instrumentation	√	√
Sample collection by prepared canister (holding time = 30d)	√	√
Sample collection by sorbent tube	√	√
Water management techniques (avoid loss of polar compds)		√
Enhanced provisions for quality control		√
Method performance criteria		√
Selected Ion Monitoring (SIM)		√
Specific Cleaning procedures		√
Air sample concentrated onto solid sorbent trap	√	√
Use of other detectors for GC (e.g., GC/MD)	√	
Detection Limit 0.2 – 25 ppbV	√	√

¹ The *Compendium of Methods for the Determination of Toxic Organic Compounds in Ambient Air* and updates can be obtained from the EPA's OAQPS Technology Transfer Network website at <http://www.epa.gov/ttn/amtic/airtox.html>.

Method TO-15

Method TO-15 is applicable to a subset of the 97 volatile organic compounds (VOCs) that are from the list of 189 hazardous air pollutants (HAPs) included in Title III of the Clean Air Act Amendments. By employing water management steps, this method reduces the losses of water-soluble VOCs. Because of this, the analyte list includes polar and non-polar VOCs. Table 1 of the compendium provides the list of VOCs with a cross-reference to Method TO-14a and the CLP-SOW.

The analytical instrumentation for Method TO-15 is the high resolution gas chromatograph (GC) coupled to a mass spectrometer (MS). This instrument can be operated in a continuously scanning mode (SCAN mode) or using select-ion monitoring mode (SIM). Method TO-15 does not provide the option for a different detector such as FID, ECD, PID, or a multi-detector arrangement as specified in TO-14a. However, the GC/MS provides a more definitive identification technique and, with the SIM option, the lower detection limits can be achieved.

As written, Method TO-15 provides specific method performance criteria to allow for alternative techniques and provides inherent quality control with criteria for internal standards, system verifications, method blanks, and compound identification. With the quality control tools more assurances of good data are possible.

In summary, Method TO-15 provides information on alternative water management systems that allow for the analysis of polar compounds, has a more complete quality control section, provides performance criteria for any monitoring technique employed, and provides guidance specifically directed at compound identification by mass spectrometry.

Select Method TO-15 when:

- An expanded analyte list is needed; especially polar compounds
- Need for special canisters considerations
- Need for specified quality control tools; method performance criteria for acceptance of data

Some difficulties do exist when using Method TO15. Calibration difficulties do exist for some of the compounds listed in Table 1, such as, formaldehyde and diazomethane. In addition, some compounds may be subject to interferences or contamination.

Method TO-14a

Method TO-14a is a revised method for TO-14 and the target analyte list is provided in Table 1 of the method. The revision addresses technological updates and provides more explanations where necessary. The use of pressurized and sub-atmospheric pressure canisters is described for sample collection.

The analytical instrumentation for TO-14 is the high-resolution gas chromatograph (GC) coupled to one or more appropriate GC detectors. These non-specific detectors include, but are not limited to, the nitrogen-phosphorous detector (NPD), the flame ionization detector (FID), the electron capture detector (ECD), and the photo-ionization detector (PID). The specific detector includes MS operating in the SCAN mode or SIM mode. The non-specific detectors can provide lower detection limits but mis-identification can occur due to the reliance on the peak assignment in the chromatogram. Since most toxic VOCs in urban air are of low concentration in complex mixtures the use of the specific detectors is highly recommended by Method TO-14a.

In Method TO-14a, the Nafion® permeable membrane dryer is used to remove the water vapor from the air sample but this may cause the loss of polar organic compounds. For analysis of polar organic compounds Method TO-15 is recommended.

Select Method TO-14a when:

- Screening of analytes (can use non-specific detectors which may be cheaper)
- When only interested in a sub-set of analytes (non-polar compounds only)

TABLE 1. VOLATILE ORGANIC COMPOUNDS ON THE TITLE III CLEAN AIR AMENDMENT LIST-- MEMBERSHIP IN COMPENDIUM METHOD TO-14A LIST AND THE SOW-CLP LIST OF VOCs

Compound	CAS No.	TO-14A	TO-15	CLP-SOW ²
1,1-Dimethylhydrazine; C ₂ H ₈ N ₂	57-14-7		X	
1,1,2-Trichloroethane; C ₂ H ₃ Cl ₃	79-00-5	X	X	X
1,2-Propyleneimine (2-methylaziridine); C ₃ H ₇ N	75-55-8		X	
1,2,4-Trichlorobenzene; C ₆ H ₃ Cl ₃	120-82-1	X	X	X
1,2-Dibromo-3-chloropropane; C ₃ H ₅ Br ₂ Cl	96-12-8		X	
1,2-Epoxybutane (1,2-butylene oxide); C ₄ H ₈ O	106-88-7		X	
1,3-Butadiene ; C ₄ H ₆	106-99-0		X	X
1,3-Dichloropropene; C ₃ H ₄ Cl ₂ (cis)	542-75-6	X	X	X
1,3-Propane sultone; C ₃ H ₆ O ₃ S	1120-71-4		X	
1,4-Dioxane (1,4-diethylene oxide); C ₄ H ₈ O ₂	123-91-1		X	
1,4-Dichlorobenzene (p-); C ₆ H ₄ Cl ₂	106-46-7	X	X	X
2,2,4-Trimethyl pentane C ₈ H ₁₈	540-84-1		X	
2-Nitropropane; C ₃ H ₇ NO ₂	79-46-9		X	
Acetaldehyde (ethanal); C ₂ H ₄ O	75-07-0		X	
Acetonitrile (cyanomethane); C ₂ H ₃ N	75-05-8		X	X
Acetophenone; C ₈ H ₈ O	98-86-2		X	
Acrolein (2-propenal); C ₃ H ₄ O	107-02-8		X	X
Acrylamide; C ₃ H ₅ NO	79-06-1		X	
Acrylic acid; C ₃ H ₄ O ₂	79-10-7		X	
Acrylonitrile (2-propenenitrile); C ₃ H ₃ N	107-13-1	X	X	
Allyl chloride (3-chloropropene); C ₃ H ₅ Cl	107-05-1	X	X	X
Aniline (aminobenzene); C ₆ H ₇ N	62-53-3		X	
Benzene; C ₆ H ₆	71-43-2	X	X	X
Benzyl chloride (a-chlorotoluene); C ₇ H ₇ Cl	100-44-7	X	X	X
Beta-Propiolactone; C ₃ H ₄ O ₂	57-57-8		X	
Bis(2-Chloroethyl)ether; C ₄ H ₈ Cl ₂ O	111-44-4		X	
Bis(chloromethyl) ether; C ₂ H ₄ Cl ₂ O	542-88-1		X	
Bromoform (tribromomethane); CHBr ₃	75-25-2		X	
Carbon disulfide; CS ₂	75-15-0		X	
Carbon tetrachloride; CCl ₄	56-23-5	X	X	X
Carbonyl sulfide; COS	463-58-1		X	
Catechol (o-hydroxyphenol); C ₆ H ₆ O ₂	120-80-9		X	
Catechol (o-hydroxyphenol); C ₆ H ₆ O ₂	120-80-9		X	
Chloroacetic acid; C ₂ H ₃ ClO ₂	79-11-8		X	
Chlorobenzene; C ₆ H ₅ Cl	108-90-7	X	X	X
Chloroform; CHCl ₃	67-66-3	X	X	X
Chloromethyl methyl ether; C ₂ H ₅ ClO	107-30-2		X	
Chloroprene (2-chloro-1,3-butadiene); C ₄ H ₅ Cl	126-99-8		X	
Cresylic acid (cresol isomer mixture); C ₇ H ₈ O	1319-77-3		X	
Cumene (isopropylbenzene); C ₉ H ₁₂	98-82-8		X	
Diazomethane; CH ₂ N ₂	334-88-3		X	
Diethyl sulfate; C ₄ H ₁₀ O ₄ S	64-67-5		X	
Dimethyl sulfate; C ₂ H ₆ O ₄ S	77-78-1		X	
Dimethylcarbamyl chloride; C ₃ H ₆ ClNO	79-44-7		X	
Epichlorohydrin (1-chloro-2,3-epoxy propane); C ₃ H ₅ ClO	106-89-8		X	
Ethyl acrylate; C ₅ H ₈ O ₂	140-88-5		X	
Ethyl carbamate (urethane); C ₃ H ₇ NO ₂	51-79-6		X	
Ethyl chloride (chloroethane); C ₂ H ₅ Cl	75-00-3	X	X	X

² Statement of Work (SOW) for the Analysis of Air Toxics at Superfund Sites, Winberry W.T. Jr., et.al. USEPA, OSW, CLP, Washington D.C., Draft Report June 1990.

TABLE 1. VOLATILE ORGANIC COMPOUNDS ON THE TITLE III CLEAN AIR AMENDMENT LIST-- MEMBERSHIP IN COMPENDIUM METHOD TO-14A LIST AND THE SOW-CLP LIST OF VOCs

Compound	CAS No.	TO-14A	TO-15	CLP-SOW ²
Ethylbenzene; C ₈ H ₁₀	100-41-4	X	X	X
Ethylene dibromide (1,2-dibromoethane); C ₂ H ₄ Br ₂	106-93-4	X	X	X
Ethylene dichloride (1,2-dichloroethane); C ₂ H ₄ Cl ₂	107-06-2	X	X	X
Ethylene oxide; C ₂ H ₄ O	75-21-8		X	
Ethyleneimine (aziridine); C ₂ H ₅ N	151-56-4		X	
Ethylidene dichloride (1,1-dichloroethane); C ₂ H ₄ Cl ₂	75-34-3	X	X	
Formaldehyde ; CH ₂ O	50-00-0		X	
Hexachlorobutadiene; C ₄ Cl ₆	87-68-3	X	X	X
Hexachloroethane; C ₂ Cl ₆	67-72-1		X	
Hexane; C ₆ H ₁₄	110-54-3	X	X	
Isophorone; C ₉ H ₁₄ O	78-59-1		X	
1,1,2,2-Tetrachloroethane; C ₂ H ₂ Cl ₄	79-34-5	X	X	X
Methanol; CH ₄ O	67-56-1		X	X
Methyl bromide (bromomethane); CH ₃ Br	74-83-9	X	X	X
Methyl chloride (chloromethane); CH ₃ Cl	74-87-3	X	X	X
Methyl chloroform (1,1,1-trichloroethane); C ₂ H ₃ Cl ₃	71-55-6	X	X	X
Methyl ethyl ketone (2-butanone); C ₄ H ₈ O	78-93-3		X	X
Methyl iodide (iodomethane); CH ₃ I	74-88-4		X	
Methyl isobutyl ketone (hexone); C ₆ H ₁₂ O	108-10-1		X	
Methyl isocyanate; C ₂ H ₃ NO	624-83-9		X	
Methyl methacrylate; C ₅ H ₈ O ₂	80-62-6		X	
Methyl tert-butyl ether; C ₅ H ₁₂ O	1634-04-4		X	
Methylene chloride; CH ₂ Cl ₂	75-09-2	X	X	X
Methylhydrazine; CH ₆ N ₂	60-34-4		X	
m-Xylene; C ₈ H ₁₀	108-38-3	X	X	X
N,N-Dimethylaniline; C ₈ H ₁₁ N	121-69-7		X	
N,N-Dimethylformamide; C ₃ H ₇ NO	68-12-2		X	
Nitrobenzene; C ₆ H ₅ NO ₂	98-95-3		X	
N-Nitrosodimethylamine; C ₂ H ₆ N ₂ O	62-75-9		X	
N-Nitrosomorpholine; C ₄ H ₈ N ₂ O ₂	59-89-2		X	
N-Nitroso-N-methylurea; C ₂ H ₅ N ₃ O ₂	684-93-5		X	
o-Cresol; C ₇ H ₈ O	95-48-7		X	
o-Xylene; C ₈ H ₁₀	95-47-6	X	X	X
Phenol; C ₆ H ₆ O	108-95-2		X	
Phosgene; CCl ₂ O	75-44-5		X	
Propionaldehyde; C ₂ H ₅ CHO	123-38-6		X	
Propylene dichloride (1,2-dichloropropane); C ₃ H ₆ Cl ₂	78-87-5	X	X	X
Propylene oxide; C ₃ H ₆ O	75-56-9		X	
p-Xylene; C ₈ H ₁₀	106-42-3	X	X	X
Styrene oxide; C ₈ H ₈ O	96-09-3		X	
Styrene; C ₈ H ₈	100-42-5	X	X	X
Tetrachloroethylene; C ₂ Cl ₄	127-18-4	X	X	X
Toluene; C ₇ H ₈	108-88-3	X	X	X
Trichloroethylene; C ₂ HCl ₃	79-01-6	X	X	X
Triethylamine; C ₆ H ₁₅ N	121-44-8		X	
Vinyl acetate; C ₄ H ₆ O ₂	108-05-4		X	X
Vinyl bromide (bromoethene); C ₂ H ₃ Br	593-60-2		X	
Vinyl chloride (chloroethene); C ₂ H ₃ Cl	75-01-4	X	X	X
Vinylidene chloride (1,1-dichloroethylene); C ₂ H ₂ Cl ₂	75-35-4	X	X	X
Xylenes (isomer & mixtures); C ₈ H ₁₀	1330-20-7	X	X	X