An aerial photograph of a reservoir showing a significant green algal bloom. The water is a deep green color, and several black cylindrical buoys are visible in the water. The shoreline is lined with dense green trees and vegetation. A road or highway is visible in the bottom right corner.


Microcystins Analysis Methods

April 28, 2016

Heather Raymond
Ohio Harmful Algal Bloom Coordinator

Microcystins Testing

No "Perfect" Analytical Method for Detecting TOTAL Microcystins



- 1 - D-Alanine
- 2 - Variable L-amino acid
- 3 - D-Methylaspartic acid
- 4 - Variable L-amino acid
- 5 - 3-amino-9-methoxy-2,6,8-trimethyl-10-phenyldeca-4,6-dienoic acid (Adda)
- 6 - D-Glutamic acid
- 7 - N-Methyldehydroalanine

- Over 140 Microcystin Variants
- Standards Not Available for Majority

Microcystins Testing - ELISA

- **Enzyme-Linked ImmunoSorbent Assay (ELISA) Microcystin-ADDA Method** (detection of antigen using an antibody)
 - Measures Total Microcystins (all variants/congeners, based on ADDA)
 - Highly Selective/Specific (for ADDA)
 - Certified by USEPA (ETV Program)
 - Moderately sensitive (RL: 0.30ug/L)
 - Suitable for raw & finished water (complex matrices)
 - Quick (four hours), useful for operational adjustments
 - Relatively inexpensive
 - Does not require high end equipment or expertise to run (can be used in water system lab)
 - Does not require pre-concentration solid phase extraction (SPE) step

 - Does not provide concentrations of specific microcystin variants
 - Is an indirect measure of the toxin



Ohio EPA DES Method 701.0

Total Microcystins – ADDA by ELISA

- Developed in consultation with USEPA, PWSs, and National Experts.
- Helps ensure consistent sample handling, preparation, and application of analytical method.
 - Finished water samples and treatment train samples that are subjected to an oxidant must be quenched upon collection with 10 mg of sodium thiosulfate per 100 ml of sample.
 - At Lab, sample pH must be adjusted within the range of 5-11.
 - At Lab, all samples are subjected to three freeze/thaw cycles to lyse (break apart) cyanobacteria cells and release toxins.
- Labs must demonstrate they can achieve an acceptable level of precision and accuracy.
- Ohio EPA conducts site visits at labs performing analysis; provides acceptance letters and certification
- Compliance microcystins monitoring under 3745-90-03 must use this method or “another method accepted by the director in writing.”

A Standard Method Increases Consistency & Confidence
in Results



High Performance Liquid Chromatography (HPLC) – Ultraviolet (UV) or Photodiode Array (PDA)

- HPLC-UV/PDA
 - Liquid Chromatography separates components
 - Microcystins have UV absorption maxima at 238 nm
 - Non-selective detector; co-eluting interferences prevent accurate identification of components and quantitation (potential for false positives and false negatives)
 - Less expensive than mass spectrometry.
 - Less sensitive than mass spec (average LOQ ~ 0.3 µg/L)

Adapted From USEPA



Liquid Chromatography(LC) –Tandem Mass Spectrometry (MS/MS)

- LC/MS/MS
 - Highly specific identification of components (based on standards)
 - MS can identify a component in the presence of co-eluting interferences but quantitation may be compromised
 - Presence of co-eluting interferences can act to suppress or enhance response resulting in analytical bias
 - Sensitive (LOQ ~ 0.02 µg/L)
 - Expensive and requires highly skilled analysts
 - USEPA Method 544
 - Standard Method- includes QA/QC protocols and reduces variability in results between labs
 - Limited to 6 microcystin variants and finished water only

Adapted From USEPA



LC/MS/MS

Additional Considerations

- **Microcystins**
 - “Weak” product ion abundance limits sensitivity
 - Limited based on available standards
 - Requires preconcentration with SPE to augment sensitivity (LOQs $\leq 0.02 \mu\text{g/L}$)
 - Preconcentrates NOM too
- **Anatoxin-a and Cylindrospermopsin**
 - Abundant product ion responses
 - Direct injection (no SPE required to improve sensitivity – LOQs: anatoxin-a $\sim 0.02 \mu\text{g/L}$; cylindrospermopsin $\sim 0.06 \mu\text{g/L}$)

Source: USEPA



LC-MS/MS MMPB Method


- MMPB (2-methyl-3(methoxy)-4-phenylbutyric acid) method analyzes the chemically cleaved Adda group common to all microcystin variants
 - Measures total microcystins (all variants)
 - Quick (~2 hours, does not require freeze/thaw or sonication)
 - Sensitive (0.05 ug/L)
 - Suitable for raw and finished water
 - Does not require standards for individual variants
 - Utilizes 4 PB internal standard
-
- Does not provide data on individual congeners
 - Requires oxidation step
 - Potential for detection of microcystins disinfection byproducts

Toxicon 104 (2015) 91-101 (Foss & Aabel): Using the MMPB technique to confirm microcystin concentrations in water measured by ELISA and HPLC (UV, MS, MS/MS)



LC-UV/PDA & LC-MS Scan


- Uses two LC-based methods in tandem to independently confirm presence of microcystins
 - Can detect microcystin variants without standards
 - No standard methods, expensive, requires complex data-interpretation, time-consuming



Analytical Method Comparison & Microcystin Variant Evaluation


- 11 Sites: 4 Up-ground Reservoirs, 2 In-stream Reservoirs, 2 Lake Erie locations, 2 Canal-feeder Lakes, and 1 River Source.
- 22 Samples from 2014 Selected to Help Evaluate Spatial and Temporal Variability Within Source Waters
- Variety of Cyanobacteria Genera Represented
- Each Sample Analyzed Using 5 Separate Analytical Methods

Results of Method Comparison



* LC-UV data presented does not include false-positives that were eliminated from total (Based on lack of confirmation with LC-MS methods).
Sample # 14 was non-detect using LC-UV.

Results of LC-MS/MS MMPB and Individual Variant Analysis Compared to ELISA



Spatial and Temporal Variability in Microcystin Variants

Lake Erie Microcystin Variants					
MC-Variant	Site 1 8/25/14	Site 2 8/4/14	Site 2 8/18/14	Site 2 9/29/14	Site 2 10/14/14
MC-RR	2.1	20	10	5.5	8.5
MC-YR	0.6	6	5	1.2	2.5
MC-LR	2.9	16	10	5.5	6.1
MC-WR	0.6		3-9	0.2-0.6	0.2-0.6
MC-LY			2-6		
8.7 min 1043.5 m/z			10-30	3.6	

Inland Lake Microcystin Variants (Planktothrix)			
MC-Variant	Site 1	Site 2	Site 2
	6/16/14	6/16/14	9/2/14
[DAsp3] MC-RR	5.3	6.1	17.5
[Dha7] MC-LR	1.1	1.4	1.5
MC-YR	0.2-0.6	0.2-0.6	1.2
MC-RR		0.1-0.3	

Inland Lake Microcystin Variants (Mixed Bloom)				
MC-Variant	Site 1	Site 2	Site 2	Site 3
	6/18/14	6/18/14	7/9/14	6/30/14
[Dha7] MC-RR	2.9	3-9	1.0	0.08
MC-RR	1.4	39	1.0	0.01-0.03
MC-YR	1.1	15	1.0	
MC-LR	4.0	67	2.4	0.55
[DAsp3] MC-LR	0.6	18	0.4	0.03
[Dha7] MC-LR	3.6		1.0	0.05
MC-WR	0.2-0.6		0.2-0.6	
MC-LA	0.2-0.6			
MC-LY	0.2-0.6	6	0.2-0.6	0.10

Key Findings

- LC-based Methods Confirmed ELISA Results
- 16 Different MC-variants were detected
- MC-LR was only detected at 5 of 11 sites (45%)
- Most common variants were: MC-YR, [Dha7] MC-LR and [DAsp3] MC-RR
- HPLC-PDA Methods Prone to Interference, Especially at Lower Concentrations
- 91% of samples had MC-variants not detectable by USEPA Method 544 (including dominant MC-variant in some samples)
- LC-MS/MS individual variant analysis under-reported total microcystins, based on MMPB and LC-UV/MS Scan Data
- Generally, the Dominant/Co-Dominant MC-Variants Did Not Vary Spatially nor Temporally
- Secondary & Minor MC-Variants Did Vary.

ELISA MC-ADDA Matrix Interference Studies

Treatment Chemical	Microcystins – ADDA ELISA Assay Tolerance (< / =)
Sodium Carbonate (Soda Ash)	≤25 gpg
Sodium Hexametaphosphate	≤250 ppm
Sodium Silicofluoride	≤10 ppm
Aluminum Sulfate ¹	≤100 gpg (with pH adjustment within assay tolerance)
Calcium Oxide (Lime) ¹	≤2000 gpg (with pH adjustment to within assay tolerance)
Potassium Permanganate ²	≤10 ppm (with quenching using 1 mg sodium thiosulfate per 1 ml sample)
Sodium Chlorite ²	≤10 ppm (with quenching using 1 mg sodium thiosulfate per 1 ml sample)
Carbon ³	≤2 ppm with filtering at time of sampling

¹ Natural pH of solution outside assay tolerance, Chemical tolerance levels determined after pH adjustment.

² Oxidizers degrade microcystins, tolerance determined after quenching.

³ Tolerance level due to effect of carbon on toxin, not assay performance.



ELISA MC-ADDA Matrix Interference Studies

Lisa Kamp, et. al, 2016. *The effects of water sample treatment, preparation, and storage prior to cyanotoxin analysis for cylindrospermopsin, microcystin and Saxitoxin.* Chemico-Biological Interactions.

Studies by USEPA as part of ELISA MC-ADDA Method Development for UCMR 4:

- Storage Stability – Holding Times
- Sample Preservation and Container Studies
- Matrix Interference Studies
 - Microcystins Variant Fortified Sample Studies (finished water, raw water, reagent water with chemical addition, etc.)
 - Dilution Experiments (real world raw/finished water samples)
- USEPA Method Validation & Interlab Validation

LC-MS/MS MMPB Method Evaluation

- Potential concern regarding detection of microcystins disinfection byproducts
- ELISA MC-ADDA does not detect microcystins disinfection byproducts



Ohio EPA Method Comparison Study- Round 2

- 17 Raw and treatment train samples from 2015 submitted for method comparison
- Additional MC variant standards available
- Preliminary LC-MS/MS results indicate MC-HtYR, MC-HiLR, MC-LF, MC-LW, MC-RR, MC-YR, MC-LR, MC-WR, and MC-LY present in Lake Erie Samples
- 100% of samples had MC variants not included in USEPA Method 544
- Next Steps: Split sample interlab method comparison study



2015 PWS HAB Response Strategy

Analytical Methods

	Microcystins ($\mu\text{g/L}$)	Cylindro- spermopsin ($\mu\text{g/L}$)	Saxitoxins ($\mu\text{g/L}$)	Anatoxin-a ($\mu\text{g/L}$)
Surveillance sampling	ELISA (MC-ADDA)	ELISA	ELISA	LC-MS/MS
Repeat sampling in response to a finished water detection	ELISA (MC-ADDA)	LC-MS/MS	LC-MS/MS	LC-MS/MS

ELISA: Enzyme-Linked Immunosorbent Assay

LC-MS/MS: Liquid Chromatography followed by
tandem Mass Spectrometry



Data Reporting Considerations

- Specify Analysis Method
- If reporting “total microcystins” using LC-MS/MS, include number of variants analyzing for in comment field
- For ELISA results greater than calibration curve, dilute and reanalyze sample or qualify result as estimate (or report using “>” symbol)

Questions?

Heather.Raymond@epa.ohio.gov

(614) 644-2752

<http://epa.ohio.gov/ddagw/HAB.aspx>

